


Surgical Saw Blades

CATALOGUE

Accessible and Reliable







Saw Blades
Compatible with
DEPUY SYNTHES[®]



Brand Code **220**




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 130 090 089	S 220 130 090 089	13 mm	90 mm	0,89 mm	0,89 mm	Straight
220 130 090 090	S 220 130 090 090	13 mm	90 mm	0,90 mm	0,90 mm	Straight
220 130 090 100	S 220 130 090 100	13 mm	90 mm	1,00 mm	1,00 mm	Straight
220 130 090 119	S 220 130 090 119	13 mm	90 mm	1,19 mm	1,19 mm	Straight
220 130 090 120	S 220 130 090 120	13 mm	90 mm	1,20 mm	1,20 mm	Straight
220 130 090 127	S 220 130 090 127	13 mm	90 mm	1,27 mm	1,27 mm	Straight
220 130 090 130	S 220 130 090 130	13 mm	90 mm	1,30 mm	1,30 mm	Straight
220 130 090 137	S 220 130 090 137	13 mm	90 mm	1,37 mm	1,37 mm	Straight
220 130 090 140	S 220 130 090 140	13 mm	90 mm	1,40 mm	1,40 mm	Straight
220 130 090 147	S 220 130 090 147	13 mm	90 mm	1,47 mm	1,47 mm	Straight
220 130 090 150	S 220 130 090 150	13 mm	90 mm	1,50 mm	1,50 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 190 090 089	S 220 190 090 089	19 mm	90 mm	0,89 mm	0,89 mm	Straight
220 190 090 090	S 220 190 090 090	19 mm	90 mm	0,90 mm	0,90 mm	Straight
220 190 090 100	S 220 190 090 100	19 mm	90 mm	1,00 mm	1,00 mm	Straight
220 190 090 119	S 220 190 090 119	19 mm	90 mm	1,19 mm	1,19 mm	Straight
220 190 090 120	S 220 190 090 120	19 mm	90 mm	1,20 mm	1,20 mm	Straight
220 190 090 127	S 220 190 090 127	19 mm	90 mm	1,27 mm	1,27 mm	Straight
220 190 090 130	S 220 190 090 130	19 mm	90 mm	1,30 mm	1,30 mm	Straight
220 190 090 137	S 220 190 090 137	19 mm	90 mm	1,37 mm	1,37 mm	Straight
220 190 090 140	S 220 190 090 140	19 mm	90 mm	1,40 mm	1,40 mm	Straight
220 190 090 147	S 220 190 090 147	19 mm	90 mm	1,47 mm	1,47 mm	Straight
220 190 090 150	S 220 190 090 150	19 mm	90 mm	1,50 mm	1,50 mm	Straight

Compatible Connections;


Compact AirDrive II + Power Drive 511.800, 511.801, Air Power Line 511.610, Battery Power Line 530.610, TR5 05.001.224, 05.001.223, 05.001.240, Battery Power Line II 530.710




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 190 100 089	S 220 190 100 089	19 mm	100 mm	0,89 mm	0,89 mm	Straight
220 190 100 090	S 220 190 100 090	19 mm	100 mm	0,90 mm	0,90 mm	Straight
220 190 100 100	S 220 190 100 100	19 mm	100 mm	1,00 mm	1,00 mm	Straight
220 190 100 119	S 220 190 100 119	19 mm	100 mm	1,19 mm	1,19 mm	Straight
220 190 100 120	S 220 190 100 120	19 mm	100 mm	1,20 mm	1,20 mm	Straight
220 190 100 127	S 220 190 100 127	19 mm	100 mm	1,27 mm	1,27 mm	Straight
220 190 100 130	S 220 190 100 130	19 mm	100 mm	1,30 mm	1,30 mm	Straight
220 190 100 137	S 220 190 100 137	19 mm	100 mm	1,37 mm	1,37 mm	Straight
220 190 100 140	S 220 190 100 140	19 mm	100 mm	1,40 mm	1,40 mm	Straight
220 190 100 147	S 220 190 100 147	19 mm	100 mm	1,47 mm	1,47 mm	Straight
220 190 100 150	S 220 190 100 150	19 mm	100 mm	1,50 mm	1,50 mm	Straight



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 195 090 089	S 220 195 090 089	19,5 mm	90 mm	0,89 mm	0,89 mm	Straight
220 195 090 090	S 220 195 090 090	19,5 mm	90 mm	0,90 mm	0,90 mm	Straight
220 195 090 100	S 220 195 090 100	19,5 mm	90 mm	1,00 mm	1,00 mm	Straight
220 195 090 119	S 220 195 090 119	19,5 mm	90 mm	1,19 mm	1,19 mm	Straight
220 195 090 120	S 220 195 090 120	19,5 mm	90 mm	1,20 mm	1,20 mm	Straight
220 195 090 127	S 220 195 090 127	19,5 mm	90 mm	1,27 mm	1,27 mm	Straight
220 195 090 130	S 220 195 090 130	19,5 mm	90 mm	1,30 mm	1,30 mm	Straight
220 195 090 137	S 220 195 090 137	19,5 mm	90 mm	1,37 mm	1,37 mm	Straight
220 195 090 140	S 220 195 090 140	19,5 mm	90 mm	1,40 mm	1,40 mm	Straight
220 195 090 147	S 220 195 090 147	19,5 mm	90 mm	1,47 mm	1,47 mm	Straight
220 195 090 150	S 220 195 090 150	19,5 mm	90 mm	1,50 mm	1,50 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 222 090 089	S 220 222 090 089	22,2 mm	90 mm	0,89 mm	0,89 mm	Straight
220 222 090 090	S 220 222 090 090	22,2 mm	90 mm	0,90 mm	0,90 mm	Straight
220 222 090 100	S 220 222 090 100	22,2 mm	90 mm	1,00 mm	1,00 mm	Straight
220 222 090 119	S 220 222 090 119	22,2 mm	90 mm	1,19 mm	1,19 mm	Straight
220 222 090 120	S 220 222 090 120	22,2 mm	90 mm	1,20 mm	1,20 mm	Straight
220 222 090 127	S 220 222 090 127	22,2 mm	90 mm	1,27 mm	1,27 mm	Straight
220 222 090 130	S 220 222 090 130	22,2 mm	90 mm	1,30 mm	1,30 mm	Straight
220 222 090 137	S 220 222 090 137	22,2 mm	90 mm	1,37 mm	1,37 mm	Straight
220 222 090 140	S 220 222 090 140	22,2 mm	90 mm	1,40 mm	1,40 mm	Straight
220 222 090 147	S 220 222 090 147	22,2 mm	90 mm	1,47 mm	1,47 mm	Straight
220 222 090 150	S 220 222 090 150	22,2 mm	90 mm	1,50 mm	1,50 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 250 090 089	S 220 250 090 089	25 mm	90 mm	0,89 mm	0,89 mm	Straight
220 250 090 090	S 220 250 090 090	25 mm	90 mm	0,90 mm	0,90 mm	Straight
220 250 090 100	S 220 250 090 100	25 mm	90 mm	1,00 mm	1,00 mm	Straight
220 250 090 119	S 220 250 090 119	25 mm	90 mm	1,19 mm	1,19 mm	Straight
220 250 090 120	S 220 250 090 120	25 mm	90 mm	1,20 mm	1,20 mm	Straight
220 250 090 127	S 220 250 090 127	25 mm	90 mm	1,27 mm	1,27 mm	Straight
220 250 090 130	S 220 250 090 130	25 mm	90 mm	1,30 mm	1,30 mm	Straight
220 250 090 137	S 220 250 090 137	25 mm	90 mm	1,37 mm	1,37 mm	Straight
220 250 090 140	S 220 250 090 140	25 mm	90 mm	1,40 mm	1,40 mm	Straight
220 250 090 147	S 220 250 090 147	25 mm	90 mm	1,47 mm	1,47 mm	Straight
220 250 090 150	S 220 250 090 150	25 mm	90 mm	1,50 mm	1,50 mm	Straight

Compatible Connections:


Compact AirDrive II + Power Drive 511.800, 511.801, Air Power Line 511.610, Battery Power Line 530.610, TR5 05.001.224, 05.001.223, 05.001.240, Battery Power Line II 530.710




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 254 090 089	S 220 254 090 089	25,4 mm	90 mm	0,89 mm	0,89 mm	Straight
220 254 090 090	S 220 254 090 090	25,4 mm	90 mm	0,90 mm	0,90 mm	Straight
220 254 090 100	S 220 254 090 100	25,4 mm	90 mm	1,00 mm	1,00 mm	Straight
220 254 090 119	S 220 254 090 119	25,4 mm	90 mm	1,19 mm	1,19 mm	Straight
220 254 090 120	S 220 254 090 120	25,4 mm	90 mm	1,20 mm	1,20 mm	Straight
220 254 090 127	S 220 254 090 127	25,4 mm	90 mm	1,27 mm	1,27 mm	Straight
220 254 090 130	S 220 254 090 130	25,4 mm	90 mm	1,30 mm	1,30 mm	Straight
220 254 090 137	S 220 254 090 137	25,4 mm	90 mm	1,37 mm	1,37 mm	Straight
220 254 090 140	S 220 254 090 140	25,4 mm	90 mm	1,40 mm	1,40 mm	Straight
220 254 090 147	S 220 254 090 147	25,4 mm	90 mm	1,47 mm	1,47 mm	Straight
220 254 090 150	S 220 254 090 150	25,4 mm	90 mm	1,50 mm	1,50 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 125 090 127	S 220 125 090 127	12,5 mm	90 mm	1,27 mm	1,27 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 190 095 120	S 220 190 095 120	19 mm	95 mm	1,20 mm	1,20 mm	Straight
220 190 095 127	S 220 190 095 127	19 mm	95 mm	1,27 mm	1,27 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 230 090 127	S 220 230 090 127	23 mm	90 mm	1,27 mm	1,27 mm	Straight



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 250 095 119	S 220 250 095 119	25 mm	95 mm	1,19 mm	1,19 mm	Straight
220 250 095 127	S 220 250 095 127	25 mm	95 mm	1,27 mm	1,27 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 090 025 060	S 220 090 025 060	9 mm	25 mm	0,60 mm	0,60 mm	Straight
X 220 090 025 060	SX 220 090 025 060	9 mm	25 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 070 025 060	S 220 070 025 060	7 mm	25 mm	0,60 mm	0,60 mm	Straight

Compatible Connections:


Compact AirDrive II + Power Drive 511.800, 511.801, Air Power Line 511.610, Battery Power Line 530.610, TRS 05.001.224, 05.001.223, 05.001.240, Battery Power Line II 530.710




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 100 026 060	S 220 100 026 060	10 mm	26 mm	0,60 mm	0,60 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 140 026 060	S 220 140 026 060	14 mm	26 mm	0,60 mm	0,60 mm	Straight
X 220 140 026 060	SX 220 140 026 060	14 mm	26 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 140 051 060	S 220 140 051 060	14 mm	51 mm	0,60 mm	0,60 mm	Straight
X 220 140 051 060	SX 220 140 051 060	14 mm	51 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 141 050 060	S 220 141 050 060	14,1 mm	50 mm	0,60 mm	0,60 mm	Straight



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 130 045 127	S 220 130 045 127	13 mm	45 mm	1,27 mm	1,27 m	Straight



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 130 060 090	S 220 130 060 090	13 mm	60 mm	0,90 mm	0,90 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 150 060 090	S 220 150 060 090	15 mm	60 mm	0,90 mm	0,90 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 270 051 060	S 220 270 051 060	27 mm	51 mm	0,60 mm	0,60 mm	Straight
X 220 270 051 060	SX 220 270 051 060	27 mm	51 mm	0,60 mm	0,80 mm	Cross

Compatible Connections:


Compact AirDrive II + Power Drive 511.800, 511.801, Air Power Line 511.610, Battery Power Line 530.610, TRS 05.001.224, 05.001.223, 05.001.240, Battery Power Line II 530.710




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 140 020 060	S 220 140 020 060	14 mm	20 mm	0,60 mm	0,60 mm	Straight
X 220 140 020 060	SX 220 140 020 060	14 mm	20 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 480 070 060	S 220 480 070 060	48 mm	70 mm	0,60 mm	0,60 mm	Straight
X 220 480 070 060	SX 220 480 070 060	48 mm	70 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 490 070 090	S 220 490 070 090	49 mm	70 mm	0,60 mm	0,60 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 125 077 090	S 220 125 077 090	12,5 mm	77 mm	0,90 mm	0,90 mm	Straight
220 125 077 100	S 220 125 077 100	12,5 mm	77 mm	1,00 mm	1,00 mm	Straight
220 125 077 120	S 220 125 077 120	12,5 mm	77 mm	1,20 mm	1,20 mm	Straight



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 125 089 090	S 220 125 089 090	12,5 mm	89 mm	0,90 mm	0,90 mm	Straight
220 125 089 100	S 220 125 089 100	12,5 mm	89 mm	1,00 mm	1,00 mm	Straight
220 125 089 120	S 220 125 089 120	12,5 mm	89 mm	1,20 mm	1,20 mm	Straight



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 120 067 090	S 220 120 067 090	12 mm	67 mm	0,90 mm	0,90 mm	Straight
220 120 067 100	S 220 120 067 100	12 mm	67 mm	1,00 mm	1,00 mm	Straight
220 120 067 120	S 220 120 067 120	12 mm	67 mm	1,20 mm	1,20 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 130 068 090	S 220 130 068 090	13 mm	68 mm	0,90 mm	0,90 mm	Straight
220 130 068 100	S 220 130 068 100	13 mm	68 mm	1,00 mm	1,00 mm	Straight
220 130 068 120	S 220 130 068 120	13 mm	68 mm	1,20 mm	1,20 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 050 018 060	S 220 050 018 060	5 mm	18 mm	0,60 mm	0,60 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 060 015 040	S 220 060 015 040	6 mm	15 mm	0,40 mm	0,40 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 060 030 060	S 220 060 030 060	6 mm	30 mm	0,60 mm	0,60 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 070 015 060	S 220 070 015 060	7 mm	15 mm	0,60 mm	0,60 mm	Straight
X 220 070 015 060	SX 220 070 015 060	7 mm	15 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 080 022 060	S 220 080 022 060	8 mm	22 mm	0,60 mm	0,60 mm	Straight
X 220 080 022 060	SX 220 080 022 060	8 mm	22 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 100 025 060	S 220 100 025 060	10 mm	25 mm	0,60 mm	0,60 mm	Straight
X 220 100 025 060	SX 220 100 025 060	10 mm	25 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
X 220 100 027 040	SX 220 100 027 040	10 mm	27 mm	0,40 mm	0,60 mm	Cross
220 100 027 060	S 220 100 027 060	10 mm	27 mm	0,60 mm	0,60 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 100 030 060	S 220 100 030 060	10 mm	30 mm	0,60 mm	0,60 mm	Straight
X 220 100 030 060	SX 220 100 030 060	10 mm	30 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 100 031 060	S 220 100 031 060	10 mm	31 mm	0,60 mm	0,60 mm	Straight
X 220 100 031 060	SX 220 100 031 060	10 mm	31 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 100 049 060	S 220 100 049 060	10 mm	49 mm	0,60 mm	0,60 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
X 220 100 050 040	SX 220 100 050 040	10 mm	50 mm	0,40 mm	0,60 mm	Cross
220 100 050 060	S 220 100 050 060	10 mm	50 mm	0,60 mm	0,60 mm	Straight
X 220 100 050 060	SX 220 100 050 060	10 mm	50 mm	0,60 mm	0,80 mm	Cross



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 120 022 060	S 220 120 022 060	12 mm	22 mm	0,60 mm	0,60 mm	Straight
X 220 120 022 060	SX 220 120 022 060	12 mm	22 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 120 040 060	S 220 120 040 060	12 mm	40 mm	0,60 mm	0,60 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 140 025 060	S 220 140 025 060	14 mm	25 mm	0,60 mm	0,60 mm	Straight
X 220 140 025 060	SX 220 140 025 060	14 mm	25 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 140 027 060	S 220 140 027 060	14 mm	27 mm	0,60 mm	0,60 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 140 050 060	S 220 140 050 060	14 mm	50 mm	0,60 mm	0,60 mm	Straight



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 145 050 060	S 220 145 050 060	14,5 mm	50 mm	0,60 mm	0,60 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 198 050 060	S 220 198 050 060	19,8 mm	50 mm	0,60 mm	0,60 mm	Straight



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 200 049 060	S 220 200 049 060	20 mm	49 mm	0,60 mm	0,60 mm	Straight



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
X 220 200 050 040	SX 220 200 050 040	20 mm	50 mm	0,40 mm	0,60 mm	Cross



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 270 050 060	S 220 270 050 060	27 mm	50 mm	0,60 mm	0,60 mm	Straight
X 220 270 050 060	SX 220 270 050 060	27 mm	50 mm	0,60 mm	0,80 mm	Cross



ACF Medikal Urunler Mak. San. ve Tic. Ltd. Sti.

Guzelyurt Mah. 5754 Sok. No: 3/A
Yunusemre / Manisa, Turkey, 45030

+90 236 236 28 20

info@acf.com.tr

www.acf.com.tr