


# Surgical Saw Blades

## CATALOGUE

Accessible and Reliable






**Saw Blades  
Compatible with  
DE SOUTTER®**




**Brand Code**


**200 / 205**



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 130 085 127	S 200 130 085 127	13 mm	85 mm	1,27 mm	1,27 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 130 090 089	S 200 130 090 089	13 mm	90 mm	0,89 mm	0,89 mm	Straight
200 130 090 090	S 200 130 090 090	13 mm	90 mm	0,90 mm	0,90 mm	Straight
200 130 090 100	S 200 130 090 100	13 mm	90 mm	1,00 mm	1,00 mm	Straight
200 130 090 119	S 200 130 090 119	13 mm	90 mm	1,19 mm	1,19 mm	Straight
200 130 090 120	S 200 130 090 120	13 mm	90 mm	1,20 mm	1,20 mm	Straight
200 130 090 127	S 200 130 090 127	13 mm	90 mm	1,27 mm	1,27 mm	Straight
200 130 090 130	S 200 130 090 130	13 mm	90 mm	1,30 mm	1,30 mm	Straight
200 130 090 137	S 200 130 090 137	13 mm	90 mm	1,37 mm	1,37 mm	Straight
200 130 090 140	S 200 130 090 140	13 mm	90 mm	1,40 mm	1,40 mm	Straight
200 130 090 147	S 200 130 090 147	13 mm	90 mm	1,47 mm	1,47 mm	Straight
200 130 090 150	S 200 130 090 150	13 mm	90 mm	1,50 mm	1,50 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 190 085 090	S 200 190 085 090	19 mm	85 mm	0,90 mm	0,90 mm	Straight

## Compatible Connections;


MDX Multi Drive KX 600, MDX Multi Drive KDX-600,  
MPX Multi Drive SX 500, MPX Multi Drive DPX-676,  
MCX Multi Drive DCX 500, Orthodrive DBK-700




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 190 090 089	S 200 190 090 089	19 mm	90 mm	0,89 mm	0,89 mm	Straight
200 190 090 090	S 200 190 090 090	19 mm	90 mm	0,90 mm	0,90 mm	Straight
200 190 090 100	S 200 190 090 100	19 mm	90 mm	1,00 mm	1,00 mm	Straight
200 190 090 119	S 200 190 090 119	19 mm	90 mm	1,19 mm	1,19 mm	Straight
200 190 090 120	S 200 190 090 120	19 mm	90 mm	1,20 mm	1,20 mm	Straight
200 190 090 127	S 200 190 090 127	19 mm	90 mm	1,27 mm	1,27 mm	Straight
200 190 090 130	S 200 190 090 130	19 mm	90 mm	1,30 mm	1,30 mm	Straight
200 190 090 137	S 200 190 090 137	19 mm	90 mm	1,37 mm	1,37 mm	Straight
200 190 090 140	S 200 190 090 140	19 mm	90 mm	1,40 mm	1,40 mm	Straight
200 190 090 147	S 200 190 090 147	19 mm	90 mm	1,47 mm	1,47 mm	Straight
200 190 090 150	S 200 190 090 150	19 mm	90 mm	1,50 mm	1,50 mm	Straight



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 190 100 089	S 200 190 100 089	19 mm	100 mm	0,89 mm	0,89 mm	Straight
200 190 100 090	S 200 190 100 090	19 mm	100 mm	0,90 mm	0,90 mm	Straight
200 190 100 100	S 200 190 100 100	19 mm	100 mm	1,00 mm	1,00 mm	Straight
200 190 100 119	S 200 190 100 119	19 mm	100 mm	1,19 mm	1,19 mm	Straight
200 190 100 120	S 200 190 100 120	19 mm	100 mm	1,20 mm	1,20 mm	Straight
200 190 100 127	S 200 190 100 127	19 mm	100 mm	1,27 mm	1,27 mm	Straight
200 190 100 130	S 200 190 100 130	19 mm	100 mm	1,30 mm	1,30 mm	Straight
200 190 100 137	S 200 190 100 137	19 mm	100 mm	1,37 mm	1,37 mm	Straight
200 190 100 140	S 200 190 100 140	19 mm	100 mm	1,40 mm	1,40 mm	Straight
200 190 100 147	S 200 190 100 147	19 mm	100 mm	1,47 mm	1,47 mm	Straight
200 190 100 150	S 200 190 100 150	19 mm	100 mm	1,50 mm	1,50 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 195 090 089	S 200 195 090 089	19,5 mm	90 mm	0,89 mm	0,89 mm	Straight
200 195 090 090	S 200 195 090 090	19,5 mm	90 mm	0,90 mm	0,90 mm	Straight
200 195 090 100	S 200 195 090 100	19,5 mm	90 mm	1,00 mm	1,00 mm	Straight
200 195 090 119	S 200 195 090 119	19,5 mm	90 mm	1,19 mm	1,19 mm	Straight
200 195 090 120	S 200 195 090 120	19,5 mm	90 mm	1,20 mm	1,20 mm	Straight
200 195 090 127	S 200 195 090 127	19,5 mm	90 mm	1,27 mm	1,27 mm	Straight
200 195 090 130	S 200 195 090 130	19,5 mm	90 mm	1,30 mm	1,30 mm	Straight
200 195 090 137	S 200 195 090 137	19,5 mm	90 mm	1,37 mm	1,37 mm	Straight
200 195 090 140	S 200 195 090 140	19,5 mm	90 mm	1,40 mm	1,40 mm	Straight
200 195 090 147	S 200 195 090 147	19,5 mm	90 mm	1,47 mm	1,47 mm	Straight
200 195 090 150	S 200 195 090 150	19,5 mm	90 mm	1,50 mm	1,50 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 220 090 089	S 200 220 090 089	22 mm	90 mm	0,89 mm	0,89 mm	Straight
200 220 090 090	S 200 220 090 090	22 mm	90 mm	0,90 mm	0,90 mm	Straight
200 220 090 100	S 200 220 090 100	22 mm	90 mm	1,00 mm	1,00 mm	Straight
200 220 090 119	S 200 220 090 119	22 mm	90 mm	1,19 mm	1,19 mm	Straight
200 220 090 120	S 200 220 090 120	22 mm	90 mm	1,20 mm	1,20 mm	Straight
200 220 090 127	S 200 220 090 127	22 mm	90 mm	1,27 mm	1,27 mm	Straight
200 220 090 130	S 200 220 090 130	22 mm	90 mm	1,30 mm	1,30 mm	Straight
200 220 090 137	S 200 220 090 137	22 mm	90 mm	1,37 mm	1,37 mm	Straight
200 220 090 140	S 200 220 090 140	22 mm	90 mm	1,40 mm	1,40 mm	Straight
200 220 090 147	S 200 220 090 147	22 mm	90 mm	1,47 mm	1,47 mm	Straight
200 220 090 150	S 200 220 090 150	22 mm	90 mm	1,50 mm	1,50 mm	Straight

**Compatible Connections;**


MDX Multi Drive KX 600, MDX Multi Drive KDX-600,  
MPX Multi Drive SX 500, MPX Multi Drive DPX-676,  
MCX Multi Drive DCX 500, Orthodrive DBK-700




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 250 090 089	S 200 250 090 089	25 mm	90 mm	0,89 mm	0,89 mm	Straight
200 250 090 090	S 200 250 090 090	25 mm	90 mm	0,90 mm	0,90 mm	Straight
200 250 090 100	S 200 250 090 100	25 mm	90 mm	1,00 mm	1,00 mm	Straight
200 250 090 119	S 200 250 090 119	25 mm	90 mm	1,19 mm	1,19 mm	Straight
200 250 090 120	S 200 250 090 120	25 mm	90 mm	1,20 mm	1,20 mm	Straight
200 250 090 127	S 200 250 090 127	25 mm	90 mm	1,27 mm	1,27 mm	Straight
200 250 090 130	S 200 250 090 130	25 mm	90 mm	1,30 mm	1,30 mm	Straight
200 250 090 137	S 200 250 090 137	25 mm	90 mm	1,37 mm	1,37 mm	Straight
200 250 090 140	S 200 250 090 140	25 mm	90 mm	1,40 mm	1,40 mm	Straight
200 250 090 147	S 200 250 090 147	25 mm	90 mm	1,47 mm	1,47 mm	Straight
200 250 090 150	S 200 250 090 150	25 mm	90 mm	1,50 mm	1,50 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 254 090 089	S 200 254 090 089	25,4 mm	90 mm	0,89 mm	0,89 mm	Straight
200 254 090 090	S 200 254 090 090	25,4 mm	90 mm	0,90 mm	0,90 mm	Straight
200 254 090 100	S 200 254 090 100	25,4 mm	90 mm	1,00 mm	1,00 mm	Straight
200 254 090 119	S 200 254 090 119	25,4 mm	90 mm	1,19 mm	1,19 mm	Straight
200 254 090 120	S 200 254 090 120	25,4 mm	90 mm	1,20 mm	1,20 mm	Straight
200 254 090 127	S 200 254 090 127	25,4 mm	90 mm	1,27 mm	1,27 mm	Straight
200 254 090 130	S 200 254 090 130	25,4 mm	90 mm	1,30 mm	1,30 mm	Straight
200 254 090 137	S 200 254 090 137	25,4 mm	90 mm	1,37 mm	1,37 mm	Straight
200 254 090 140	S 200 254 090 140	25,4 mm	90 mm	1,40 mm	1,40 mm	Straight
200 254 090 147	S 200 254 090 147	25,4 mm	90 mm	1,47 mm	1,47 mm	Straight
200 254 090 150	S 200 254 090 150	25,4 mm	90 mm	1,50 mm	1,50 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 250 085 127	S 200 250 085 127	25 mm	85 mm	1,27 mm	1,27 mm	Straight



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 070 020 080	S 200 070 020 080	7 mm	20 mm	0,80 mm	0,80 mm	Straight



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 090 025 060	S 200 090 025 060	9 mm	25 mm	0,60 mm	0,60 mm	Straight
X 200 090 025 060	SX 200 090 025 060	9 mm	25 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 120 030 060	S 200 120 030 060	12 mm	30 mm	0,60 mm	0,60 mm	Straight

## Compatible Connections;


MDX Multi Drive KX 600, MDX Multi Drive KDX-600,  
MPX Multi Drive SX 500, MPX Multi Drive DPX-676,  
MCX Multi Drive DCX 500, Orthodrive DBK-700



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 120 025 060	S 200 120 025 060	12 mm	25 mm	0,60 mm	0,60 mm	Straight
X 200 120 025 060	SX 200 120 025 060	12 mm	25 mm	0,60 mm	0,80 mm	Cross



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 120 060 060	S 200 120 060 060	12 mm	60 mm	0,60 mm	0,60 mm	Straight
X 200 120 060 060	SX 200 120 060 060	12 mm	60 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 150 025 060	S 200 150 025 060	15 mm	25 mm	0,60 mm	0,60 mm	Straight
X 200 150 025 060	SX 200 150 025 060	15 mm	25 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 150 048 060	S 200 150 048 060	15 mm	48 mm	0,60 mm	0,60 mm	Straight
X 200 150 048 060	SX 200 150 048 060	15 mm	48 mm	0,60 mm	0,80 mm	Cross







NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 200 055 060	S 200 200 055 060	20 mm	55 mm	0,60 mm	0,60 mm	Straight
X 200 200 055 060	SX 200 200 055 060	20 mm	55 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 270 048 060	S 200 270 048 060	27 mm	48 mm	0,60 mm	0,60 mm	Straight
X 200 270 048 060	SX 200 270 048 060	27 mm	48 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 270 068 060	S 200 270 068 060	27 mm	68 mm	0,60 mm	0,60 mm	Straight
X 200 270 068 060	SX 200 270 068 060	27 mm	68 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 090 055 060	S 200 090 055 060	9 mm	55 mm	0,60 mm	0,60 mm	Straight
X 200 090 055 060	SX 200 090 055 060	9 mm	55 mm	0,60 mm	0,80 mm	Cross
200 090 055 090	S 200 090 055 090	9 mm	55 mm	0,90 mm	0,90 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 090 080 060	S 200 090 080 060	9 mm	80 mm	0,60 mm	0,60 mm	Straight
X 200 090 080 060	SX 200 090 080 060	9 mm	80 mm	0,60 mm	0,80 mm	Cross
200 090 080 090	S 200 090 080 090	9 mm	80 mm	0,90 mm	0,90 mm	Straight
200 090 080 100	S 200 090 080 100	9 mm	80 mm	1,00 mm	1,00 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 090 089 100	S 200 090 089 100	9 mm	89 mm	1,00 mm	1,00 mm	Straight
200 090 089 120	S 200 090 089 120	9 mm	89 mm	1,20 mm	1,20 mm	Straight
200 090 089 130	S 200 090 089 130	9 mm	89 mm	1,30 mm	1,30 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 100 089 120	S 200 100 089 120	10 mm	89 mm	1,20 mm	1,20 mm	Straight
200 100 089 130	S 200 100 089 130	10 mm	89 mm	1,30 mm	1,30 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 120 070 100	S 200 120 070 100	12 mm	70 mm	1,00 mm	1,00 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 120 082 120	S 200 120 082 120	12 mm	82 mm	1,20 mm	1,20 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 055 019 060	S 200 055 019 060	5,5 mm	19 mm	0,60 mm	0,60 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 055 023 060	S 200 055 023 060	5,5 mm	23 mm	0,60 mm	0,60 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 055 026 060	S 200 055 026 060	5,5 mm	26 mm	0,60 mm	0,60 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 055 034 060	S 200 055 034 060	5,5 mm	34 mm	0,60 mm	0,60 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 060 030 060	S 200 060 030 060	6 mm	30 mm	0,60 mm	0,60 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 070 016 060	S 200 070 016 060	7 mm	16 mm	0,60 mm	0,60 mm	Straight
X 200 070 016 060	SX 200 070 016 060	7 mm	16 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 070 020 060	S 200 070 020 060	7 mm	20 mm	0,60 mm	0,60 mm	Straight
X 200 070 020 060	SX 200 070 020 060	7 mm	20 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 070 022 060	S 200 070 022 060	7 mm	22 mm	0,60 mm	0,60 mm	Straight
X 200 070 022 060	SX 200 070 022 060	7 mm	22 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 070 031 060	S 200 070 031 060	7 mm	31 mm	0,60 mm	0,60 mm	Straight
X 200 070 031 060	SX 200 070 031 060	7 mm	31 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 080 030 060	S 200 080 030 060	8 mm	30 mm	0,60 mm	0,60 mm	Straight
X 200 080 030 060	SX 200 080 030 060	8 mm	30 mm	0,60 mm	0,80 mm	Cross



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 090 022 060	S 200 090 022 060	9 mm	22 mm	0,60 mm	0,60 mm	Straight
X 200 090 022 060	SX 200 090 022 060	9 mm	22 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 090 027 060	S 200 090 027 060	9 mm	27 mm	0,60 mm	0,60 mm	Straight
X 200 090 027 060	SX 200 090 027 060	9 mm	27 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 090 030 060	S 200 090 030 060	9 mm	30 mm	0,60 mm	0,60 mm	Straight
X 200 090 030 060	SX 200 090 030 060	9 mm	30 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 090 032 060	S 200 090 032 060	9 mm	32 mm	0,60 mm	0,60 mm	Straight
X 200 090 032 060	SX 200 090 032 060	9 mm	32 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 090 036 060	S 200 090 036 060	9 mm	36 mm	0,60 mm	0,60 mm	Straight
X 200 090 036 060	SX 200 090 036 060	9 mm	36 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 090 039 060	S 200 090 039 060	9 mm	39 mm	0,60 mm	0,60 mm	Straight
X 200 090 039 060	SX 200 090 039 060	9 mm	39 mm	0,60 mm	0,80 mm	Cross



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 130 026 060	S 200 130 026 060	13 mm	26 mm	0,60 mm	0,60 mm	Straight
X 200 130 026 060	SX 200 130 026 060	13 mm	26 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 130 039 060	S 200 130 039 060	13 mm	39 mm	0,60 mm	0,60 mm	Straight
X 200 130 039 060	SX 200 130 039 060	13 mm	39 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 190 039 060	S 200 190 039 060	19 mm	39 mm	0,60 mm	0,60 mm	Straight
X 200 190 039 060	SX 200 190 039 060	19 mm	39 mm	0,60 mm	0,80 mm	Cross







NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
205 130 090 089	S 205 130 090 089	13 mm	90 mm	0,89 mm	0,89 mm	Straight
205 130 090 090	S 205 130 090 090	13 mm	90 mm	0,90 mm	0,90 mm	Straight
205 130 090 100	S 205 130 090 100	13 mm	90 mm	1,00 mm	1,00 mm	Straight
205 130 090 119	S 205 130 090 119	13 mm	90 mm	1,19 mm	1,19 mm	Straight
205 130 090 120	S 205 130 090 120	13 mm	90 mm	1,20 mm	1,20 mm	Straight
205 130 090 127	S 205 130 090 127	13 mm	90 mm	1,27 mm	1,27 mm	Straight
205 130 090 130	S 205 130 090 130	13 mm	90 mm	1,30 mm	1,30 mm	Straight
205 130 090 137	S 205 130 090 137	13 mm	90 mm	1,37 mm	1,37 mm	Straight
205 130 090 140	S 205 130 090 140	13 mm	90 mm	1,40 mm	1,40 mm	Straight
205 130 090 147	S 205 130 090 147	13 mm	90 mm	1,47 mm	1,47 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
205 190 090 089	S 205 190 090 089	19 mm	90 mm	0,89 mm	0,89 mm	Straight
205 190 090 090	S 205 190 090 090	19 mm	90 mm	0,90 mm	0,90 mm	Straight
205 190 090 100	S 205 190 090 100	19 mm	90 mm	1,00 mm	1,00 mm	Straight
205 190 090 119	S 205 190 090 119	19 mm	90 mm	1,19 mm	1,19 mm	Straight
205 190 090 120	S 205 190 090 120	19 mm	90 mm	1,20 mm	1,20 mm	Straight
205 190 090 127	S 205 190 090 127	19 mm	90 mm	1,27 mm	1,27 mm	Straight
205 190 090 130	S 205 190 090 130	19 mm	90 mm	1,30 mm	1,30 mm	Straight
205 190 090 137	S 205 190 090 137	19 mm	90 mm	1,37 mm	1,37 mm	Straight
205 190 090 140	S 205 190 090 140	19 mm	90 mm	1,40 mm	1,40 mm	Straight
205 190 090 147	S 205 190 090 147	19 mm	90 mm	1,47 mm	1,47 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
205 190 100 089	S 205 190 100 089	19 mm	100 mm	0,89 mm	0,89 mm	Straight
205 190 100 090	S 205 190 100 090	19 mm	100 mm	0,90 mm	0,90 mm	Straight
205 190 100 100	S 205 190 100 100	19 mm	100 mm	1,00 mm	1,00 mm	Straight
205 190 100 119	S 205 190 100 119	19 mm	100 mm	1,19 mm	1,19 mm	Straight
205 190 100 120	S 205 190 100 120	19 mm	100 mm	1,20 mm	1,20 mm	Straight
205 190 100 127	S 205 190 100 127	19 mm	100 mm	1,27 mm	1,27 mm	Straight
205 190 100 130	S 205 190 100 130	19 mm	100 mm	1,30 mm	1,30 mm	Straight
205 190 100 137	S 205 190 100 137	19 mm	100 mm	1,37 mm	1,37 mm	Straight
205 190 100 140	S 205 190 100 140	19 mm	100 mm	1,40 mm	1,40 mm	Straight
205 190 100 147	S 205 190 100 147	19 mm	100 mm	1,47 mm	1,47 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
205 195 090 089	S 205 195 090 089	19.5 mm	90 mm	0,89 mm	0,89 mm	Straight
205 195 090 090	S 205 195 090 090	19.5 mm	90 mm	0,90 mm	0,90 mm	Straight
205 195 090 100	S 205 195 090 100	19.5 mm	90 mm	1,00 mm	1,00 mm	Straight
205 195 090 119	S 205 195 090 119	19.5 mm	90 mm	1,19 mm	1,19 mm	Straight
205 195 090 120	S 205 195 090 120	19.5 mm	90 mm	1,20 mm	1,20 mm	Straight
205 195 090 127	S 205 195 090 127	19.5 mm	90 mm	1,27 mm	1,27 mm	Straight
205 195 090 130	S 205 195 090 130	19.5 mm	90 mm	1,30 mm	1,30 mm	Straight
205 195 090 137	S 205 195 090 137	19.5 mm	90 mm	1,37 mm	1,37 mm	Straight
205 195 090 140	S 205 195 090 140	19.5 mm	90 mm	1,40 mm	1,40 mm	Straight
205 195 090 147	S 205 195 090 147	19.5 mm	90 mm	1,47 mm	1,47 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
205 220 090 089	S 205 220 090 089	22 mm	90 mm	0,89 mm	0,89 mm	Straight
205 220 090 090	S 205 220 090 090	22 mm	90 mm	0,90 mm	0,90 mm	Straight
205 220 090 100	S 205 220 090 100	22 mm	90 mm	1,00 mm	1,00 mm	Straight
205 220 090 119	S 205 220 090 119	22 mm	90 mm	1,19 mm	1,19 mm	Straight
205 220 090 120	S 205 220 090 120	22 mm	90 mm	1,20 mm	1,20 mm	Straight
205 220 090 127	S 205 220 090 127	22 mm	90 mm	1,27 mm	1,27 mm	Straight
205 220 090 130	S 205 220 090 130	22 mm	90 mm	1,30 mm	1,30 mm	Straight
205 220 090 137	S 205 220 090 137	22 mm	90 mm	1,37 mm	1,37 mm	Straight
205 220 090 140	S 205 220 090 140	22 mm	90 mm	1,40 mm	1,40 mm	Straight
205 220 090 147	S 205 220 090 147	22 mm	90 mm	1,47 mm	1,47 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
205 223 090 089	S 205 223 090 089	22,3 mm	90 mm	0,89 mm	0,89 mm	Straight
205 223 090 090	S 205 223 090 090	22,3 mm	90 mm	0,90 mm	0,90 mm	Straight
205 223 090 100	S 205 223 090 100	22,3 mm	90 mm	1,00 mm	1,00 mm	Straight
205 223 090 119	S 205 223 090 119	22,3 mm	90 mm	1,19 mm	1,19 mm	Straight
205 223 090 120	S 205 223 090 120	22,3 mm	90 mm	1,20 mm	1,20 mm	Straight
205 223 090 127	S 205 223 090 127	22,3 mm	90 mm	1,27 mm	1,27 mm	Straight
205 223 090 130	S 205 223 090 130	22,3 mm	90 mm	1,30 mm	1,30 mm	Straight
205 223 090 137	S 205 223 090 137	22,3 mm	90 mm	1,37 mm	1,37 mm	Straight
205 223 090 140	S 205 223 090 140	22,3 mm	90 mm	1,40 mm	1,40 mm	Straight
205 223 090 147	S 205 223 090 147	22,3 mm	90 mm	1,47 mm	1,47 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
205 250 090 089	S 205 250 090 089	25 mm	90 mm	0,89 mm	0,89 mm	Straight
205 250 090 090	S 205 250 090 090	25 mm	90 mm	0,90 mm	0,90 mm	Straight
205 250 090 100	S 205 250 090 100	25 mm	90 mm	1,00 mm	1,00 mm	Straight
205 250 090 119	S 205 250 090 119	25 mm	90 mm	1,19 mm	1,19 mm	Straight
205 250 090 120	S 205 250 090 120	25 mm	90 mm	1,20 mm	1,20 mm	Straight
205 250 090 127	S 205 250 090 127	25 mm	90 mm	1,27 mm	1,27 mm	Straight
205 250 090 130	S 205 250 090 130	25 mm	90 mm	1,30 mm	1,30 mm	Straight
205 250 090 137	S 205 250 090 137	25 mm	90 mm	1,37 mm	1,37 mm	Straight
205 250 090 140	S 205 250 090 140	25 mm	90 mm	1,40 mm	1,40 mm	Straight
205 250 090 147	S 205 250 090 147	25 mm	90 mm	1,47 mm	1,47 mm	Straight



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
205 251 090 089	S 205 251 090 089	25,1 mm	90 mm	0,89 mm	0,89 mm	Straight
205 251 090 090	S 205 251 090 090	25,1 mm	90 mm	0,90 mm	0,90 mm	Straight
205 251 090 100	S 205 251 090 100	25,1 mm	90 mm	1,00 mm	1,00 mm	Straight
205 251 090 119	S 205 251 090 119	25,1 mm	90 mm	1,19 mm	1,19 mm	Straight
205 251 090 120	S 205 251 090 120	25,1 mm	90 mm	1,20 mm	1,20 mm	Straight
205 251 090 127	S 205 251 090 127	25,1 mm	90 mm	1,27 mm	1,27 mm	Straight
205 251 090 130	S 205 251 090 130	25,1 mm	90 mm	1,30 mm	1,30 mm	Straight
205 251 090 137	S 205 251 090 137	25,1 mm	90 mm	1,37 mm	1,37 mm	Straight
205 251 090 140	S 205 251 090 140	25,1 mm	90 mm	1,40 mm	1,40 mm	Straight
205 251 090 147	S 205 251 090 147	25,1 mm	90 mm	1,47 mm	1,47 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
205 254 090 089	S 205 254 090 089	25,4 mm	90 mm	0,89 mm	0,89 mm	Straight
205 254 090 090	S 205 254 090 090	25,4 mm	90 mm	0,90 mm	0,90 mm	Straight
205 254 090 100	S 205 254 090 100	25,4 mm	90 mm	1,00 mm	1,00 mm	Straight
205 254 090 119	S 205 254 090 119	25,4 mm	90 mm	1,19 mm	1,19 mm	Straight
205 254 090 120	S 205 254 090 120	25,4 mm	90 mm	1,20 mm	1,20 mm	Straight
205 254 090 127	S 205 254 090 127	25,4 mm	90 mm	1,27 mm	1,27 mm	Straight
205 254 090 130	S 205 254 090 130	25,4 mm	90 mm	1,30 mm	1,30 mm	Straight
205 254 090 137	S 205 254 090 137	25,4 mm	90 mm	1,37 mm	1,37 mm	Straight
205 254 090 140	S 205 254 090 140	25,4 mm	90 mm	1,40 mm	1,40 mm	Straight
205 254 090 147	S 205 254 090 147	25,4 mm	90 mm	1,47 mm	1,47 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
205 070 020 060	S 205 070 020 060	7 mm	20 mm	0,60 mm	0,60 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
205 090 025 060	S 205 090 025 060	9 mm	25 mm	0,60 mm	0,60 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
205 120 025 060	S 205 120 025 060	12 mm	25 mm	0,60 mm	0,60 mm	Straight



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
205 150 025 060	S 205 150 025 060	15 mm	25 mm	0,60 mm	0,60 mm	Straight



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
205 130 045 089	S 205 130 045 089	13 mm	45 mm	0,89 mm	0,89 mm	Straight
205 130 045 090	S 205 130 045 090	13 mm	45 mm	0,90 mm	0,90 mm	Straight
205 130 045 100	S 205 130 045 100	13 mm	45 mm	1,00 mm	1,00 mm	Straight
205 130 045 119	S 205 130 045 119	13 mm	45 mm	1,19 mm	1,19 mm	Straight
205 130 045 120	S 205 130 045 120	13 mm	45 mm	1,20 mm	1,20 mm	Straight
205 130 045 127	S 205 130 045 127	13 mm	45 mm	1,27 mm	1,27 mm	Straight
205 130 045 130	S 205 130 045 130	13 mm	45 mm	1,30 mm	1,30 mm	Straight
205 130 045 137	S 205 130 045 137	13 mm	45 mm	1,37 mm	1,37 mm	Straight
205 130 045 140	S 205 130 045 140	13 mm	45 mm	1,40 mm	1,40 mm	Straight
205 130 045 147	S 205 130 045 147	13 mm	45 mm	1,47 mm	1,47 mm	Straight



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
205 150 048 060	S 205 150 048 060	15 mm	48 mm	0,60 mm	0,60 mm	Straight



ACF Medikal Urunler Mak. San. ve Tic. Ltd. Sti.

Guzelyurt Mah. 5754 Sok. No: 3/A  
Yunusemre / Manisa, Turkey, 45030

+90 236 236 28 20

info@acf.com.tr

[www.acf.com.tr](http://www.acf.com.tr)