


Surgical Saw Blades

CATALOGUE

Accessible and Reliable






**Saw Blades
Compatible with
ZIMMER[®]**




Brand Code 220


For Reciprocating Saws
223




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 130 090 089	S 220 130 090 089	13 mm	90 mm	0,89 mm	0,89 mm	Straight
220 130 090 090	S 220 130 090 090	13 mm	90 mm	0,90 mm	0,90 mm	Straight
220 130 090 100	S 220 130 090 100	13 mm	90 mm	1,00 mm	1,00 mm	Straight
220 130 090 119	S 220 130 090 119	13 mm	90 mm	1,19 mm	1,19 mm	Straight
220 130 090 120	S 220 130 090 120	13 mm	90 mm	1,20 mm	1,20 mm	Straight
220 130 090 127	S 220 130 090 127	13 mm	90 mm	1,27 mm	1,27 mm	Straight
220 130 090 130	S 220 130 090 130	13 mm	90 mm	1,30 mm	1,30 mm	Straight
220 130 090 137	S 220 130 090 137	13 mm	90 mm	1,37 mm	1,37 mm	Straight
220 130 090 140	S 220 130 090 140	13 mm	90 mm	1,40 mm	1,40 mm	Straight
220 130 090 147	S 220 130 090 147	13 mm	90 mm	1,47 mm	1,47 mm	Straight
220 130 090 150	S 220 130 090 150	13 mm	90 mm	1,50 mm	1,50 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 190 090 089	S 220 190 090 089	19 mm	90 mm	0,89 mm	0,89 mm	Straight
220 190 090 090	S 220 190 090 090	19 mm	90 mm	0,90 mm	0,90 mm	Straight
220 190 090 100	S 220 190 090 100	19 mm	90 mm	1,00 mm	1,00 mm	Straight
220 190 090 119	S 220 190 090 119	19 mm	90 mm	1,19 mm	1,19 mm	Straight
220 190 090 120	S 220 190 090 120	19 mm	90 mm	1,20 mm	1,20 mm	Straight
220 190 090 127	S 220 190 090 127	19 mm	90 mm	1,27 mm	1,27 mm	Straight
220 190 090 130	S 220 190 090 130	19 mm	90 mm	1,30 mm	1,30 mm	Straight
220 190 090 137	S 220 190 090 137	19 mm	90 mm	1,37 mm	1,37 mm	Straight
220 190 090 140	S 220 190 090 140	19 mm	90 mm	1,40 mm	1,40 mm	Straight
220 190 090 147	S 220 190 090 147	19 mm	90 mm	1,47 mm	1,47 mm	Straight
220 190 090 150	S 220 190 090 150	19 mm	90 mm	1,50 mm	1,50 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 190 100 089	S 220 190 100 089	19 mm	100 mm	0,89 mm	0,89 mm	Straight
220 190 100 090	S 220 190 100 090	19 mm	100 mm	0,90 mm	0,90 mm	Straight
220 190 100 100	S 220 190 100 100	19 mm	100 mm	1,00 mm	1,00 mm	Straight
220 190 100 119	S 220 190 100 119	19 mm	100 mm	1,19 mm	1,19 mm	Straight
220 190 100 120	S 220 190 100 120	19 mm	100 mm	1,20 mm	1,20 mm	Straight
220 190 100 127	S 220 190 100 127	19 mm	100 mm	1,27 mm	1,27 mm	Straight
220 190 100 130	S 220 190 100 130	19 mm	100 mm	1,30 mm	1,30 mm	Straight
220 190 100 137	S 220 190 100 137	19 mm	100 mm	1,37 mm	1,37 mm	Straight
220 190 100 140	S 220 190 100 140	19 mm	100 mm	1,40 mm	1,40 mm	Straight
220 190 100 147	S 220 190 100 147	19 mm	100 mm	1,47 mm	1,47 mm	Straight
220 190 100 150	S 220 190 100 150	19 mm	100 mm	1,50 mm	1,50 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 195 090 089	S 220 195 090 089	19,5 mm	90 mm	0,89 mm	0,89 mm	Straight
220 195 090 090	S 220 195 090 090	19,5 mm	90 mm	0,90 mm	0,90 mm	Straight
220 195 090 100	S 220 195 090 100	19,5 mm	90 mm	1,00 mm	1,00 mm	Straight
220 195 090 119	S 220 195 090 119	19,5 mm	90 mm	1,19 mm	1,19 mm	Straight
220 195 090 120	S 220 195 090 120	19,5 mm	90 mm	1,20 mm	1,20 mm	Straight
220 195 090 127	S 220 195 090 127	19,5 mm	90 mm	1,27 mm	1,27 mm	Straight
220 195 090 130	S 220 195 090 130	19,5 mm	90 mm	1,30 mm	1,30 mm	Straight
220 195 090 137	S 220 195 090 137	19,5 mm	90 mm	1,37 mm	1,37 mm	Straight
220 195 090 140	S 220 195 090 140	19,5 mm	90 mm	1,40 mm	1,40 mm	Straight
220 195 090 147	S 220 195 090 147	19,5 mm	90 mm	1,47 mm	1,47 mm	Straight
220 195 090 150	S 220 195 090 150	19,5 mm	90 mm	1,50 mm	1,50 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 222 090 089	S 220 222 090 089	22,2 mm	90 mm	0,89 mm	0,89 mm	Straight
220 222 090 090	S 220 222 090 090	22,2 mm	90 mm	0,90 mm	0,90 mm	Straight
220 222 090 100	S 220 222 090 100	22,2 mm	90 mm	1,00 mm	1,00 mm	Straight
220 222 090 119	S 220 222 090 119	22,2 mm	90 mm	1,19 mm	1,19 mm	Straight
220 222 090 120	S 220 222 090 120	22,2 mm	90 mm	1,20 mm	1,20 mm	Straight
220 222 090 127	S 220 222 090 127	22,2 mm	90 mm	1,27 mm	1,27 mm	Straight
220 222 090 130	S 220 222 090 130	22,2 mm	90 mm	1,30 mm	1,30 mm	Straight
220 222 090 137	S 220 222 090 137	22,2 mm	90 mm	1,37 mm	1,37 mm	Straight
220 222 090 140	S 220 222 090 140	22,2 mm	90 mm	1,40 mm	1,40 mm	Straight
220 222 090 147	S 220 222 090 147	22,2 mm	90 mm	1,47 mm	1,47 mm	Straight
220 222 090 150	S 220 222 090 150	22,2 mm	90 mm	1,50 mm	1,50 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 250 090 089	S 220 250 090 089	25 mm	90 mm	0,89 mm	0,89 mm	Straight
220 250 090 090	S 220 250 090 090	25 mm	90 mm	0,90 mm	0,90 mm	Straight
220 250 090 100	S 220 250 090 100	25 mm	90 mm	1,00 mm	1,00 mm	Straight
220 250 090 119	S 220 250 090 119	25 mm	90 mm	1,19 mm	1,19 mm	Straight
220 250 090 120	S 220 250 090 120	25 mm	90 mm	1,20 mm	1,20 mm	Straight
220 250 090 127	S 220 250 090 127	25 mm	90 mm	1,27 mm	1,27 mm	Straight
220 250 090 130	S 220 250 090 130	25 mm	90 mm	1,30 mm	1,30 mm	Straight
220 250 090 137	S 220 250 090 137	25 mm	90 mm	1,37 mm	1,37 mm	Straight
220 250 090 140	S 220 250 090 140	25 mm	90 mm	1,40 mm	1,40 mm	Straight
220 250 090 147	S 220 250 090 147	25 mm	90 mm	1,47 mm	1,47 mm	Straight
220 250 090 150	S 220 250 090 150	25 mm	90 mm	1,50 mm	1,50 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 254 090 089	S 220 254 090 089	25,4 mm	90 mm	0,89 mm	0,89 mm	Straight
220 254 090 090	S 220 254 090 090	25,4 mm	90 mm	0,90 mm	0,90 mm	Straight
220 254 090 100	S 220 254 090 100	25,4 mm	90 mm	1,00 mm	1,00 mm	Straight
220 254 090 119	S 220 254 090 119	25,4 mm	90 mm	1,19 mm	1,19 mm	Straight
220 254 090 120	S 220 254 090 120	25,4 mm	90 mm	1,20 mm	1,20 mm	Straight
220 254 090 127	S 220 254 090 127	25,4 mm	90 mm	1,27 mm	1,27 mm	Straight
220 254 090 130	S 220 254 090 130	25,4 mm	90 mm	1,30 mm	1,30 mm	Straight
220 254 090 137	S 220 254 090 137	25,4 mm	90 mm	1,37 mm	1,37 mm	Straight
220 254 090 140	S 220 254 090 140	25,4 mm	90 mm	1,40 mm	1,40 mm	Straight
220 254 090 147	S 220 254 090 147	25,4 mm	90 mm	1,47 mm	1,47 mm	Straight
220 254 090 150	S 220 254 090 150	25,4 mm	90 mm	1,50 mm	1,50 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 125 090 127	S 220 125 090 127	12,5 mm	90 mm	1,27 mm	1,27 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 190 095 120	S 220 190 095 120	19 mm	95 mm	1,20 mm	1,20 mm	Straight
220 190 095 127	S 220 190 095 127	19 mm	95 mm	1,27 mm	1,27 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 230 090 127	S 220 230 090 127	23 mm	90 mm	1,27 mm	1,27 mm	Straight



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 250 095 119	S 220 250 095 119	25 mm	95 mm	1,19 mm	1,19 mm	Straight
220 250 095 127	S 220 250 095 127	25 mm	95 mm	1,27 mm	1,27 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 090 025 060	S 220 090 025 060	9 mm	25 mm	0,60 mm	0,60 mm	Straight
X 220 090 025 060	SX 220 090 025 060	9 mm	25 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 070 025 060	S 220 070 025 060	7 mm	25 mm	0,60 mm	0,60 mm	Straight

Zimmer®


Compatible Connections:
Zimmer Universal Power System,
89-8509-450-60



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 100 026 060	S 220 100 026 060	10 mm	26 mm	0,60 mm	0,60 mm	Straight



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 140 026 060	S 220 140 026 060	14 mm	26 mm	0,60 mm	0,60 mm	Straight
X 220 140 026 060	SX 220 140 026 060	14 mm	26 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 140 051 060	S 220 140 051 060	14 mm	51 mm	0,60 mm	0,60 mm	Straight
X 220 140 051 060	SX 220 140 051 060	14 mm	51 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 141 050 060	S 220 141 050 060	14,1 mm	50 mm	0,60 mm	0,60 mm	Straight



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 130 045 127	S 220 130 045 127	13 mm	45 mm	1,27 mm	1,27 m	Straight



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 130 060 090	S 220 130 060 090	13 mm	60 mm	0,90 mm	0,90 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 150 060 090	S 220 150 060 090	15 mm	60 mm	0,90 mm	0,90 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 270 051 060	S 220 270 051 060	27 mm	51 mm	0,60 mm	0,60 mm	Straight
X 220 270 051 060	SX 220 270 051 060	27 mm	51 mm	0,60 mm	0,80 mm	Cross

Zimmer®


Compatible Connections:
Zimmer Universal Power System,
89-8509-450-60




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 140 020 060	S 220 140 020 060	14 mm	20 mm	0,60 mm	0,60 mm	Straight
X 220 140 020 060	SX 220 140 020 060	14 mm	20 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 480 070 060	S 220 480 070 060	48 mm	70 mm	0,60 mm	0,60 mm	Straight
X 220 480 070 060	SX 220 480 070 060	48 mm	70 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 490 070 090	S 220 490 070 090	49 mm	70 mm	0,90 mm	0,90 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
223 125 062 090	S 223 125 062 090	12,5 mm	62 mm	0,90 mm	0,90 mm	Straight
223 125 062 100	S 223 125 062 100	12,5 mm	62 mm	1,00 mm	1,00 mm	Straight
223 125 062 120	S 223 125 062 120	12,5 mm	62 mm	1,20 mm	1,20 mm	Straight



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
223 125 073 090	S 223 125 073 090	12,5 mm	73 mm	0,90 mm	0,90 mm	Straight
223 125 073 100	S 223 125 073 100	12,5 mm	73 mm	1,00 mm	1,00 mm	Straight
223 125 073 120	S 223 125 073 120	12,5 mm	73 mm	1,20 mm	1,20 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
223 125 075 090	S 223 125 075 090	12,5 mm	75 mm	0,90 mm	0,90 mm	Straight
223 125 075 100	S 223 125 075 100	12,5 mm	75 mm	1,00 mm	1,00 mm	Straight
223 125 075 120	S 223 125 075 120	12,5 mm	75 mm	1,20 mm	1,20 mm	Straight



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
223 125 089 090	S 223 125 089 090	12,5 mm	89 mm	0,90 mm	0,90 mm	Straight
223 125 089 100	S 223 125 089 100	12,5 mm	89 mm	1,00 mm	1,00 mm	Straight
223 125 089 120	S 223 125 089 120	12,5 mm	89 mm	1,20 mm	1,20 mm	Straight

Zimmer®

Compatible Connections:
Zimmer Universal Power System,
89-8509-451-20



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
223 135 076 090	S 223 135 076 090	13,5 mm	76 mm	0,90 mm	0,90 mm	Straight
223 135 076 100	S 223 135 076 100	13,5 mm	76 mm	1,00 mm	1,00 mm	Straight
223 135 076 120	S 223 135 076 120	13,5 mm	76 mm	1,20 mm	1,20 mm	Straight



ACF Medikal Urunler Mak. San. ve Tic. Ltd. Sti.

Guzelyurt Mah. 5754 Sok. No: 3/A
Yunusemre / Manisa, Turkey, 45030

+90 236 236 28 20

info@acf.com.tr

www.acf.com.tr