


Surgical Saw Blades

CATALOGUE

Accessible and Reliable



Contents

Where Are Surgical Saw Blades Used?	3
Choosing Saw Blades	4
Warning & Cautions / Copyright Information	5
Saw Blades Compatible with ACF [®] (Brand Code 360)	6
Saw Blades Compatible with Microaire [®] (Brand Code 100)	18
Saw Blades Compatible with Stryker [®] (Brand Code 120)	33
Saw Blades Compatible with Conmed Linvatec [®] (Brand Code 140 / 100)	52
Saw Blades Compatible with De Soutter [®] (Brand Code 200/205)	68
Saw Blades Compatible with Smith & Nephew Dyonics [®] (Brand Code 210)	90
Saw Blades Compatible with Depuy Synthes [®] (Brand Code 220)	93
Saw Blades Compatible with Zimmer [®] (Brand Code 220 / 223)	109
Saw Blades Compatible with Brasseler [®] / TAVA [®] (Brand Code 260)	120
Saw Blades Compatible with B Braun Aesculap [®] (Brand Code 320 / 330)	127
Saw Blades Compatible with Bojin [®] (Brand Code 180)	148
Saw Blades Compatible with Chinese Power Tools [®] (Brand Code 300)	153

Where Are Surgical Saw Blades Used?

1 Arthroplasty

Knee, Hip, Shoulder Replacement Surgeries



2 Trauma

Plate, Nails



3 ACL(Anterior Cruciate Ligament), Hand, Foot and Wrist Surgeries

Finger, Foot Implants, Ankle Arthroplasty, Wrist Joint Replacement




4 Plastic Bone Surgery

Maxillofacial Trauma



5 Amputation Surgery

Surgical removal of part of the body, such as an arm or leg






Choosing Saw Blades

Speaking in Reference Code

Most simple reference code system. Understandable and ready to order.

REF 360 190 090 120



- Brand Code** : This code shows the power tool brand that the blades are compatible.
- Cut Edge** : The width of the blade measured from the outside base of the two extreme teeth.
- Cut Depth** : The actual depth attainable, not the overall blade length, while making a cut.
- Thickness** : The distance between the top and the bottom of the stainless steel.

360 150 048 060 : This means the saw blade is **non-sterile**.

S 360 150 048 060 : This means the saw blade is **sterile**. (with gamma irradiation)

X 360 150 048 060 : This means the saw blade is **non-sterile** and it has **cross teeth**. (material and cutting thickness are different)

SX 360 150 048 060 : This means the saw blade is **sterile** and it has **cross teeth**.

Warning & Cautions / Copyright Information

Warning & Cautions

- *Read instructions for use accompany saw blades.*
- *Reusing is not recommended. Multiple use will cause wear on the blade may cause breakage.*
- *Never resharpen the blade. This may negatively affect the tool and the performance.*
- *Ensure that the saw blade is properly inserted and securely locked before activating the tool.*
- *Do not exceed the manufacturer's recommended operating limits for the tool.*
- *Monitor the temperature of the saw blade periodically. Using irrigation while operating blades will help reduce the possibility of thermal necrosis.*
- *The blades must be able to move freely within the guide.*
- *Forceful some other side rather than oscillation may cause the fracture of the saw blade.*
- *Measurements on the blade are only for information purposes.*
- *Use care when handling sharps to prevent accidental cutting.*

Copyright Information

- *Microaire® is a registered trademark of MicroAire Surgical Instruments LLC.*
- *Stryker® is a registered trademark of Stryker Corporation.*
- *Conmed Linvatec® is a registered trademark of ConMed Corporation.*
- *De Soutter® is a registered trademark of De Soutter Medical Ltd.*
- *Smith & Nephew Dyonics® is a registered trademark of Smith & Nephew, Inc.*
- *Depuy Synthes® is a registered trademark of DePuy Synthes, Inc.*
- *Zimmer® is a registered trademark of Zimmer Inc.*
- *Brasseler® is a registered trademark of Brasseler U.S.A., Inc.*
- *TAVA® is a registered trademark of TAVA Surgical.*
- *B Braun Aesculap® is a registered trademark of B. Braun Melsungen AG and Aesculap AG.*


ACF Medical are not associated with the above listed trademark owners.




**Saw Blades
Compatible with
ACF[®]**



Brand Code 360




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
360 130 090 089	S 360 130 090 089	13 mm	90 mm	0,89 mm	0,89 mm	Straight
360 130 090 090	S 360 130 090 090	13 mm	90 mm	0,90 mm	0,90 mm	Straight
360 130 090 100	S 360 130 090 100	13 mm	90 mm	1,00 mm	1,00 mm	Straight
360 130 090 119	S 360 130 090 119	13 mm	90 mm	1,19 mm	1,19 mm	Straight
360 130 090 120	S 360 130 090 120	13 mm	90 mm	1,20 mm	1,20 mm	Straight
360 130 090 127	S 360 130 090 127	13 mm	90 mm	1,27 mm	1,27 mm	Straight
360 130 090 130	S 360 130 090 130	13 mm	90 mm	1,30 mm	1,30 mm	Straight
360 130 090 137	S 360 130 090 137	13 mm	90 mm	1,37 mm	1,37 mm	Straight
360 130 090 140	S 360 130 090 140	13 mm	90 mm	1,40 mm	1,40 mm	Straight
360 130 090 147	S 360 130 090 147	13 mm	90 mm	1,47 mm	1,47 mm	Straight
360 130 090 150	S 360 130 090 150	13 mm	90 mm	1,50 mm	1,50 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
360 190 080 090	S 360 190 080 090	19 mm	80 mm	0,90 mm	0,90 mm	Straight
360 190 080 100	S 360 190 080 100	19 mm	80 mm	1,00 mm	1,00 mm	Straight
360 190 080 120	S 360 190 080 120	19 mm	80 mm	1,20 mm	1,20 mm	Straight

ACF Medical®


Compatible Connections:
Sign Modular Systems,
A00 120 Cable System,
A00 140 Battery System




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
360 190 090 089	S 360 190 090 089	19 mm	90 mm	0,89 mm	0,89 mm	Straight
360 190 090 090	S 360 190 090 090	19 mm	90 mm	0,90 mm	0,90 mm	Straight
360 190 090 100	S 360 190 090 100	19 mm	90 mm	1,00 mm	1,00 mm	Straight
360 190 090 119	S 360 190 090 119	19 mm	90 mm	1,19 mm	1,19 mm	Straight
360 190 090 120	S 360 190 090 120	19 mm	90 mm	1,20 mm	1,20 mm	Straight
360 190 090 127	S 360 190 090 127	19 mm	90 mm	1,27 mm	1,27 mm	Straight
360 190 090 130	S 360 190 090 130	19 mm	90 mm	1,30 mm	1,30 mm	Straight
360 190 090 137	S 360 190 090 137	19 mm	90 mm	1,37 mm	1,37 mm	Straight
360 190 090 140	S 360 190 090 140	19 mm	90 mm	1,40 mm	1,40 mm	Straight
360 190 090 147	S 360 190 090 147	19 mm	90 mm	1,47 mm	1,47 mm	Straight
360 190 090 150	S 360 190 090 150	19 mm	90 mm	1,50 mm	1,50 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
360 195 090 089	S 360 195 090 089	19,5 mm	90 mm	0,89 mm	0,89 mm	Straight
360 195 090 090	S 360 195 090 090	19,5 mm	90 mm	0,90 mm	0,90 mm	Straight
360 195 090 100	S 360 195 090 100	19,5 mm	90 mm	1,00 mm	1,00 mm	Straight
360 195 090 119	S 360 195 090 119	19,5 mm	90 mm	1,19 mm	1,19 mm	Straight
360 195 090 120	S 360 195 090 120	19,5 mm	90 mm	1,20 mm	1,20 mm	Straight
360 195 090 127	S 360 195 090 127	19,5 mm	90 mm	1,27 mm	1,27 mm	Straight
360 195 090 130	S 360 195 090 130	19,5 mm	90 mm	1,30 mm	1,30 mm	Straight
360 195 090 137	S 360 195 090 137	19,5 mm	90 mm	1,37 mm	1,37 mm	Straight
360 195 090 140	S 360 195 090 140	19,5 mm	90 mm	1,40 mm	1,40 mm	Straight
360 195 090 147	S 360 195 090 147	19,5 mm	90 mm	1,47 mm	1,47 mm	Straight
360 195 090 150	S 360 195 090 150	19,5 mm	90 mm	1,50 mm	1,50 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
360 220 090 089	S 360 220 090 089	22 mm	90 mm	0,89 mm	0,89 mm	Straight
360 220 090 090	S 360 220 090 090	22 mm	90 mm	0,90 mm	0,90 mm	Straight
360 220 090 100	S 360 220 090 100	22 mm	90 mm	1,00 mm	1,00 mm	Straight
360 220 090 119	S 360 220 090 119	22 mm	90 mm	1,19 mm	1,19 mm	Straight
360 220 090 120	S 360 220 090 120	22 mm	90 mm	1,20 mm	1,20 mm	Straight
360 220 090 127	S 360 220 090 127	22 mm	90 mm	1,27 mm	1,27 mm	Straight
360 220 090 130	S 360 220 090 130	22 mm	90 mm	1,30 mm	1,30 mm	Straight
360 220 090 137	S 360 220 090 137	22 mm	90 mm	1,37 mm	1,37 mm	Straight
360 220 090 140	S 360 220 090 140	22 mm	90 mm	1,40 mm	1,40 mm	Straight
360 220 090 147	S 360 220 090 147	22 mm	90 mm	1,47 mm	1,47 mm	Straight
360 220 090 150	S 360 220 090 150	22 mm	90 mm	1,50 mm	1,50 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
360 250 090 089	S 360 250 090 089	25 mm	90 mm	0,89 mm	0,89 mm	Straight
360 250 090 090	S 360 250 090 090	25 mm	90 mm	0,90 mm	0,90 mm	Straight
360 250 090 100	S 360 250 090 100	25 mm	90 mm	1,00 mm	1,00 mm	Straight
360 250 090 119	S 360 250 090 119	25 mm	90 mm	1,19 mm	1,19 mm	Straight
360 250 090 120	S 360 250 090 120	25 mm	90 mm	1,20 mm	1,20 mm	Straight
360 250 090 127	S 360 250 090 127	25 mm	90 mm	1,27 mm	1,27 mm	Straight
360 250 090 130	S 360 250 090 130	25 mm	90 mm	1,30 mm	1,30 mm	Straight
360 250 090 137	S 360 250 090 137	25 mm	90 mm	1,37 mm	1,37 mm	Straight
360 250 090 140	S 360 250 090 140	25 mm	90 mm	1,40 mm	1,40 mm	Straight
360 250 090 147	S 360 250 090 147	25 mm	90 mm	1,47 mm	1,47 mm	Straight
360 250 090 150	S 360 250 090 150	25 mm	90 mm	1,50 mm	1,50 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
360 254 090 089	S 360 254 090 089	25,4 mm	90 mm	0,89 mm	0,89 mm	Straight
360 254 090 090	S 360 254 090 090	25,4 mm	90 mm	0,90 mm	0,90 mm	Straight
360 254 090 100	S 360 254 090 100	25,4 mm	90 mm	1,00 mm	1,00 mm	Straight
360 254 090 119	S 360 254 090 119	25,4 mm	90 mm	1,19 mm	1,19 mm	Straight
360 254 090 120	S 360 254 090 120	25,4 mm	90 mm	1,20 mm	1,20 mm	Straight
360 254 090 127	S 360 254 090 127	25,4 mm	90 mm	1,27 mm	1,27 mm	Straight
360 254 090 130	S 360 254 090 130	25,4 mm	90 mm	1,30 mm	1,30 mm	Straight
360 254 090 137	S 360 254 090 137	25,4 mm	90 mm	1,37 mm	1,37 mm	Straight
360 254 090 140	S 360 254 090 140	25,4 mm	90 mm	1,40 mm	1,40 mm	Straight
360 254 090 147	S 360 254 090 147	25,4 mm	90 mm	1,47 mm	1,47 mm	Straight
360 254 090 150	S 360 254 090 150	25,4 mm	90 mm	1,50 mm	1,50 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
360 130 064 090	S 360 130 064 090	13 mm	64 mm	0,90 mm	0,90 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
360 190 064 090	S 360 190 064 090	19 mm	64 mm	0,90 mm	0,90 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
360 250 064 090	S 360 250 064 090	25 mm	64 mm	0,90 mm	0,90 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
360 320 030 090	S 360 320 030 090	32 mm	30 mm	0,90 mm	0,90 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
360 320 064 090	S 360 320 064 090	32 mm	64 mm	0,90 mm	0,90 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
360 330 064 090	S 360 330 064 090	33 mm	64 mm	0,90 mm	0,90 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
360 055 025 060	S 360 055 025 060	5,5 mm	25 mm	0,60 mm	0,60 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
360 070 025 060	S 360 070 025 060	7 mm	25 mm	0,60 mm	0,60 mm	Straight
X 360 070 025 060	SX 360 070 025 060	7 mm	25 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
360 090 025 060	S 360 090 025 060	9 mm	25 mm	0,60	0,60 mm	Straight
X 360 090 025 060	SX 360 090 025 060	9 mm	25 mm	0,60	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
360 090 026 060	S 360 090 026 060	9 mm	26 mm	0,60 mm	0,60 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
360 090 040 060	S 360 090 040 060	9 mm	40 mm	0,60 mm	0,60 mm	Straight
X 360 090 040 060	SX 360 090 040 060	9 mm	40 mm	0,60 mm	0,80 mm	Cross



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
360 100 090 060	S 360 100 090 060	10 mm	90 mm	0,60 mm	0,60 mm	Straight
360 100 090 120	S 360 100 090 120	10 mm	90 mm	1,20 mm	1,20 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
360 195 040 060	S 360 195 040 060	19,5 mm	40 mm	0,60 mm	0,60 mm	Straight
X 360 195 040 060	SX 360 195 040 060	19,5 mm	40 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
360 150 048 060	S 360 150 048 060	15 mm	48 mm	0,60 mm	0,60 mm	Straight
X 360 150 048 060	SX 360 150 048 060	15 mm	48 mm	0,60 mm	0,80 mm	Cross

ACF Medical®


Compatible Connections;
Sign Modular Systems,
A00 120 Cable System,
A00 140 Battery System




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
360 280 048 060	S 360 280 048 060	28 mm	48 mm	0,60 mm	0,60 mm	Straight
X 360 280 048 060	SX 360 280 048 060	28 mm	48 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
360 520 030 060	S 360 520 030 060	52 mm	30 mm	0,60 mm	0,60 mm	Straight
360 520 030 090	S 360 520 030 090	52 mm	30 mm	0,90 mm	0,90 mm	Straight
X 360 520 030 060	SX 360 520 030 060	52 mm	30 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
360 520 040 060	S 360 520 040 060	52 mm	40 mm	0,60 mm	0,60 mm	Straight
360 520 040 090	S 360 520 040 090	52 mm	40 mm	0,90 mm	0,90 mm	Straight
X 360 520 040 060	SX 360 520 040 060	52 mm	40 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
360 090 050 060	S 360 090 050 060	9 mm	50 mm	0,60 mm	0,60 mm	Straight
360 090 050 100	S 360 090 050 100	9 mm	50 mm	1,00 mm	1,00 mm	Straight
X 360 090 050 060	SX 360 090 050 060	9 mm	50 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
360 090 060 060	S 360 090 060 060	9 mm	60 mm	0,60 mm	0,60 mm	Straight
360 090 060 100	S 360 090 060 100	9 mm	60 mm	1,00 mm	1,00 mm	Straight
X 360 090 060 060	SX 360 090 060 060	9 mm	60 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
360 100 070 060	S 360 100 070 060	10 mm	70 mm	0,60 mm	0,60 mm	Straight
360 100 070 100	S 360 100 070 100	10 mm	70 mm	1,00 mm	1,00 mm	Straight



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
360 100 080 060	S 360 100 080 060	10 mm	80 mm	0,60 mm	0,60 mm	Straight
360 100 080 100	S 360 100 080 100	10 mm	80 mm	1,00 mm	1,00 mm	Straight



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
360 120 090 060	S 360 120 090 060	12 mm	90 mm	0,60 mm	0,60 mm	Straight
360 120 090 100	S 360 120 090 100	12 mm	90 mm	1,00 mm	1,00 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
360 125 076 060	S 360 125 076 060	12,5 mm	76 mm	0,60 mm	0,60 mm	Straight
360 125 076 100	S 360 125 076 100	12,5 mm	76 mm	1,00 mm	1,00 mm	Straight




Saw Blades Compatible with MICROAIRE[®]




Brand Code **100**




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
100 110 090 089	S 100 110 090 089	11 mm	90 mm	0,89 mm	0,89 mm	Straight
100 110 090 119	S 100 110 090 119	11 mm	90 mm	1,19 mm	1,19 mm	Straight
100 110 090 127	S 100 110 090 127	11 mm	90 mm	1,27 mm	1,27 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
100 130 090 089	S 100 130 090 089	13 mm	90 mm	0,89 mm	0,89 mm	Straight
100 130 090 090	S 100 130 090 090	13 mm	90 mm	0,90 mm	0,90 mm	Straight
100 130 090 100	S 100 130 090 100	13 mm	90 mm	1,00 mm	1,00 mm	Straight
100 130 090 119	S 100 130 090 119	13 mm	90 mm	1,19 mm	1,19 mm	Straight
100 130 090 120	S 100 130 090 120	13 mm	90 mm	1,20 mm	1,20 mm	Straight
100 130 090 127	S 100 130 090 127	13 mm	90 mm	1,27 mm	1,27 mm	Straight
100 130 090 130	S 100 130 090 130	13 mm	90 mm	1,30 mm	1,30 mm	Straight
100 130 090 137	S 100 130 090 137	13 mm	90 mm	1,37 mm	1,37 mm	Straight
100 130 090 140	S 100 130 090 140	13 mm	90 mm	1,40 mm	1,40 mm	Straight
100 130 090 147	S 100 130 090 147	13 mm	90 mm	1,47 mm	1,47 mm	Straight
100 130 090 150	S 100 130 090 150	13 mm	90 mm	1,50 mm	1,50 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
100 190 090 089	S 100 190 090 089	19 mm	90 mm	0,89 mm	0,89 mm	Straight
100 190 090 090	S 100 190 090 090	19 mm	90 mm	0,90 mm	0,90 mm	Straight
100 190 090 100	S 100 190 090 100	19 mm	90 mm	1,00 mm	1,00 mm	Straight
100 190 090 119	S 100 190 090 119	19 mm	90 mm	1,19 mm	1,19 mm	Straight
100 190 090 120	S 100 190 090 120	19 mm	90 mm	1,20 mm	1,20 mm	Straight
100 190 090 127	S 100 190 090 127	19 mm	90 mm	1,27 mm	1,27 mm	Straight
100 190 090 130	S 100 190 090 130	19 mm	90 mm	1,30 mm	1,30 mm	Straight
100 190 090 137	S 100 190 090 137	19 mm	90 mm	1,37 mm	1,37 mm	Straight
100 190 090 140	S 100 190 090 140	19 mm	90 mm	1,40 mm	1,40 mm	Straight
100 190 090 147	S 100 190 090 147	19 mm	90 mm	1,47 mm	1,47 mm	Straight
100 190 090 150	S 100 190 090 150	19 mm	90 mm	1,50 mm	1,50 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
100 190 100 089	S 100 190 100 089	19 mm	100 mm	0,89 mm	0,89 mm	Straight
100 190 100 090	S 100 190 100 090	19 mm	100 mm	0,90 mm	0,90 mm	Straight
100 190 100 100	S 100 190 100 100	19 mm	100 mm	1,00 mm	1,00 mm	Straight
100 190 100 119	S 100 190 100 119	19 mm	100 mm	1,19 mm	1,19 mm	Straight
100 190 100 120	S 100 190 100 120	19 mm	100 mm	1,20 mm	1,20 mm	Straight
100 190 100 127	S 100 190 100 127	19 mm	100 mm	1,27 mm	1,27 mm	Straight
100 190 100 130	S 100 190 100 130	19 mm	100 mm	1,30 mm	1,30 mm	Straight
100 190 100 137	S 100 190 100 137	19 mm	100 mm	1,37 mm	1,37 mm	Straight
100 190 100 140	S 100 190 100 140	19 mm	100 mm	1,40 mm	1,40 mm	Straight
100 190 100 147	S 100 190 100 147	19 mm	100 mm	1,47 mm	1,47 mm	Straight
100 190 100 150	S 100 190 100 150	19 mm	100 mm	1,50 mm	1,50 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
100 195 090 089	S 100 195 090 089	19,5 mm	90 mm	0,89 mm	0,89 mm	Straight
100 195 090 090	S 100 195 090 090	19,5 mm	90 mm	0,90 mm	0,90 mm	Straight
100 195 090 100	S 100 195 090 100	19,5 mm	90 mm	1,00 mm	1,00 mm	Straight
100 195 090 119	S 100 195 090 119	19,5 mm	90 mm	1,19 mm	1,19 mm	Straight
100 195 090 120	S 100 195 090 120	19,5 mm	90 mm	1,20 mm	1,20 mm	Straight
100 195 090 127	S 100 195 090 127	19,5 mm	90 mm	1,27 mm	1,27 mm	Straight
100 195 090 130	S 100 195 090 130	19,5 mm	90 mm	1,30 mm	1,30 mm	Straight
100 195 090 137	S 100 195 090 137	19,5 mm	90 mm	1,37 mm	1,37 mm	Straight
100 195 090 140	S 100 195 090 140	19,5 mm	90 mm	1,40 mm	1,40 mm	Straight
100 195 090 147	S 100 195 090 147	19,5 mm	90 mm	1,47 mm	1,47 mm	Straight
100 195 090 150	S 100 195 090 150	19,5 mm	90 mm	1,50 mm	1,50 mm	Straight



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
100 220 090 089	S 100 220 090 089	22 mm	90 mm	0,89 mm	0,89 mm	Straight
100 220 090 090	S 100 220 090 090	22 mm	90 mm	0,90 mm	0,90 mm	Straight
100 220 090 100	S 100 220 090 100	22 mm	90 mm	1,00 mm	1,00 mm	Straight
100 220 090 119	S 100 220 090 119	22 mm	90 mm	1,19 mm	1,19 mm	Straight
100 220 090 120	S 100 220 090 120	22 mm	90 mm	1,20 mm	1,20 mm	Straight
100 220 090 127	S 100 220 090 127	22 mm	90 mm	1,27 mm	1,27 mm	Straight
100 220 090 130	S 100 220 090 130	22 mm	90 mm	1,30 mm	1,30 mm	Straight
100 220 090 137	S 100 220 090 137	22 mm	90 mm	1,37 mm	1,37 mm	Straight
100 220 090 140	S 100 220 090 140	22 mm	90 mm	1,40 mm	1,40 mm	Straight
100 220 090 147	S 100 220 090 147	22 mm	90 mm	1,47 mm	1,47 mm	Straight
100 220 090 150	S 100 220 090 150	22 mm	90 mm	1,50 mm	1,50 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
100 250 090 089	S 100 250 090 089	25 mm	90 mm	0,89 mm	0,89 mm	Straight
100 250 090 090	S 100 250 090 090	25 mm	90 mm	0,90 mm	0,90 mm	Straight
100 250 090 100	S 100 250 090 100	25 mm	90 mm	1,00 mm	1,00 mm	Straight
100 250 090 119	S 100 250 090 119	25 mm	90 mm	1,19 mm	1,19 mm	Straight
100 250 090 120	S 100 250 090 120	25 mm	90 mm	1,20 mm	1,20 mm	Straight
100 250 090 127	S 100 250 090 127	25 mm	90 mm	1,27 mm	1,27 mm	Straight
100 250 090 130	S 100 250 090 130	25 mm	90 mm	1,30 mm	1,30 mm	Straight
100 250 090 137	S 100 250 090 137	25 mm	90 mm	1,37 mm	1,37 mm	Straight
100 250 090 140	S 100 250 090 140	25 mm	90 mm	1,40 mm	1,40 mm	Straight
100 250 090 147	S 100 250 090 147	25 mm	90 mm	1,47 mm	1,47 mm	Straight
100 250 090 150	S 100 250 090 150	25 mm	90 mm	1,50 mm	1,50 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
100 254 100 089	S 100 254 100 089	25,4 mm	100 mm	0,89 mm	0,89 mm	Straight
100 254 100 090	S 100 254 100 090	25,4 mm	100 mm	0,90 mm	0,90 mm	Straight
100 254 100 100	S 100 254 100 100	25,4 mm	100 mm	1,00 mm	1,00 mm	Straight
100 254 100 119	S 100 254 100 119	25,4 mm	100 mm	1,19 mm	1,19 mm	Straight
100 254 100 120	S 100 254 100 120	25,4 mm	100 mm	1,20 mm	1,20 mm	Straight
100 254 100 127	S 100 254 100 127	25,4 mm	100 mm	1,27 mm	1,27 mm	Straight
100 254 100 130	S 100 254 100 130	25,4 mm	100 mm	1,30 mm	1,30 mm	Straight
100 254 100 137	S 100 254 100 137	25,4 mm	100 mm	1,37 mm	1,37 mm	Straight
100 254 100 140	S 100 254 100 140	25,4 mm	100 mm	1,40 mm	1,40 mm	Straight
100 254 100 147	S 100 254 100 147	25,4 mm	100 mm	1,47 mm	1,47 mm	Straight
100 254 100 150	S 100 254 100 150	25,4 mm	100 mm	1,50 mm	1,50 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
100 300 050 090	S 100 300 050 090	30 mm	50 mm	0,90 mm	0,90 mm	Straight
100 300 050 100	S 100 300 050 100	30 mm	50 mm	1,00 mm	1,00 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
100 480 064 090	S 100 480 064 090	48 mm	64 mm	0,90 mm	0,90 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
100 130 065 090	S 100 130 065 090	13 mm	65 mm	0,90 mm	0,90 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
100 195 062 090	S 100 195 062 090	19,5 mm	62 mm	0,90 mm	0,90 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
100 195 070 090	S 100 195 070 090	19,5 mm	70 mm	0,90 mm	0,90 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
100 320 067 090	S 100 320 067 090	32 mm	67 mm	0,90 mm	0,90 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
100 320 064 090	S 100 320 064 090	32 mm	64 mm	0,90 mm	0,90 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
100 420 050 090	S 100 420 050 090	42 mm	50 mm	0,90 mm	0,90 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
100 070 020 080	S 100 070 020 080	7 mm	20 mm	0,80 mm	0,80 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
100 090 025 060	S 100 090 025 060	9 mm	25 mm	0,60 mm	0,60 mm	Straight
X 100 090 025 060	SX 100 090 025 060	9 mm	25 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
100 090 035 060	S 100 090 035 060	9 mm	35 mm	0,60 mm	0,60 mm	Straight
X 100 090 035 060	SX 100 090 035 060	9 mm	35 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
100 150 030 060	S 100 150 030 060	15 mm	30 mm	0,60 mm	0,60 mm	Straight
X 100 150 030 060	SX 100 150 030 060	15 mm	30 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
100 195 063 060	S 100 195 063 060	19,5 mm	63 mm	0,60 mm	0,60 mm	Straight
X 100 195 063 060	SX 100 195 063 060	19,5 mm	63 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
100 090 055 060	S 100 090 055 060	9 mm	55 mm	0,60 mm	0,60 mm	Straight
100 090 055 090	S 100 090 055 090	9 mm	55 mm	0,90 mm	0,90 mm	Straight
X 100 090 055 060	SX 100 090 055 060	9 mm	55 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
100 090 080 060	S 100 090 080 060	9 mm	80 mm	0,60 mm	0,60 mm	Straight
X 100 090 080 060	SX 100 090 080 060	9 mm	80 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
100 120 090 060	S 100 120 090 060	12 mm	90 mm	0,60 mm	0,60 mm	Straight
100 120 090 090	S 100 120 090 090	12 mm	90 mm	0,90 mm	0,90 mm	Straight
X 100 120 090 060	SX 100 120 090 060	12 mm	90 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
100 120 070 060	S 100 120 070 060	12 mm	70 mm	0,60 mm	0,60 mm	Straight
100 120 070 090	S 100 120 070 090	12 mm	70 mm	0,90 mm	0,90 mm	Straight
100 120 070 100	S 100 120 070 100	12 mm	70 mm	1,00 mm	1,00 mm	Straight
100 120 070 120	S 100 120 070 120	12 mm	70 mm	1,20 mm	1,20 mm	Straight
X 100 120 070 060	SX 100 120 070 060	12 mm	70 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
100 120 076 060	S 100 120 076 060	12 mm	76 mm	0,60 mm	0,60 mm	Straight
100 120 076 090	S 100 120 076 090	12 mm	76 mm	0,90 mm	0,90 mm	Straight
100 120 076 100	S 100 120 076 100	12 mm	76 mm	1,00 mm	1,00 mm	Straight
100 120 076 120	S 100 120 076 120	12 mm	76 mm	1,20 mm	1,20 mm	Straight
X 100 120 076 060	SX 100 120 076 060	12 mm	76 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
100 090 036 060	S 100 090 036 060	9 mm	36 mm	0,60 mm	0,60 mm	Straight
X 100 090 036 060	SX 100 090 036 060	9 mm	36 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
100 040 006 060	S 100 040 006 060	4 mm	6 mm	0,60 mm	0,60 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
100 045 026 060	S 100 045 026 060	4,5 mm	26 mm	0,60 mm	0,60 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
100 055 018 060	S 100 055 018 060	5,5 mm	18 mm	0,60 mm	0,60 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
100 055 026 060	S 100 055 026 060	5,5 mm	26 mm	0,60 mm	0,60 mm	Straight



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
100 060 033 060	S 100 060 033 060	6 mm	33 mm	0,60 mm	0,60 mm	Straight



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
100 095 025 060	S 100 095 025 060	9,5 mm	25 mm	0,60 mm	0,60 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
100 100 018 060	S 100 100 018 060	10 mm	18 mm	0,60 mm	0,60 mm	Straight
X 100 100 018 060	SX 100 100 018 060	10 mm	18 mm	0,60 mm	0,80 mm	Cross



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
100 100 030 060	S 100 100 030 060	10 mm	30 mm	0,60 mm	0,60 mm	Straight
X 100 100 030 060	SX 100 100 030 060	10 mm	30 mm	0,60 mm	0,80 mm	Cross

Microaire®


Compatible Connections:
1950, 1955, 2250, 6670,
6672, 6673 sagittal saws




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
100 100 032 060	S 100 100 032 060	10 mm	32 mm	0,60 mm	0,60 mm	Straight
X 100 100 032 060	SX 100 100 032 060	10 mm	32 mm	0,60 mm	0,80 mm	Cross



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
100 100 033 060	S 100 100 033 060	10 mm	33 mm	0,60 mm	0,60 mm	Straight
X 100 100 033 060	SX 100 100 033 060	10 mm	33 mm	0,60 mm	0,80 mm	Cross



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
100 100 040 060	S 100 100 040 060	10 mm	40mm	0,60 mm	0,60 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
100 150 040 060	S 100 150 040 060	15 mm	40 mm	0,60 mm	0,60 mm	Straight
X 100 150 040 060	SX 100 150 040 060	15 mm	40 mm	0,60 mm	0,80 mm	Cross




Saw Blades Compatible with STRYKER®



Brand Code **120**




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
120 130 085 089	S 120 130 085 089	13 mm	85 mm	0,89 mm	0,89 mm	Straight
120 130 085 090	S 120 130 085 090	13 mm	85 mm	0,90 mm	0,90 mm	Straight
120 130 085 100	S 120 130 085 100	13 mm	85 mm	1,00 mm	1,00 mm	Straight
120 130 085 119	S 120 130 085 119	13 mm	85 mm	1,19 mm	1,19 mm	Straight
120 130 085 120	S 120 130 085 120	13 mm	85 mm	1,20 mm	1,20 mm	Straight
120 130 085 127	S 120 130 085 127	13 mm	85 mm	1,27 mm	1,27 mm	Straight
120 130 085 130	S 120 130 085 130	13 mm	85 mm	1,30 mm	1,30 mm	Straight
120 130 085 137	S 120 130 085 137	13 mm	85 mm	1,37 mm	1,37 mm	Straight
120 130 085 140	S 120 130 085 140	13 mm	85 mm	1,40 mm	1,40 mm	Straight
120 130 085 147	S 120 130 085 147	13 mm	85 mm	1,47 mm	1,47 mm	Straight
120 130 085 150	S 120 130 085 150	13 mm	85 mm	1,50 mm	1,50 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
120 130 090 089	S 120 130 090 089	13 mm	90 mm	0,89 mm	0,89 mm	Straight
120 130 090 090	S 120 130 090 090	13 mm	90 mm	0,90 mm	0,90 mm	Straight
120 130 090 100	S 120 130 090 100	13 mm	90 mm	1,00 mm	1,00 mm	Straight
120 130 090 119	S 120 130 090 119	13 mm	90 mm	1,19 mm	1,19 mm	Straight
120 130 090 120	S 120 130 090 120	13 mm	90 mm	1,20 mm	1,20 mm	Straight
120 130 090 127	S 120 130 090 127	13 mm	90 mm	1,27 mm	1,27 mm	Straight
120 130 090 130	S 120 130 090 130	13 mm	90 mm	1,30 mm	1,30 mm	Straight
120 130 090 137	S 120 130 090 137	13 mm	90 mm	1,37 mm	1,37 mm	Straight
120 130 090 140	S 120 130 090 140	13 mm	90 mm	1,40 mm	1,40 mm	Straight
120 130 090 147	S 120 130 090 147	13 mm	90 mm	1,47 mm	1,47 mm	Straight
120 130 090 150	S 120 130 090 150	13 mm	90 mm	1,50 mm	1,50 mm	Straight

Stryker®


Compatible Connections:
Stryker System EHD(296-98),
2000(2108), 54(4108), 55(4208),
56(6208), 57(7208), 58(8208)




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
120 190 085 089	S 120 190 085 089	19 mm	85 mm	0,89 mm	0,89 mm	Straight
120 190 085 090	S 120 190 085 090	19 mm	85 mm	0,90 mm	0,90 mm	Straight
120 190 085 100	S 120 190 085 100	19 mm	85 mm	1,00 mm	1,00 mm	Straight
120 190 085 119	S 120 190 085 119	19 mm	85 mm	1,19 mm	1,19 mm	Straight
120 190 085 120	S 120 190 085 120	19 mm	85 mm	1,20 mm	1,20 mm	Straight
120 190 085 127	S 120 190 085 127	19 mm	85 mm	1,27 mm	1,27 mm	Straight
120 190 085 130	S 120 190 085 130	19 mm	85 mm	1,30 mm	1,30 mm	Straight
120 190 085 137	S 120 190 085 137	19 mm	85 mm	1,37 mm	1,37 mm	Straight
120 190 085 140	S 120 190 085 140	19 mm	85 mm	1,40 mm	1,40 mm	Straight
120 190 085 147	S 120 190 085 147	19 mm	85 mm	1,47 mm	1,47 mm	Straight
120 190 085 150	S 120 190 085 150	19 mm	85 mm	1,50 mm	1,50 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
120 190 090 089	S 120 190 090 089	19 mm	90 mm	0,89 mm	0,89 mm	Straight
120 190 090 090	S 120 190 090 090	19 mm	90 mm	0,90 mm	0,90 mm	Straight
120 190 090 100	S 120 190 090 100	19 mm	90 mm	1,00 mm	1,00 mm	Straight
120 190 090 119	S 120 190 090 119	19 mm	90 mm	1,19 mm	1,19 mm	Straight
120 190 090 120	S 120 190 090 120	19 mm	90 mm	1,20 mm	1,20 mm	Straight
120 190 090 127	S 120 190 090 127	19 mm	90 mm	1,27 mm	1,27 mm	Straight
120 190 090 130	S 120 190 090 130	19 mm	90 mm	1,30 mm	1,30 mm	Straight
120 190 090 137	S 120 190 090 137	19 mm	90 mm	1,37 mm	1,37 mm	Straight
120 190 090 140	S 120 190 090 140	19 mm	90 mm	1,40 mm	1,40 mm	Straight
120 190 090 147	S 120 190 090 147	19 mm	90 mm	1,47 mm	1,47 mm	Straight
120 190 090 150	S 120 190 090 150	19 mm	90 mm	1,50 mm	1,50 mm	Straight



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
120 190 070 120	S 120 190 070 120	19 mm	70 mm	1,20 mm	1,20 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
120 190 095 100	S 120 190 095 100	19 mm	95 mm	1,00 mm	1,00 mm	Straight
120 190 095 120	S 120 190 095 120	19 mm	95 mm	1,20 mm	1,20 mm	Straight
120 190 095 127	S 120 190 095 127	19 mm	95 mm	1,27 mm	1,27 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
120 190 100 089	S 120 190 100 089	19 mm	100 mm	0,89 mm	0,89 mm	Straight
120 190 100 090	S 120 190 100 090	19 mm	100 mm	0,90 mm	0,90 mm	Straight
120 190 100 100	S 120 190 100 100	19 mm	100 mm	1,00 mm	1,00 mm	Straight
120 190 100 119	S 120 190 100 119	19 mm	100 mm	1,19 mm	1,19 mm	Straight
120 190 100 120	S 120 190 100 120	19 mm	100 mm	1,20 mm	1,20 mm	Straight
120 190 100 127	S 120 190 100 127	19 mm	100 mm	1,27 mm	1,27 mm	Straight
120 190 100 130	S 120 190 100 130	19 mm	100 mm	1,30 mm	1,30 mm	Straight
120 190 100 137	S 120 190 100 137	19 mm	100 mm	1,37 mm	1,37 mm	Straight
120 190 100 140	S 120 190 100 140	19 mm	100 mm	1,40 mm	1,40 mm	Straight
120 190 100 147	S 120 190 100 147	19 mm	100 mm	1,47 mm	1,47 mm	Straight
120 190 100 150	S 120 190 100 150	19 mm	100 mm	1,50 mm	1,50 mm	Straight

Stryker®
Compatible Connections:
Stryker System EHD(296-98),
2000(2108), S4(4108), S5(4208),
S6(6208), S7(7208), S8(8208)




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
120 195 090 089	S 120 195 090 089	19,5 mm	90 mm	0,89 mm	0,89 mm	Straight
120 195 090 090	S 120 195 090 090	19,5 mm	90 mm	0,90 mm	0,90 mm	Straight
120 195 090 100	S 120 195 090 100	19,5 mm	90 mm	1,00 mm	1,00 mm	Straight
120 195 090 119	S 120 195 090 119	19,5 mm	90 mm	1,19 mm	1,19 mm	Straight
120 195 090 120	S 120 195 090 120	19,5 mm	90 mm	1,20 mm	1,20 mm	Straight
120 195 090 127	S 120 195 090 127	19,5 mm	90 mm	1,27 mm	1,27 mm	Straight
120 195 090 130	S 120 195 090 130	19,5 mm	90 mm	1,30 mm	1,30 mm	Straight
120 195 090 137	S 120 195 090 137	19,5 mm	90 mm	1,37 mm	1,37 mm	Straight
120 195 090 140	S 120 195 090 140	19,5 mm	90 mm	1,40 mm	1,40 mm	Straight
120 195 090 147	S 120 195 090 147	19,5 mm	90 mm	1,47 mm	1,47 mm	Straight
120 195 090 150	S 120 195 090 150	19,5 mm	90 mm	1,50 mm	1,50 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
120 220 090 089	S 120 220 090 089	22 mm	90 mm	0,89 mm	0,89 mm	Straight
120 220 090 090	S 120 220 090 090	22 mm	90 mm	0,90 mm	0,90 mm	Straight
120 220 090 100	S 120 220 090 100	22 mm	90 mm	1,00 mm	1,00 mm	Straight
120 220 090 119	S 120 220 090 119	22 mm	90 mm	1,19 mm	1,19 mm	Straight
120 220 090 120	S 120 220 090 120	22 mm	90 mm	1,20 mm	1,20 mm	Straight
120 220 090 127	S 120 220 090 127	22 mm	90 mm	1,27 mm	1,27 mm	Straight
120 220 090 130	S 120 220 090 130	22 mm	90 mm	1,30 mm	1,30 mm	Straight
120 220 090 137	S 120 220 090 137	22 mm	90 mm	1,37 mm	1,37 mm	Straight
120 220 090 140	S 120 220 090 140	22 mm	90 mm	1,40 mm	1,40 mm	Straight
120 220 090 147	S 120 220 090 147	22 mm	90 mm	1,47 mm	1,47 mm	Straight
120 220 090 150	S 120 220 090 150	22 mm	90 mm	1,50 mm	1,50 mm	Straight

STRYKER®




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
120 250 085 089	S 120 250 085 089	25 mm	85 mm	0,89 mm	0,89 mm	Straight
120 250 085 090	S 120 250 085 090	25 mm	85 mm	0,90 mm	0,90 mm	Straight
120 250 085 100	S 120 250 085 100	25 mm	85 mm	1,00 mm	1,00 mm	Straight
120 250 085 119	S 120 250 085 119	25 mm	85 mm	1,19 mm	1,19 mm	Straight
120 250 085 120	S 120 250 085 120	25 mm	85 mm	1,20 mm	1,20 mm	Straight
120 250 085 127	S 120 250 085 127	25 mm	85 mm	1,27 mm	1,27 mm	Straight
120 250 085 130	S 120 250 085 130	25 mm	85 mm	1,30 mm	1,30 mm	Straight
120 250 085 137	S 120 250 085 137	25 mm	85 mm	1,37 mm	1,37 mm	Straight
120 250 085 140	S 120 250 085 140	25 mm	85 mm	1,40 mm	1,40 mm	Straight
120 250 085 147	S 120 250 085 147	25 mm	85 mm	1,47 mm	1,47 mm	Straight
120 250 085 150	S 120 250 085 150	25 mm	85 mm	1,50 mm	1,50 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
120 250 090 089	S 120 250 090 089	25 mm	90 mm	0,89 mm	0,89 mm	Straight
120 250 090 090	S 120 250 090 090	25 mm	90 mm	0,90 mm	0,90 mm	Straight
120 250 090 100	S 120 250 090 100	25 mm	90 mm	1,00 mm	1,00 mm	Straight
120 250 090 119	S 120 250 090 119	25 mm	90 mm	1,19 mm	1,19 mm	Straight
120 250 090 120	S 120 250 090 120	25 mm	90 mm	1,20 mm	1,20 mm	Straight
120 250 090 127	S 120 250 090 127	25 mm	90 mm	1,27 mm	1,27 mm	Straight
120 250 090 130	S 120 250 090 130	25 mm	90 mm	1,30 mm	1,30 mm	Straight
120 250 090 137	S 120 250 090 137	25 mm	90 mm	1,37 mm	1,37 mm	Straight
120 250 090 140	S 120 250 090 140	25 mm	90 mm	1,40 mm	1,40 mm	Straight
120 250 090 147	S 120 250 090 147	25 mm	90 mm	1,47 mm	1,47 mm	Straight
120 250 090 150	S 120 250 090 150	25 mm	90 mm	1,50 mm	1,50 mm	Straight


Stryker®
Compatible Connections:
Stryker System EHD(296-98),
2000(2108), 54(4108), 55(4208),
56(6208), 57(7208), 58(8208)




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
120 254 090 089	S 120 254 090 089	25,4 mm	90 mm	0,89 mm	0,89 mm	Straight
120 254 090 090	S 120 254 090 090	25,4 mm	90 mm	0,90 mm	0,90 mm	Straight
120 254 090 100	S 120 254 090 100	25,4 mm	90 mm	1,00 mm	1,00 mm	Straight
120 254 090 119	S 120 254 090 119	25,4 mm	90 mm	1,19 mm	1,19 mm	Straight
120 254 090 120	S 120 254 090 120	25,4 mm	90 mm	1,20 mm	1,20 mm	Straight
120 254 090 127	S 120 254 090 127	25,4 mm	90 mm	1,27 mm	1,27 mm	Straight
120 254 090 130	S 120 254 090 130	25,4 mm	90 mm	1,30 mm	1,30 mm	Straight
120 254 090 137	S 120 254 090 137	25,4 mm	90 mm	1,37 mm	1,37 mm	Straight
120 254 090 140	S 120 254 090 140	25,4 mm	90 mm	1,40 mm	1,40 mm	Straight
120 254 090 147	S 120 254 090 147	25,4 mm	90 mm	1,47 mm	1,47 mm	Straight
120 254 090 150	S 120 254 090 150	25,4 mm	90 mm	1,50 mm	1,50 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
120 250 075 089	S 120 250 075 089	25 mm	75 mm	0,89 mm	0,89 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
120 250 095 100	S 120 250 095 100	25 mm	95 mm	1,00 mm	1,00 mm	Straight
120 250 095 120	S 120 250 095 120	25 mm	95 mm	1,20 mm	1,20 mm	Straight
120 250 095 127	S 120 250 095 127	25 mm	95 mm	1,27 mm	1,27 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
120 070 020 080	S 120 070 020 080	7 mm	20 mm	0,80 mm	0,80 mm	Straight



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
120 090 020 060	S 120 090 020 060	9 mm	20 mm	0,60 mm	0,60 mm	Straight
120 090 020 090	S 120 090 020 090	9 mm	20 mm	0,90 mm	0,90 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
120 090 025 060	S 120 090 025 060	9 mm	25 mm	0,60 mm	0,60 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
120 090 035 060	S 120 090 035 060	9 mm	35 mm	0,60 mm	0,60 mm	Straight
X 120 090 035 060	SX 120 090 035 060	9 mm	35 mm	0,60 mm	0,80 mm	Cross

Stryker®


Compatible Connections:
 Stryker System EHD(296-98),
 2000(2108), S4(4108), S5(4208),
 S6(6208), S7(7208), S8(8208)



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
120 150 035 060	S 120 150 035 060	15 mm	35 mm	0,60 mm	0,60 mm	Straight
X 120 150 035 060	SX 120 150 035 060	15 mm	35 mm	0,60 mm	0,80 mm	Cross
120 150 035 100	S 120 150 035 100	15 mm	35 mm	1,00 mm	1,00 mm	Straight
120 150 035 120	S 120 150 035 120	15 mm	35 mm	1,20 mm	1,20 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
120 150 030 120	S 120 150 030 120	15 mm	30 mm	1,20 mm	1,20 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
120 150 045 100	S 120 150 045 100	15 mm	45 mm	1,00 mm	1,00 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
120 270 058 060	S 120 270 058 060	27 mm	58 mm	0,60 mm	0,60 mm	Straight
X 120 270 058 060	SX 120 270 058 060	27 mm	58 mm	0,60 mm	0,80 mm	Cross



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
120 300 032 060	S 120 300 032 060	30 mm	32 mm	0,60 mm	0,60 mm	Straight
X 120 300 032 060	SX 120 300 032 060	30 mm	32 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
120 405 475 060	S 120 405 475 060	40,5 mm	47,5 mm	0,60 mm	0,60 mm	Straight
X 120 405 475 060	SX 120 405 475 060	40,5 mm	47,5 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
120 530 032 060	S 120 530 032 060	53 mm	32 mm	0,60 mm	0,60 mm	Straight
X 120 530 032 060	SX 120 530 032 060	53 mm	32 mm	0,60 mm	0,80 mm	Cross


Stryker®
Compatible Connections:
Stryker System EHD(296-98),
2000(2108), S4(4108), S5(4208),
S6(6208), S7(7208), S8(8208)




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
120 600 030 120	S 120 600 030 120	60 mm	30 mm	1,20 mm	1,20 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
120 800 025 060	S 120 800 025 060	80 mm	25 mm	0,60 mm	0,60 mm	Straight
X 120 800 025 060	SX 120 800 025 060	80 mm	25 mm	0,60 mm	0,80 mm	Cross



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
120 060 032 060	S 120 060 032 060	6 mm	32 mm	0,60 mm	0,60 mm	Straight
X 120 060 032 060	SX 120 060 032 060	6 mm	32 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
120 060 060 060	S 120 060 060 060	6 mm	60 mm	0,60 mm	0,60 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
120 105 072 060	S 120 105 072 060	10,5 mm	72 mm	0,60 mm	0,60 mm	Straight
X 120 105 072 060	SX 120 105 072 060	10,5 mm	72 mm	0,60 mm	0,80 mm	Cross
120 105 072 090	S 120 105 072 090	10,5 mm	72 mm	0,90 mm	0,90 mm	Straight
120 105 072 100	S 120 105 072 100	10,5 mm	72 mm	1,00 mm	1,00 mm	Straight
120 105 072 120	S 120 105 072 120	10,5 mm	72 mm	1,20 mm	1,20 mm	Straight


Stryker®
Compatible Connections;
4206, 4106, 2106, 296-96, 6206,
7206, 8206 Reciprocating Saws



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
120 105 077 060	S 120 105 077 060	10,5 mm	77 mm	0,60 mm	0,60 mm	Straight
X 120 105 077 060	SX 120 105 077 060	10,5 mm	77 mm	0,60 mm	0,80 mm	Cross
120 105 077 090	S 120 105 077 090	10,5 mm	77 mm	0,90 mm	0,90 mm	Straight
120 105 077 100	S 120 105 077 100	10,5 mm	77 mm	1,00 mm	1,00 mm	Straight
120 105 077 120	S 120 105 077 120	10,5 mm	77 mm	1,20 mm	1,20 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
120 125 068 090	S 120 125 068 090	12,5 mm	68 mm	0,90 mm	0,90 mm	Straight
120 125 068 100	S 120 125 068 100	12,5 mm	68 mm	1,00 mm	1,00 mm	Straight
120 125 068 120	S 120 125 068 120	12,5 mm	68 mm	1,20 mm	1,20 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
120 135 068 090	S 120 135 068 090	13,5 mm	68 mm	0,90 mm	0,90 mm	Straight
120 135 068 100	S 120 135 068 100	13,5 mm	68 mm	1,00 mm	1,00 mm	Straight
120 135 068 120	S 120 135 068 120	13,5 mm	68 mm	1,20 mm	1,20 mm	Straight

TPS5100-34 Sagittal Saw, TPS5100-31 Oscillating Saw, Command2 2296-31, 2296-234, 2296-34, 2296-231, 2296-080-134, Cordless Driver 2+3+4 Att. 4100-400, Sabo Sagittal Saw 4300-034, Sabo 2 4408-000, Micro Bone Saw, Core 5400-34, 5400-31, NSK Primado, RemB Micro Power 6400-034, 6400-031, 4100-400




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
120 025 018 060	S 120 025 018 060	2,5 mm	18 mm	0,60 mm	0,60 mm	Straight



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
120 050 012 060	S 120 050 012 060	5 mm	12 mm	0,60 mm	0,60 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
120 050 018 060	S 120 050 018 060	5 mm	18 mm	0,60 mm	0,60 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
120 050 025 060	S 120 050 025 060	5 mm	25 mm	0,60 mm	0,60 mm	Straight

Compatible Connections;


TP55100-34 Sagittal Saw, TP55100-31 Oscillating Saw, Command2 2296-31, 2296-234, 2296-34, 2296-231, 2296-080-134, Cordless Driver 2+3+4 Att. 4100-400, Sabo Sagittal Saw 4300-034, Sabo 2 4408-000, Micro Bone Saw, Core 5400-34, 5400-31, NSK Primado, RemB Micro Power 6400-034, 6400-031, 4100-400




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
120 070 015 040	S 120 070 015 040	7 mm	15 mm	0,40 mm	0,40 mm	Straight



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
120 070 018 060	S 120 070 018 060	7 mm	18 mm	0,60 mm	0,60 mm	Straight
X 120 070 018 060	SX 120 070 018 060	7 mm	18 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
120 070 025 060	S 120 070 025 060	7 mm	25 mm	0,60 mm	0,60 mm	Straight
X 120 070 025 060	SX 120 070 025 060	7 mm	25 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
X 120 080 022 030	SX 120 080 022 030	8 mm	22 mm	0,30 mm	0,40 mm	Cross


TPS5100-34 Sagittal Saw, TPS5100-31 Oscillating Saw, Command2 2296-31, 2296-234, 2296-34, 2296-231, 2296-080-134, Cordless Driver 2+3+4 Att. 4100-400, Sabo Sagittal Saw 4300-034, Sabo 2 4408-000, Micro Bone Saw, Core 5400-34, 5400-31, NSK Primado, RemB Micro Power 6400-034, 6400-031, 4100-400




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
120 080 025 060	S 120 080 025 060	8 mm	25 mm	0,60 mm	0,60 mm	Straight
X 120 080 025 060	SX 120 080 025 060	8 mm	25 mm	0,60 mm	0,80 mm	Cross



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
120 080 035 060	S 120 080 035 060	8 mm	35 mm	0,60 mm	0,60 mm	Straight
X 120 080 035 060	SX 120 080 035 060	8 mm	35 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
120 090 012 060	S 120 090 012 060	9 mm	12 mm	0,60 mm	0,60 mm	Straight
X 120 090 012 060	SX 120 090 012 060	9 mm	12 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
120 090 018 060	S 120 090 018 060	9 mm	18 mm	0,60 mm	0,60 mm	Straight
X 120 090 018 060	SX 120 090 018 060	9 mm	18 mm	0,60 mm	0,80 mm	Cross

Compatible Connections;


TP55100-34 Sagittal Saw, TP55100-31 Oscillating Saw, Command2 2296-31, 2296-234, 2296-34, 2296-231, 2296-080-134, Cordless Driver 2+3+4 Att. 4100-400, Sabo Sagittal Saw 4300-034, Sabo 2 4408-000, Micro Bone Saw, Core 5400-34, 5400-31, NSK Primado, RemB Micro Power 6400-034, 6400-031, 4100-400




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
120 090 024 040	S 120 090 024 040	9 mm	24 mm	0,40 mm	0,40 mm	Straight
X 120 090 024 040	SX 120 090 024 040	9 mm	24 mm	0,40 mm	0,60 mm	Cross
120 090 024 060	S 120 090 024 060	9 mm	24 mm	0,60 mm	0,60 mm	Straight
X 120 090 024 060	SX 120 090 024 060	9 mm	24 mm	0,60 mm	0,80 mm	Cross



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
X 120 090 030 040	SX 120 090 030 040	9 mm	30 mm	0,40 mm	0,60 mm	Cross
120 090 030 060	S 120 090 030 060	9 mm	30 mm	0,60 mm	0,60 mm	Straight
X 120 090 030 060	SX 120 090 030 060	9 mm	30 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
120 100 025 060	S 120 100 025 060	10 mm	25 mm	0,60 mm	0,60 mm	Straight
X 120 100 025 060	SX 120 100 025 060	10 mm	25 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
120 100 031 030	S 120 100 031 030	10 mm	31 mm	0,30 mm	0,30 mm	Straight
X 120 100 031 030	SX 120 100 031 030	10 mm	31 mm	0,30 mm	0,40 mm	Cross


TPS5100-34 Sagittal Saw, TPS5100-31 Oscillating Saw, Command2 2296-31, 2296-234, 2296-34, 2296-231, 2296-080-134, Cordless Driver 2+3+4 Att. 4100-400, Sabo Sagittal Saw 4300-034, Sabo 2 4408-000, Micro Bone Saw, Core 5400-34, 5400-31, NSK Primado, RemB Micro Power 6400-034, 6400-031, 4100-400



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
120 100 035 060	S 120 100 035 060	10 mm	35 mm	0,60 mm	0,60 mm	Straight
X 120 100 035 060	SX 120 100 035 060	10 mm	35 mm	0,60 mm	0,80 mm	Cross



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
120 100 050 040	S 120 100 050 040	10 mm	50mm	0,40 mm	0,40 mm	Straight
X 120 100 050 040	SX 120 100 050 040	10 mm	50mm	0,40 mm	0,60 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
120 120 031 060	S 120 120 031 060	12 mm	31 mm	0,60 mm	0,60 mm	Straight
X 120 120 031 060	SX 120 120 031 060	12 mm	31 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
120 130 035 040	S 120 130 035 040	13 mm	35 mm	0,40 mm	0,40 mm	Straight
120 130 035 060	S 120 130 035 060	13 mm	35 mm	0,60 mm	0,60 mm	Straight
X 120 130 035 060	SX 120 130 035 060	13 mm	35 mm	0,60 mm	0,80 mm	Cross

Compatible Connections;


TP55100-34 Sagittal Saw, TP55100-31 Oscillating Saw, Command2 2296-31, 2296-234, 2296-34, 2296-231, 2296-080-134, Cordless Driver 2+3+4 Att. 4100-400, Sabo Sagittal Saw 4300-034, Sabo 2 4408-000, Micro Bone Saw, Core 5400-34, 5400-31, NSK Primado, RemB Micro Power 6400-034, 6400-031, 4100-400



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
120 130 040 060	S 120 130 040 060	13 mm	40 mm	0,60 mm	0,60 mm	Straight
X 120 130 040 060	SX 120 130 040 060	13 mm	40 mm	0,60 mm	0,80 mm	Cross



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
120 140 031 060	S 120 140 031 060	14 mm	31 mm	0,60 mm	0,60 mm	Straight
X 120 140 031 060	SX 120 140 031 060	14 mm	31 mm	0,60 mm	0,80 mm	Cross



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
X 120 200 050 040	SX 120 200 050 040	20 mm	50 mm	0,40 mm	0,60 mm	Cross




Saw Blades Compatible with CONMED LINVATEC®




Brand Code **140**


For Reciprocating Saws and Some Micro Blades
100




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
140 130 090 089	S 140 130 090 089	13 mm	90 mm	0,89 mm	0,89 mm	Straight
140 130 090 090	S 140 130 090 090	13 mm	90 mm	0,90 mm	0,90 mm	Straight
140 130 090 100	S 140 130 090 100	13 mm	90 mm	1,00 mm	1,00 mm	Straight
140 130 090 119	S 140 130 090 119	13 mm	90 mm	1,19 mm	1,19 mm	Straight
140 130 090 120	S 140 130 090 120	13 mm	90 mm	1,20 mm	1,20 mm	Straight
140 130 090 127	S 140 130 090 127	13 mm	90 mm	1,27 mm	1,27 mm	Straight
140 130 090 130	S 140 130 090 130	13 mm	90 mm	1,30 mm	1,30 mm	Straight
140 130 090 137	S 140 130 090 137	13 mm	90 mm	1,37 mm	1,37 mm	Straight
140 130 090 140	S 140 130 090 140	13 mm	90 mm	1,40 mm	1,40 mm	Straight
140 130 090 147	S 140 130 090 147	13 mm	90 mm	1,47 mm	1,47 mm	Straight
140 130 090 150	S 140 130 090 150	13 mm	90 mm	1,50 mm	1,50 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
140 190 090 089	S 140 190 090 089	19 mm	90 mm	0,89 mm	0,89 mm	Straight
140 190 090 090	S 140 190 090 090	19 mm	90 mm	0,90 mm	0,90 mm	Straight
140 190 090 100	S 140 190 090 100	19 mm	90 mm	1,00 mm	1,00 mm	Straight
140 190 090 119	S 140 190 090 119	19 mm	90 mm	1,19 mm	1,19 mm	Straight
140 190 090 120	S 140 190 090 120	19 mm	90 mm	1,20 mm	1,20 mm	Straight
140 190 090 127	S 140 190 090 127	19 mm	90 mm	1,27 mm	1,27 mm	Straight
140 190 090 130	S 140 190 090 130	19 mm	90 mm	1,30 mm	1,30 mm	Straight
140 190 090 137	S 140 190 090 137	19 mm	90 mm	1,37 mm	1,37 mm	Straight
140 190 090 140	S 140 190 090 140	19 mm	90 mm	1,40 mm	1,40 mm	Straight
140 190 090 147	S 140 190 090 147	19 mm	90 mm	1,47 mm	1,47 mm	Straight
140 190 090 150	S 140 190 090 150	19 mm	90 mm	1,50 mm	1,50 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
140 190 100 089	S 140 190 100 089	19 mm	100 mm	0,89 mm	0,89 mm	Straight
140 190 100 090	S 140 190 100 090	19 mm	100 mm	0,90 mm	0,90 mm	Straight
140 190 100 100	S 140 190 100 100	19 mm	100 mm	1,00 mm	1,00 mm	Straight
140 190 100 119	S 140 190 100 119	19 mm	100 mm	1,19 mm	1,19 mm	Straight
140 190 100 120	S 140 190 100 120	19 mm	100 mm	1,20 mm	1,20 mm	Straight
140 190 100 127	S 140 190 100 127	19 mm	100 mm	1,27 mm	1,27 mm	Straight
140 190 100 130	S 140 190 100 130	19 mm	100 mm	1,30 mm	1,30 mm	Straight
140 190 100 137	S 140 190 100 137	19 mm	100 mm	1,37 mm	1,37 mm	Straight
140 190 100 140	S 140 190 100 140	19 mm	100 mm	1,40 mm	1,40 mm	Straight
140 190 100 147	S 140 190 100 147	19 mm	100 mm	1,47 mm	1,47 mm	Straight
140 190 100 150	S 140 190 100 150	19 mm	100 mm	1,50 mm	1,50 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
140 195 090 089	S 140 195 090 089	19,5 mm	90 mm	0,89 mm	0,89 mm	Straight
140 195 090 090	S 140 195 090 090	19,5 mm	90 mm	0,90 mm	0,90 mm	Straight
140 195 090 100	S 140 195 090 100	19,5 mm	90 mm	1,00 mm	1,00 mm	Straight
140 195 090 119	S 140 195 090 119	19,5 mm	90 mm	1,19 mm	1,19 mm	Straight
140 195 090 120	S 140 195 090 120	19,5 mm	90 mm	1,20 mm	1,20 mm	Straight
140 195 090 127	S 140 195 090 127	19,5 mm	90 mm	1,27 mm	1,27 mm	Straight
140 195 090 130	S 140 195 090 130	19,5 mm	90 mm	1,30 mm	1,30 mm	Straight
140 195 090 137	S 140 195 090 137	19,5 mm	90 mm	1,37 mm	1,37 mm	Straight
140 195 090 140	S 140 195 090 140	19,5 mm	90 mm	1,40 mm	1,40 mm	Straight
140 195 090 147	S 140 195 090 147	19,5 mm	90 mm	1,47 mm	1,47 mm	Straight
140 195 090 150	S 140 195 090 150	19,5 mm	90 mm	1,50 mm	1,50 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
140 220 090 089	S 140 220 090 089	22 mm	90 mm	0,89 mm	0,89 mm	Straight
140 220 090 090	S 140 220 090 090	22 mm	90 mm	0,90 mm	0,90 mm	Straight
140 220 090 100	S 140 220 090 100	22 mm	90 mm	1,00 mm	1,00 mm	Straight
140 220 090 119	S 140 220 090 119	22 mm	90 mm	1,19 mm	1,19 mm	Straight
140 220 090 120	S 140 220 090 120	22 mm	90 mm	1,20 mm	1,20 mm	Straight
140 220 090 127	S 140 220 090 127	22 mm	90 mm	1,27 mm	1,27 mm	Straight
140 220 090 130	S 140 220 090 130	22 mm	90 mm	1,30 mm	1,30 mm	Straight
140 220 090 137	S 140 220 090 137	22 mm	90 mm	1,37 mm	1,37 mm	Straight
140 220 090 140	S 140 220 090 140	22 mm	90 mm	1,40 mm	1,40 mm	Straight
140 220 090 147	S 140 220 090 147	22 mm	90 mm	1,47 mm	1,47 mm	Straight
140 220 090 150	S 140 220 090 150	22 mm	90 mm	1,50 mm	1,50 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
140 250 090 089	S 140 250 090 089	25 mm	90 mm	0,89 mm	0,89 mm	Straight
140 250 090 090	S 140 250 090 090	25 mm	90 mm	0,90 mm	0,90 mm	Straight
140 250 090 100	S 140 250 090 100	25 mm	90 mm	1,00 mm	1,00 mm	Straight
140 250 090 119	S 140 250 090 119	25 mm	90 mm	1,19 mm	1,19 mm	Straight
140 250 090 120	S 140 250 090 120	25 mm	90 mm	1,20 mm	1,20 mm	Straight
140 250 090 127	S 140 250 090 127	25 mm	90 mm	1,27 mm	1,27 mm	Straight
140 250 090 130	S 140 250 090 130	25 mm	90 mm	1,30 mm	1,30 mm	Straight
140 250 090 137	S 140 250 090 137	25 mm	90 mm	1,37 mm	1,37 mm	Straight
140 250 090 140	S 140 250 090 140	25 mm	90 mm	1,40 mm	1,40 mm	Straight
140 250 090 147	S 140 250 090 147	25 mm	90 mm	1,47 mm	1,47 mm	Straight
140 250 090 150	S 140 250 090 150	25 mm	90 mm	1,50 mm	1,50 mm	Straight



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
140 254 090 089	S 140 254 090 089	25,4 mm	90 mm	0,89 mm	0,89 mm	Straight
140 254 090 090	S 140 254 090 090	25,4 mm	90 mm	0,90 mm	0,90 mm	Straight
140 254 090 100	S 140 254 090 100	25,4 mm	90 mm	1,00 mm	1,00 mm	Straight
140 254 090 119	S 140 254 090 119	25,4 mm	90 mm	1,19 mm	1,19 mm	Straight
140 254 090 120	S 140 254 090 120	25,4 mm	90 mm	1,20 mm	1,20 mm	Straight
140 254 090 127	S 140 254 090 127	25,4 mm	90 mm	1,27 mm	1,27 mm	Straight
140 254 090 130	S 140 254 090 130	25,4 mm	90 mm	1,30 mm	1,30 mm	Straight
140 254 090 137	S 140 254 090 137	25,4 mm	90 mm	1,37 mm	1,37 mm	Straight
140 254 090 140	S 140 254 090 140	25,4 mm	90 mm	1,40 mm	1,40 mm	Straight
140 254 090 147	S 140 254 090 147	25,4 mm	90 mm	1,47 mm	1,47 mm	Straight
140 254 090 150	S 140 254 090 150	25,4 mm	90 mm	1,50 mm	1,50 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
140 055 025 060	S 140 055 025 060	5,5 mm	25 mm	0,60 mm	0,60 mm	Straight
140 055 025 080	S 140 055 025 080	5,5 mm	25 mm	0,80 mm	0,80 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
140 070 020 080	S 140 070 020 080	7 mm	20 mm	0,80 mm	0,80 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
140 070 025 060	S 140 070 025 060	7 mm	25 mm	0,60 mm	0,60 mm	Straight
X 140 070 025 060	SX 140 070 025 060	7 mm	25 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
140 090 025 060	S 140 090 025 060	9 mm	25 mm	0,60 mm	0,60 mm	Straight
X 140 090 025 060	SX 140 090 025 060	9 mm	25 mm	0,60 mm	0,80 mm	Cross
140 090 025 090	S 140 090 025 090	9 mm	25 mm	0,90 mm	0,90 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
140 090 040 060	S 140 090 040 060	9 mm	40 mm	0,60 mm	0,60 mm	Straight
X 140 090 040 060	SX 140 090 040 060	9 mm	40 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
140 095 025 060	S 140 095 025 060	9,5 mm	25 mm	0,60 mm	0,60 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
140 150 048 060	S 140 150 048 060	15 mm	48 mm	0,60 mm	0,60 mm	Straight
X 140 150 048 060	SX 140 150 048 060	15 mm	48 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
140 195 040 060	S 140 195 040 060	19,5 mm	40 mm	0,60 mm	0,60 mm	Straight
X 140 195 040 060	SX 140 195 040 060	19,5 mm	40 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
140 280 048 060	S 140 280 048 060	28 mm	48 mm	0,60 mm	0,60 mm	Straight
X 140 280 048 060	SX 140 280 048 060	28 mm	48 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
140 130 064 090	S 140 130 064 090	13 mm	64 mm	0,90 mm	0,90 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
140 190 064 090	S 140 190 064 090	19 mm	64 mm	0,90 mm	0,90 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
140 250 064 090	S 140 250 064 090	25 mm	64 mm	0,90 mm	0,90 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
140 320 064 090	S 140 320 064 090	32 mm	64 mm	0,90 mm	0,90 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
140 320 030 090	S 140 320 030 090	32 mm	30 mm	0,90 mm	0,90 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
140 330 064 090	S 140 330 064 090	33 mm	64 mm	0,90 mm	0,90 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
140 520 030 090	S 140 520 030 090	52 mm	30 mm	0,90 mm	0,90 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
140 520 040 090	S 140 520 040 090	52 mm	40 mm	0,90 mm	0,90 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
100 090 055 060	S 100 090 055 060	9 mm	55 mm	0,60 mm	0,60 mm	Straight
100 090 055 090	S 100 090 055 090	9 mm	55 mm	0,90 mm	0,90 mm	Straight
X 100 090 055 060	SX 100 090 055 060	9 mm	55 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
100 090 080 060	S 100 090 080 060	9 mm	80 mm	0,60 mm	0,60 mm	Straight
X 100 090 080 060	SX 100 090 080 060	9 mm	80 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
100 120 090 060	S 100 120 090 060	12 mm	90 mm	0,60 mm	0,60 mm	Straight
100 120 090 090	S 100 120 090 090	12 mm	90 mm	0,90 mm	0,90 mm	Straight
X 100 120 090 060	SX 100 120 090 060	12 mm	90 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
100 120 070 060	S 100 120 070 060	12 mm	70 mm	0,60 mm	0,60 mm	Straight
100 120 070 090	S 100 120 070 090	12 mm	70 mm	0,90 mm	0,90 mm	Straight
100 120 070 100	S 100 120 070 100	12 mm	70 mm	1,00 mm	1,00 mm	Straight
100 120 070 120	S 100 120 070 120	12 mm	70 mm	1,20 mm	1,20 mm	Straight
X 100 120 070 060	SX 100 120 070 060	12 mm	70 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
100 120 076 060	S 100 120 076 060	12 mm	76 mm	0,60 mm	0,60 mm	Straight
100 120 076 090	S 100 120 076 090	12 mm	76 mm	0,90 mm	0,90 mm	Straight
100 120 076 100	S 100 120 076 100	12 mm	76 mm	1,00 mm	1,00 mm	Straight
100 120 076 120	S 100 120 076 120	12 mm	76 mm	1,20 mm	1,20 mm	Straight
X 100 120 076 060	SX 100 120 076 060	12 mm	76 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
100 076 025 060	S 100 076 025 060	7,6 mm	25 mm	0,60 mm	0,60 mm	Straight
X 100 076 025 060	SX 100 076 025 060	7,6 mm	25 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
100 076 060 060	S 100 076 060 060	7,6 mm	60 mm	0,60 mm	0,60 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
100 040 006 060	S 100 040 006 060	4 mm	6 mm	0,60 mm	0,60 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
100 045 026 060	S 100 045 026 060	4,5 mm	26 mm	0,60 mm	0,60 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
100 055 018 060	S 100 055 018 060	5,5 mm	18 mm	0,60 mm	0,60 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
100 055 026 060	S 100 055 026 060	5,5 mm	26 mm	0,60 mm	0,60 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
100 060 033 060	S 100 060 033 060	6 mm	33 mm	0,60 mm	0,60 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
100 095 025 060	S 100 095 025 060	9,5 mm	25 mm	0,60 mm	0,60 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
100 100 018 060	S 100 100 018 060	10 mm	18 mm	0,60 mm	0,60 mm	Straight
X 100 100 018 060	SX 100 100 018 060	10 mm	18 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
100 100 030 060	S 100 100 030 060	10 mm	30 mm	0,60 mm	0,60 mm	Straight
X 100 100 030 060	SX 100 100 030 060	10 mm	30 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
100 100 032 060	S 100 100 032 060	10 mm	32 mm	0,60 mm	0,60 mm	Straight
X 100 100 032 060	SX 100 100 032 060	10 mm	32 mm	0,60 mm	0,80 mm	Cross



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
100 100 033 060	S 100 100 033 060	10 mm	33 mm	0,60 mm	0,60 mm	Straight
X 100 100 033 060	SX 100 100 033 060	10 mm	33 mm	0,60 mm	0,80 mm	Cross



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
100 100 040 060	S 100 100 040 060	10 mm	40mm	0,60 mm	0,60 mm	Straight



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
100 150 040 060	S 100 150 040 060	15 mm	40 mm	0,60 mm	0,60 mm	Straight
X 100 150 040 060	SX 100 150 040 060	15 mm	40 mm	0,60 mm	0,80 mm	Cross




**Saw Blades
Compatible with
DE SOUTTER®**




Brand Code


200 / 205



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 130 085 127	S 200 130 085 127	13 mm	85 mm	1,27 mm	1,27 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 130 090 089	S 200 130 090 089	13 mm	90 mm	0,89 mm	0,89 mm	Straight
200 130 090 090	S 200 130 090 090	13 mm	90 mm	0,90 mm	0,90 mm	Straight
200 130 090 100	S 200 130 090 100	13 mm	90 mm	1,00 mm	1,00 mm	Straight
200 130 090 119	S 200 130 090 119	13 mm	90 mm	1,19 mm	1,19 mm	Straight
200 130 090 120	S 200 130 090 120	13 mm	90 mm	1,20 mm	1,20 mm	Straight
200 130 090 127	S 200 130 090 127	13 mm	90 mm	1,27 mm	1,27 mm	Straight
200 130 090 130	S 200 130 090 130	13 mm	90 mm	1,30 mm	1,30 mm	Straight
200 130 090 137	S 200 130 090 137	13 mm	90 mm	1,37 mm	1,37 mm	Straight
200 130 090 140	S 200 130 090 140	13 mm	90 mm	1,40 mm	1,40 mm	Straight
200 130 090 147	S 200 130 090 147	13 mm	90 mm	1,47 mm	1,47 mm	Straight
200 130 090 150	S 200 130 090 150	13 mm	90 mm	1,50 mm	1,50 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 190 085 090	S 200 190 085 090	19 mm	85 mm	0,90 mm	0,90 mm	Straight

De Soutter®


Compatible Connections;
MDX Multi Drive KX 600, MDX Multi Drive KDX-600,
MPX Multi Drive SX 500, MPX Multi Drive DPX-676,
MCX Multi Drive DCX 500, Orthodrive DBK-700




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 190 090 089	S 200 190 090 089	19 mm	90 mm	0,89 mm	0,89 mm	Straight
200 190 090 090	S 200 190 090 090	19 mm	90 mm	0,90 mm	0,90 mm	Straight
200 190 090 100	S 200 190 090 100	19 mm	90 mm	1,00 mm	1,00 mm	Straight
200 190 090 119	S 200 190 090 119	19 mm	90 mm	1,19 mm	1,19 mm	Straight
200 190 090 120	S 200 190 090 120	19 mm	90 mm	1,20 mm	1,20 mm	Straight
200 190 090 127	S 200 190 090 127	19 mm	90 mm	1,27 mm	1,27 mm	Straight
200 190 090 130	S 200 190 090 130	19 mm	90 mm	1,30 mm	1,30 mm	Straight
200 190 090 137	S 200 190 090 137	19 mm	90 mm	1,37 mm	1,37 mm	Straight
200 190 090 140	S 200 190 090 140	19 mm	90 mm	1,40 mm	1,40 mm	Straight
200 190 090 147	S 200 190 090 147	19 mm	90 mm	1,47 mm	1,47 mm	Straight
200 190 090 150	S 200 190 090 150	19 mm	90 mm	1,50 mm	1,50 mm	Straight



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 190 100 089	S 200 190 100 089	19 mm	100 mm	0,89 mm	0,89 mm	Straight
200 190 100 090	S 200 190 100 090	19 mm	100 mm	0,90 mm	0,90 mm	Straight
200 190 100 100	S 200 190 100 100	19 mm	100 mm	1,00 mm	1,00 mm	Straight
200 190 100 119	S 200 190 100 119	19 mm	100 mm	1,19 mm	1,19 mm	Straight
200 190 100 120	S 200 190 100 120	19 mm	100 mm	1,20 mm	1,20 mm	Straight
200 190 100 127	S 200 190 100 127	19 mm	100 mm	1,27 mm	1,27 mm	Straight
200 190 100 130	S 200 190 100 130	19 mm	100 mm	1,30 mm	1,30 mm	Straight
200 190 100 137	S 200 190 100 137	19 mm	100 mm	1,37 mm	1,37 mm	Straight
200 190 100 140	S 200 190 100 140	19 mm	100 mm	1,40 mm	1,40 mm	Straight
200 190 100 147	S 200 190 100 147	19 mm	100 mm	1,47 mm	1,47 mm	Straight
200 190 100 150	S 200 190 100 150	19 mm	100 mm	1,50 mm	1,50 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 195 090 089	S 200 195 090 089	19,5 mm	90 mm	0,89 mm	0,89 mm	Straight
200 195 090 090	S 200 195 090 090	19,5 mm	90 mm	0,90 mm	0,90 mm	Straight
200 195 090 100	S 200 195 090 100	19,5 mm	90 mm	1,00 mm	1,00 mm	Straight
200 195 090 119	S 200 195 090 119	19,5 mm	90 mm	1,19 mm	1,19 mm	Straight
200 195 090 120	S 200 195 090 120	19,5 mm	90 mm	1,20 mm	1,20 mm	Straight
200 195 090 127	S 200 195 090 127	19,5 mm	90 mm	1,27 mm	1,27 mm	Straight
200 195 090 130	S 200 195 090 130	19,5 mm	90 mm	1,30 mm	1,30 mm	Straight
200 195 090 137	S 200 195 090 137	19,5 mm	90 mm	1,37 mm	1,37 mm	Straight
200 195 090 140	S 200 195 090 140	19,5 mm	90 mm	1,40 mm	1,40 mm	Straight
200 195 090 147	S 200 195 090 147	19,5 mm	90 mm	1,47 mm	1,47 mm	Straight
200 195 090 150	S 200 195 090 150	19,5 mm	90 mm	1,50 mm	1,50 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 220 090 089	S 200 220 090 089	22 mm	90 mm	0,89 mm	0,89 mm	Straight
200 220 090 090	S 200 220 090 090	22 mm	90 mm	0,90 mm	0,90 mm	Straight
200 220 090 100	S 200 220 090 100	22 mm	90 mm	1,00 mm	1,00 mm	Straight
200 220 090 119	S 200 220 090 119	22 mm	90 mm	1,19 mm	1,19 mm	Straight
200 220 090 120	S 200 220 090 120	22 mm	90 mm	1,20 mm	1,20 mm	Straight
200 220 090 127	S 200 220 090 127	22 mm	90 mm	1,27 mm	1,27 mm	Straight
200 220 090 130	S 200 220 090 130	22 mm	90 mm	1,30 mm	1,30 mm	Straight
200 220 090 137	S 200 220 090 137	22 mm	90 mm	1,37 mm	1,37 mm	Straight
200 220 090 140	S 200 220 090 140	22 mm	90 mm	1,40 mm	1,40 mm	Straight
200 220 090 147	S 200 220 090 147	22 mm	90 mm	1,47 mm	1,47 mm	Straight
200 220 090 150	S 200 220 090 150	22 mm	90 mm	1,50 mm	1,50 mm	Straight

Compatible Connections;


MDX Multi Drive KX 600, MDX Multi Drive KDX-600,
MPX Multi Drive SX 500, MPX Multi Drive DPX-676,
MCX Multi Drive DCX 500, Orthodrive DBK-700




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 250 090 089	S 200 250 090 089	25 mm	90 mm	0,89 mm	0,89 mm	Straight
200 250 090 090	S 200 250 090 090	25 mm	90 mm	0,90 mm	0,90 mm	Straight
200 250 090 100	S 200 250 090 100	25 mm	90 mm	1,00 mm	1,00 mm	Straight
200 250 090 119	S 200 250 090 119	25 mm	90 mm	1,19 mm	1,19 mm	Straight
200 250 090 120	S 200 250 090 120	25 mm	90 mm	1,20 mm	1,20 mm	Straight
200 250 090 127	S 200 250 090 127	25 mm	90 mm	1,27 mm	1,27 mm	Straight
200 250 090 130	S 200 250 090 130	25 mm	90 mm	1,30 mm	1,30 mm	Straight
200 250 090 137	S 200 250 090 137	25 mm	90 mm	1,37 mm	1,37 mm	Straight
200 250 090 140	S 200 250 090 140	25 mm	90 mm	1,40 mm	1,40 mm	Straight
200 250 090 147	S 200 250 090 147	25 mm	90 mm	1,47 mm	1,47 mm	Straight
200 250 090 150	S 200 250 090 150	25 mm	90 mm	1,50 mm	1,50 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 254 090 089	S 200 254 090 089	25,4 mm	90 mm	0,89 mm	0,89 mm	Straight
200 254 090 090	S 200 254 090 090	25,4 mm	90 mm	0,90 mm	0,90 mm	Straight
200 254 090 100	S 200 254 090 100	25,4 mm	90 mm	1,00 mm	1,00 mm	Straight
200 254 090 119	S 200 254 090 119	25,4 mm	90 mm	1,19 mm	1,19 mm	Straight
200 254 090 120	S 200 254 090 120	25,4 mm	90 mm	1,20 mm	1,20 mm	Straight
200 254 090 127	S 200 254 090 127	25,4 mm	90 mm	1,27 mm	1,27 mm	Straight
200 254 090 130	S 200 254 090 130	25,4 mm	90 mm	1,30 mm	1,30 mm	Straight
200 254 090 137	S 200 254 090 137	25,4 mm	90 mm	1,37 mm	1,37 mm	Straight
200 254 090 140	S 200 254 090 140	25,4 mm	90 mm	1,40 mm	1,40 mm	Straight
200 254 090 147	S 200 254 090 147	25,4 mm	90 mm	1,47 mm	1,47 mm	Straight
200 254 090 150	S 200 254 090 150	25,4 mm	90 mm	1,50 mm	1,50 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 250 085 127	S 200 250 085 127	25 mm	85 mm	1,27 mm	1,27 mm	Straight



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 070 020 080	S 200 070 020 080	7 mm	20 mm	0,80 mm	0,80 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 090 025 060	S 200 090 025 060	9 mm	25 mm	0,60 mm	0,60 mm	Straight
X 200 090 025 060	SX 200 090 025 060	9 mm	25 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 120 030 060	S 200 120 030 060	12 mm	30 mm	0,60 mm	0,60 mm	Straight

Compatible Connections;


MDX Multi Drive KX 600, MDX Multi Drive KDX-600,
MPX Multi Drive SX 500, MPX Multi Drive DPX-676,
MCX Multi Drive DCX 500, Orthodrive DBK-700



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 120 025 060	S 200 120 025 060	12 mm	25 mm	0,60 mm	0,60 mm	Straight
X 200 120 025 060	SX 200 120 025 060	12 mm	25 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 120 060 060	S 200 120 060 060	12 mm	60 mm	0,60 mm	0,60 mm	Straight
X 200 120 060 060	SX 200 120 060 060	12 mm	60 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 150 025 060	S 200 150 025 060	15 mm	25 mm	0,60 mm	0,60 mm	Straight
X 200 150 025 060	SX 200 150 025 060	15 mm	25 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 150 048 060	S 200 150 048 060	15 mm	48 mm	0,60 mm	0,60 mm	Straight
X 200 150 048 060	SX 200 150 048 060	15 mm	48 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 200 055 060	S 200 200 055 060	20 mm	55 mm	0,60 mm	0,60 mm	Straight
X 200 200 055 060	SX 200 200 055 060	20 mm	55 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 270 048 060	S 200 270 048 060	27 mm	48 mm	0,60 mm	0,60 mm	Straight
X 200 270 048 060	SX 200 270 048 060	27 mm	48 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 270 068 060	S 200 270 068 060	27 mm	68 mm	0,60 mm	0,60 mm	Straight
X 200 270 068 060	SX 200 270 068 060	27 mm	68 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 090 055 060	S 200 090 055 060	9 mm	55 mm	0,60 mm	0,60 mm	Straight
X 200 090 055 060	SX 200 090 055 060	9 mm	55 mm	0,60 mm	0,80 mm	Cross
200 090 055 090	S 200 090 055 090	9 mm	55 mm	0,90 mm	0,90 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 090 080 060	S 200 090 080 060	9 mm	80 mm	0,60 mm	0,60 mm	Straight
X 200 090 080 060	SX 200 090 080 060	9 mm	80 mm	0,60 mm	0,80 mm	Cross
200 090 080 090	S 200 090 080 090	9 mm	80 mm	0,90 mm	0,90 mm	Straight
200 090 080 100	S 200 090 080 100	9 mm	80 mm	1,00 mm	1,00 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 090 089 100	S 200 090 089 100	9 mm	89 mm	1,00 mm	1,00 mm	Straight
200 090 089 120	S 200 090 089 120	9 mm	89 mm	1,20 mm	1,20 mm	Straight
200 090 089 130	S 200 090 089 130	9 mm	89 mm	1,30 mm	1,30 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 100 089 120	S 200 100 089 120	10 mm	89 mm	1,20 mm	1,20 mm	Straight
200 100 089 130	S 200 100 089 130	10 mm	89 mm	1,30 mm	1,30 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 120 070 100	S 200 120 070 100	12 mm	70 mm	1,00 mm	1,00 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 120 082 120	S 200 120 082 120	12 mm	82 mm	1,20 mm	1,20 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 055 019 060	S 200 055 019 060	5,5 mm	19 mm	0,60 mm	0,60 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 055 023 060	S 200 055 023 060	5,5 mm	23 mm	0,60 mm	0,60 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 055 026 060	S 200 055 026 060	5,5 mm	26 mm	0,60 mm	0,60 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 055 034 060	S 200 055 034 060	5,5 mm	34 mm	0,60 mm	0,60 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 060 030 060	S 200 060 030 060	6 mm	30 mm	0,60 mm	0,60 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 070 016 060	S 200 070 016 060	7 mm	16 mm	0,60 mm	0,60 mm	Straight
X 200 070 016 060	SX 200 070 016 060	7 mm	16 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 070 020 060	S 200 070 020 060	7 mm	20 mm	0,60 mm	0,60 mm	Straight
X 200 070 020 060	SX 200 070 020 060	7 mm	20 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 070 022 060	S 200 070 022 060	7 mm	22 mm	0,60 mm	0,60 mm	Straight
X 200 070 022 060	SX 200 070 022 060	7 mm	22 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 070 031 060	S 200 070 031 060	7 mm	31 mm	0,60 mm	0,60 mm	Straight
X 200 070 031 060	SX 200 070 031 060	7 mm	31 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 080 030 060	S 200 080 030 060	8 mm	30 mm	0,60 mm	0,60 mm	Straight
X 200 080 030 060	SX 200 080 030 060	8 mm	30 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 090 022 060	S 200 090 022 060	9 mm	22 mm	0,60 mm	0,60 mm	Straight
X 200 090 022 060	SX 200 090 022 060	9 mm	22 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 090 027 060	S 200 090 027 060	9 mm	27 mm	0,60 mm	0,60 mm	Straight
X 200 090 027 060	SX 200 090 027 060	9 mm	27 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 090 030 060	S 200 090 030 060	9 mm	30 mm	0,60 mm	0,60 mm	Straight
X 200 090 030 060	SX 200 090 030 060	9 mm	30 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 090 032 060	S 200 090 032 060	9 mm	32 mm	0,60 mm	0,60 mm	Straight
X 200 090 032 060	SX 200 090 032 060	9 mm	32 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 090 036 060	S 200 090 036 060	9 mm	36 mm	0,60 mm	0,60 mm	Straight
X 200 090 036 060	SX 200 090 036 060	9 mm	36 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 090 039 060	S 200 090 039 060	9 mm	39 mm	0,60 mm	0,60 mm	Straight
X 200 090 039 060	SX 200 090 039 060	9 mm	39 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 130 026 060	S 200 130 026 060	13 mm	26 mm	0,60 mm	0,60 mm	Straight
X 200 130 026 060	SX 200 130 026 060	13 mm	26 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 130 039 060	S 200 130 039 060	13 mm	39 mm	0,60 mm	0,60 mm	Straight
X 200 130 039 060	SX 200 130 039 060	13 mm	39 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
200 190 039 060	S 200 190 039 060	19 mm	39 mm	0,60 mm	0,60 mm	Straight
X 200 190 039 060	SX 200 190 039 060	19 mm	39 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
205 130 090 089	S 205 130 090 089	13 mm	90 mm	0,89 mm	0,89 mm	Straight
205 130 090 090	S 205 130 090 090	13 mm	90 mm	0,90 mm	0,90 mm	Straight
205 130 090 100	S 205 130 090 100	13 mm	90 mm	1,00 mm	1,00 mm	Straight
205 130 090 119	S 205 130 090 119	13 mm	90 mm	1,19 mm	1,19 mm	Straight
205 130 090 120	S 205 130 090 120	13 mm	90 mm	1,20 mm	1,20 mm	Straight
205 130 090 127	S 205 130 090 127	13 mm	90 mm	1,27 mm	1,27 mm	Straight
205 130 090 130	S 205 130 090 130	13 mm	90 mm	1,30 mm	1,30 mm	Straight
205 130 090 137	S 205 130 090 137	13 mm	90 mm	1,37 mm	1,37 mm	Straight
205 130 090 140	S 205 130 090 140	13 mm	90 mm	1,40 mm	1,40 mm	Straight
205 130 090 147	S 205 130 090 147	13 mm	90 mm	1,47 mm	1,47 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
205 190 090 089	S 205 190 090 089	19 mm	90 mm	0,89 mm	0,89 mm	Straight
205 190 090 090	S 205 190 090 090	19 mm	90 mm	0,90 mm	0,90 mm	Straight
205 190 090 100	S 205 190 090 100	19 mm	90 mm	1,00 mm	1,00 mm	Straight
205 190 090 119	S 205 190 090 119	19 mm	90 mm	1,19 mm	1,19 mm	Straight
205 190 090 120	S 205 190 090 120	19 mm	90 mm	1,20 mm	1,20 mm	Straight
205 190 090 127	S 205 190 090 127	19 mm	90 mm	1,27 mm	1,27 mm	Straight
205 190 090 130	S 205 190 090 130	19 mm	90 mm	1,30 mm	1,30 mm	Straight
205 190 090 137	S 205 190 090 137	19 mm	90 mm	1,37 mm	1,37 mm	Straight
205 190 090 140	S 205 190 090 140	19 mm	90 mm	1,40 mm	1,40 mm	Straight
205 190 090 147	S 205 190 090 147	19 mm	90 mm	1,47 mm	1,47 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
205 190 100 089	S 205 190 100 089	19 mm	100 mm	0,89 mm	0,89 mm	Straight
205 190 100 090	S 205 190 100 090	19 mm	100 mm	0,90 mm	0,90 mm	Straight
205 190 100 100	S 205 190 100 100	19 mm	100 mm	1,00 mm	1,00 mm	Straight
205 190 100 119	S 205 190 100 119	19 mm	100 mm	1,19 mm	1,19 mm	Straight
205 190 100 120	S 205 190 100 120	19 mm	100 mm	1,20 mm	1,20 mm	Straight
205 190 100 127	S 205 190 100 127	19 mm	100 mm	1,27 mm	1,27 mm	Straight
205 190 100 130	S 205 190 100 130	19 mm	100 mm	1,30 mm	1,30 mm	Straight
205 190 100 137	S 205 190 100 137	19 mm	100 mm	1,37 mm	1,37 mm	Straight
205 190 100 140	S 205 190 100 140	19 mm	100 mm	1,40 mm	1,40 mm	Straight
205 190 100 147	S 205 190 100 147	19 mm	100 mm	1,47 mm	1,47 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
205 195 090 089	S 205 195 090 089	19.5 mm	90 mm	0,89 mm	0,89 mm	Straight
205 195 090 090	S 205 195 090 090	19.5 mm	90 mm	0,90 mm	0,90 mm	Straight
205 195 090 100	S 205 195 090 100	19.5 mm	90 mm	1,00 mm	1,00 mm	Straight
205 195 090 119	S 205 195 090 119	19.5 mm	90 mm	1,19 mm	1,19 mm	Straight
205 195 090 120	S 205 195 090 120	19.5 mm	90 mm	1,20 mm	1,20 mm	Straight
205 195 090 127	S 205 195 090 127	19.5 mm	90 mm	1,27 mm	1,27 mm	Straight
205 195 090 130	S 205 195 090 130	19.5 mm	90 mm	1,30 mm	1,30 mm	Straight
205 195 090 137	S 205 195 090 137	19.5 mm	90 mm	1,37 mm	1,37 mm	Straight
205 195 090 140	S 205 195 090 140	19.5 mm	90 mm	1,40 mm	1,40 mm	Straight
205 195 090 147	S 205 195 090 147	19.5 mm	90 mm	1,47 mm	1,47 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
205 220 090 089	S 205 220 090 089	22 mm	90 mm	0,89 mm	0,89 mm	Straight
205 220 090 090	S 205 220 090 090	22 mm	90 mm	0,90 mm	0,90 mm	Straight
205 220 090 100	S 205 220 090 100	22 mm	90 mm	1,00 mm	1,00 mm	Straight
205 220 090 119	S 205 220 090 119	22 mm	90 mm	1,19 mm	1,19 mm	Straight
205 220 090 120	S 205 220 090 120	22 mm	90 mm	1,20 mm	1,20 mm	Straight
205 220 090 127	S 205 220 090 127	22 mm	90 mm	1,27 mm	1,27 mm	Straight
205 220 090 130	S 205 220 090 130	22 mm	90 mm	1,30 mm	1,30 mm	Straight
205 220 090 137	S 205 220 090 137	22 mm	90 mm	1,37 mm	1,37 mm	Straight
205 220 090 140	S 205 220 090 140	22 mm	90 mm	1,40 mm	1,40 mm	Straight
205 220 090 147	S 205 220 090 147	22 mm	90 mm	1,47 mm	1,47 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
205 223 090 089	S 205 223 090 089	22,3 mm	90 mm	0,89 mm	0,89 mm	Straight
205 223 090 090	S 205 223 090 090	22,3 mm	90 mm	0,90 mm	0,90 mm	Straight
205 223 090 100	S 205 223 090 100	22,3 mm	90 mm	1,00 mm	1,00 mm	Straight
205 223 090 119	S 205 223 090 119	22,3 mm	90 mm	1,19 mm	1,19 mm	Straight
205 223 090 120	S 205 223 090 120	22,3 mm	90 mm	1,20 mm	1,20 mm	Straight
205 223 090 127	S 205 223 090 127	22,3 mm	90 mm	1,27 mm	1,27 mm	Straight
205 223 090 130	S 205 223 090 130	22,3 mm	90 mm	1,30 mm	1,30 mm	Straight
205 223 090 137	S 205 223 090 137	22,3 mm	90 mm	1,37 mm	1,37 mm	Straight
205 223 090 140	S 205 223 090 140	22,3 mm	90 mm	1,40 mm	1,40 mm	Straight
205 223 090 147	S 205 223 090 147	22,3 mm	90 mm	1,47 mm	1,47 mm	Straight



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
205 250 090 089	S 205 250 090 089	25 mm	90 mm	0,89 mm	0,89 mm	Straight
205 250 090 090	S 205 250 090 090	25 mm	90 mm	0,90 mm	0,90 mm	Straight
205 250 090 100	S 205 250 090 100	25 mm	90 mm	1,00 mm	1,00 mm	Straight
205 250 090 119	S 205 250 090 119	25 mm	90 mm	1,19 mm	1,19 mm	Straight
205 250 090 120	S 205 250 090 120	25 mm	90 mm	1,20 mm	1,20 mm	Straight
205 250 090 127	S 205 250 090 127	25 mm	90 mm	1,27 mm	1,27 mm	Straight
205 250 090 130	S 205 250 090 130	25 mm	90 mm	1,30 mm	1,30 mm	Straight
205 250 090 137	S 205 250 090 137	25 mm	90 mm	1,37 mm	1,37 mm	Straight
205 250 090 140	S 205 250 090 140	25 mm	90 mm	1,40 mm	1,40 mm	Straight
205 250 090 147	S 205 250 090 147	25 mm	90 mm	1,47 mm	1,47 mm	Straight



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
205 251 090 089	S 205 251 090 089	25,1 mm	90 mm	0,89 mm	0,89 mm	Straight
205 251 090 090	S 205 251 090 090	25,1 mm	90 mm	0,90 mm	0,90 mm	Straight
205 251 090 100	S 205 251 090 100	25,1 mm	90 mm	1,00 mm	1,00 mm	Straight
205 251 090 119	S 205 251 090 119	25,1 mm	90 mm	1,19 mm	1,19 mm	Straight
205 251 090 120	S 205 251 090 120	25,1 mm	90 mm	1,20 mm	1,20 mm	Straight
205 251 090 127	S 205 251 090 127	25,1 mm	90 mm	1,27 mm	1,27 mm	Straight
205 251 090 130	S 205 251 090 130	25,1 mm	90 mm	1,30 mm	1,30 mm	Straight
205 251 090 137	S 205 251 090 137	25,1 mm	90 mm	1,37 mm	1,37 mm	Straight
205 251 090 140	S 205 251 090 140	25,1 mm	90 mm	1,40 mm	1,40 mm	Straight
205 251 090 147	S 205 251 090 147	25,1 mm	90 mm	1,47 mm	1,47 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
205 254 090 089	S 205 254 090 089	25,4 mm	90 mm	0,89 mm	0,89 mm	Straight
205 254 090 090	S 205 254 090 090	25,4 mm	90 mm	0,90 mm	0,90 mm	Straight
205 254 090 100	S 205 254 090 100	25,4 mm	90 mm	1,00 mm	1,00 mm	Straight
205 254 090 119	S 205 254 090 119	25,4 mm	90 mm	1,19 mm	1,19 mm	Straight
205 254 090 120	S 205 254 090 120	25,4 mm	90 mm	1,20 mm	1,20 mm	Straight
205 254 090 127	S 205 254 090 127	25,4 mm	90 mm	1,27 mm	1,27 mm	Straight
205 254 090 130	S 205 254 090 130	25,4 mm	90 mm	1,30 mm	1,30 mm	Straight
205 254 090 137	S 205 254 090 137	25,4 mm	90 mm	1,37 mm	1,37 mm	Straight
205 254 090 140	S 205 254 090 140	25,4 mm	90 mm	1,40 mm	1,40 mm	Straight
205 254 090 147	S 205 254 090 147	25,4 mm	90 mm	1,47 mm	1,47 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
205 070 020 060	S 205 070 020 060	7 mm	20 mm	0,60 mm	0,60 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
205 090 025 060	S 205 090 025 060	9 mm	25 mm	0,60 mm	0,60 mm	Straight



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
205 120 025 060	S 205 120 025 060	12 mm	25 mm	0,60 mm	0,60 mm	Straight



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
205 150 025 060	S 205 150 025 060	15 mm	25 mm	0,60 mm	0,60 mm	Straight



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
205 130 045 089	S 205 130 045 089	13 mm	45 mm	0,89 mm	0,89 mm	Straight
205 130 045 090	S 205 130 045 090	13 mm	45 mm	0,90 mm	0,90 mm	Straight
205 130 045 100	S 205 130 045 100	13 mm	45 mm	1,00 mm	1,00 mm	Straight
205 130 045 119	S 205 130 045 119	13 mm	45 mm	1,19 mm	1,19 mm	Straight
205 130 045 120	S 205 130 045 120	13 mm	45 mm	1,20 mm	1,20 mm	Straight
205 130 045 127	S 205 130 045 127	13 mm	45 mm	1,27 mm	1,27 mm	Straight
205 130 045 130	S 205 130 045 130	13 mm	45 mm	1,30 mm	1,30 mm	Straight
205 130 045 137	S 205 130 045 137	13 mm	45 mm	1,37 mm	1,37 mm	Straight
205 130 045 140	S 205 130 045 140	13 mm	45 mm	1,40 mm	1,40 mm	Straight
205 130 045 147	S 205 130 045 147	13 mm	45 mm	1,47 mm	1,47 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
205 150 048 060	S 205 150 048 060	15 mm	48 mm	0,60 mm	0,60 mm	Straight




**Saw Blades
Compatible with
SMITH & NEPHEW DYONICS®**



Brand Code 210




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
210 045 025 060	S 210 045 025 060	4,5 mm	25 mm	0,60 mm	0,60 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
210 050 026 060	S 210 050 026 060	5 mm	26 mm	0,60 mm	0,60 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
210 058 018 060	S 210 058 018 060	5,8 mm	18 mm	0,60 mm	0,60 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
210 058 025 060	S 210 058 025 060	5,8 mm	25 mm	0,60 mm	0,60 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
210 090 025 060	S 210 090 025 060	9 mm	25 mm	0,60 mm	0,60 mm	Straight
X 210 090 025 060	SX 210 090 025 060	9 mm	25 mm	0,60 mm	0,80 mm	Cross



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
210 094 025 060	S 210 094 025 060	9,4 mm	25 mm	0,60 mm	0,60 mm	Straight



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
210 100 026 060	S 210 100 026 060	10 mm	26 mm	0,60 mm	0,60 mm	Straight
X 210 100 026 060	SX 210 100 026 060	10 mm	26 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
210 190 025 060	S 210 190 025 060	19 mm	25 mm	0,60 mm	0,60 mm	Straight
X 210 190 025 060	SX 210 190 025 060	19 mm	25 mm	0,60 mm	0,80 mm	Cross




**Saw Blades
Compatible with
DEPUY SYNTHES[®]**




Brand Code **220**




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 130 090 089	S 220 130 090 089	13 mm	90 mm	0,89 mm	0,89 mm	Straight
220 130 090 090	S 220 130 090 090	13 mm	90 mm	0,90 mm	0,90 mm	Straight
220 130 090 100	S 220 130 090 100	13 mm	90 mm	1,00 mm	1,00 mm	Straight
220 130 090 119	S 220 130 090 119	13 mm	90 mm	1,19 mm	1,19 mm	Straight
220 130 090 120	S 220 130 090 120	13 mm	90 mm	1,20 mm	1,20 mm	Straight
220 130 090 127	S 220 130 090 127	13 mm	90 mm	1,27 mm	1,27 mm	Straight
220 130 090 130	S 220 130 090 130	13 mm	90 mm	1,30 mm	1,30 mm	Straight
220 130 090 137	S 220 130 090 137	13 mm	90 mm	1,37 mm	1,37 mm	Straight
220 130 090 140	S 220 130 090 140	13 mm	90 mm	1,40 mm	1,40 mm	Straight
220 130 090 147	S 220 130 090 147	13 mm	90 mm	1,47 mm	1,47 mm	Straight
220 130 090 150	S 220 130 090 150	13 mm	90 mm	1,50 mm	1,50 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 190 090 089	S 220 190 090 089	19 mm	90 mm	0,89 mm	0,89 mm	Straight
220 190 090 090	S 220 190 090 090	19 mm	90 mm	0,90 mm	0,90 mm	Straight
220 190 090 100	S 220 190 090 100	19 mm	90 mm	1,00 mm	1,00 mm	Straight
220 190 090 119	S 220 190 090 119	19 mm	90 mm	1,19 mm	1,19 mm	Straight
220 190 090 120	S 220 190 090 120	19 mm	90 mm	1,20 mm	1,20 mm	Straight
220 190 090 127	S 220 190 090 127	19 mm	90 mm	1,27 mm	1,27 mm	Straight
220 190 090 130	S 220 190 090 130	19 mm	90 mm	1,30 mm	1,30 mm	Straight
220 190 090 137	S 220 190 090 137	19 mm	90 mm	1,37 mm	1,37 mm	Straight
220 190 090 140	S 220 190 090 140	19 mm	90 mm	1,40 mm	1,40 mm	Straight
220 190 090 147	S 220 190 090 147	19 mm	90 mm	1,47 mm	1,47 mm	Straight
220 190 090 150	S 220 190 090 150	19 mm	90 mm	1,50 mm	1,50 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 190 100 089	S 220 190 100 089	19 mm	100 mm	0,89 mm	0,89 mm	Straight
220 190 100 090	S 220 190 100 090	19 mm	100 mm	0,90 mm	0,90 mm	Straight
220 190 100 100	S 220 190 100 100	19 mm	100 mm	1,00 mm	1,00 mm	Straight
220 190 100 119	S 220 190 100 119	19 mm	100 mm	1,19 mm	1,19 mm	Straight
220 190 100 120	S 220 190 100 120	19 mm	100 mm	1,20 mm	1,20 mm	Straight
220 190 100 127	S 220 190 100 127	19 mm	100 mm	1,27 mm	1,27 mm	Straight
220 190 100 130	S 220 190 100 130	19 mm	100 mm	1,30 mm	1,30 mm	Straight
220 190 100 137	S 220 190 100 137	19 mm	100 mm	1,37 mm	1,37 mm	Straight
220 190 100 140	S 220 190 100 140	19 mm	100 mm	1,40 mm	1,40 mm	Straight
220 190 100 147	S 220 190 100 147	19 mm	100 mm	1,47 mm	1,47 mm	Straight
220 190 100 150	S 220 190 100 150	19 mm	100 mm	1,50 mm	1,50 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 195 090 089	S 220 195 090 089	19,5 mm	90 mm	0,89 mm	0,89 mm	Straight
220 195 090 090	S 220 195 090 090	19,5 mm	90 mm	0,90 mm	0,90 mm	Straight
220 195 090 100	S 220 195 090 100	19,5 mm	90 mm	1,00 mm	1,00 mm	Straight
220 195 090 119	S 220 195 090 119	19,5 mm	90 mm	1,19 mm	1,19 mm	Straight
220 195 090 120	S 220 195 090 120	19,5 mm	90 mm	1,20 mm	1,20 mm	Straight
220 195 090 127	S 220 195 090 127	19,5 mm	90 mm	1,27 mm	1,27 mm	Straight
220 195 090 130	S 220 195 090 130	19,5 mm	90 mm	1,30 mm	1,30 mm	Straight
220 195 090 137	S 220 195 090 137	19,5 mm	90 mm	1,37 mm	1,37 mm	Straight
220 195 090 140	S 220 195 090 140	19,5 mm	90 mm	1,40 mm	1,40 mm	Straight
220 195 090 147	S 220 195 090 147	19,5 mm	90 mm	1,47 mm	1,47 mm	Straight
220 195 090 150	S 220 195 090 150	19,5 mm	90 mm	1,50 mm	1,50 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 222 090 089	S 220 222 090 089	22,2 mm	90 mm	0,89 mm	0,89 mm	Straight
220 222 090 090	S 220 222 090 090	22,2 mm	90 mm	0,90 mm	0,90 mm	Straight
220 222 090 100	S 220 222 090 100	22,2 mm	90 mm	1,00 mm	1,00 mm	Straight
220 222 090 119	S 220 222 090 119	22,2 mm	90 mm	1,19 mm	1,19 mm	Straight
220 222 090 120	S 220 222 090 120	22,2 mm	90 mm	1,20 mm	1,20 mm	Straight
220 222 090 127	S 220 222 090 127	22,2 mm	90 mm	1,27 mm	1,27 mm	Straight
220 222 090 130	S 220 222 090 130	22,2 mm	90 mm	1,30 mm	1,30 mm	Straight
220 222 090 137	S 220 222 090 137	22,2 mm	90 mm	1,37 mm	1,37 mm	Straight
220 222 090 140	S 220 222 090 140	22,2 mm	90 mm	1,40 mm	1,40 mm	Straight
220 222 090 147	S 220 222 090 147	22,2 mm	90 mm	1,47 mm	1,47 mm	Straight
220 222 090 150	S 220 222 090 150	22,2 mm	90 mm	1,50 mm	1,50 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 250 090 089	S 220 250 090 089	25 mm	90 mm	0,89 mm	0,89 mm	Straight
220 250 090 090	S 220 250 090 090	25 mm	90 mm	0,90 mm	0,90 mm	Straight
220 250 090 100	S 220 250 090 100	25 mm	90 mm	1,00 mm	1,00 mm	Straight
220 250 090 119	S 220 250 090 119	25 mm	90 mm	1,19 mm	1,19 mm	Straight
220 250 090 120	S 220 250 090 120	25 mm	90 mm	1,20 mm	1,20 mm	Straight
220 250 090 127	S 220 250 090 127	25 mm	90 mm	1,27 mm	1,27 mm	Straight
220 250 090 130	S 220 250 090 130	25 mm	90 mm	1,30 mm	1,30 mm	Straight
220 250 090 137	S 220 250 090 137	25 mm	90 mm	1,37 mm	1,37 mm	Straight
220 250 090 140	S 220 250 090 140	25 mm	90 mm	1,40 mm	1,40 mm	Straight
220 250 090 147	S 220 250 090 147	25 mm	90 mm	1,47 mm	1,47 mm	Straight
220 250 090 150	S 220 250 090 150	25 mm	90 mm	1,50 mm	1,50 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 254 090 089	S 220 254 090 089	25,4 mm	90 mm	0,89 mm	0,89 mm	Straight
220 254 090 090	S 220 254 090 090	25,4 mm	90 mm	0,90 mm	0,90 mm	Straight
220 254 090 100	S 220 254 090 100	25,4 mm	90 mm	1,00 mm	1,00 mm	Straight
220 254 090 119	S 220 254 090 119	25,4 mm	90 mm	1,19 mm	1,19 mm	Straight
220 254 090 120	S 220 254 090 120	25,4 mm	90 mm	1,20 mm	1,20 mm	Straight
220 254 090 127	S 220 254 090 127	25,4 mm	90 mm	1,27 mm	1,27 mm	Straight
220 254 090 130	S 220 254 090 130	25,4 mm	90 mm	1,30 mm	1,30 mm	Straight
220 254 090 137	S 220 254 090 137	25,4 mm	90 mm	1,37 mm	1,37 mm	Straight
220 254 090 140	S 220 254 090 140	25,4 mm	90 mm	1,40 mm	1,40 mm	Straight
220 254 090 147	S 220 254 090 147	25,4 mm	90 mm	1,47 mm	1,47 mm	Straight
220 254 090 150	S 220 254 090 150	25,4 mm	90 mm	1,50 mm	1,50 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 125 090 127	S 220 125 090 127	12,5 mm	90 mm	1,27 mm	1,27 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 190 095 120	S 220 190 095 120	19 mm	95 mm	1,20 mm	1,20 mm	Straight
220 190 095 127	S 220 190 095 127	19 mm	95 mm	1,27 mm	1,27 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 230 090 127	S 220 230 090 127	23 mm	90 mm	1,27 mm	1,27 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 250 095 119	S 220 250 095 119	25 mm	95 mm	1,19 mm	1,19 mm	Straight
220 250 095 127	S 220 250 095 127	25 mm	95 mm	1,27 mm	1,27 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 090 025 060	S 220 090 025 060	9 mm	25 mm	0,60 mm	0,60 mm	Straight
X 220 090 025 060	SX 220 090 025 060	9 mm	25 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 070 025 060	S 220 070 025 060	7 mm	25 mm	0,60 mm	0,60 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 100 026 060	S 220 100 026 060	10 mm	26 mm	0,60 mm	0,60 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 140 026 060	S 220 140 026 060	14 mm	26 mm	0,60 mm	0,60 mm	Straight
X 220 140 026 060	SX 220 140 026 060	14 mm	26 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 140 051 060	S 220 140 051 060	14 mm	51 mm	0,60 mm	0,60 mm	Straight
X 220 140 051 060	SX 220 140 051 060	14 mm	51 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 141 050 060	S 220 141 050 060	14,1 mm	50 mm	0,60 mm	0,60 mm	Straight



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 130 045 127	S 220 130 045 127	13 mm	45 mm	1,27 mm	1,27 m	Straight



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 130 060 090	S 220 130 060 090	13 mm	60 mm	0,90 mm	0,90 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 150 060 090	S 220 150 060 090	15 mm	60 mm	0,90 mm	0,90 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 270 051 060	S 220 270 051 060	27 mm	51 mm	0,60 mm	0,60 mm	Straight
X 220 270 051 060	SX 220 270 051 060	27 mm	51 mm	0,60 mm	0,80 mm	Cross

Compatible Connections:


Compact AirDrive II + Power Drive 511.800, 511.801, Air Power Line 511.610, Battery Power Line 530.610, TRS 05.001.224, 05.001.223, 05.001.240, Battery Power Line II 530.710




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 140 020 060	S 220 140 020 060	14 mm	20 mm	0,60 mm	0,60 mm	Straight
X 220 140 020 060	SX 220 140 020 060	14 mm	20 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 480 070 060	S 220 480 070 060	48 mm	70 mm	0,60 mm	0,60 mm	Straight
X 220 480 070 060	SX 220 480 070 060	48 mm	70 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 490 070 090	S 220 490 070 090	49 mm	70 mm	0,60 mm	0,60 mm	Straight



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 125 077 090	S 220 125 077 090	12,5 mm	77 mm	0,90 mm	0,90 mm	Straight
220 125 077 100	S 220 125 077 100	12,5 mm	77 mm	1,00 mm	1,00 mm	Straight
220 125 077 120	S 220 125 077 120	12,5 mm	77 mm	1,20 mm	1,20 mm	Straight



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 125 089 090	S 220 125 089 090	12,5 mm	89 mm	0,90 mm	0,90 mm	Straight
220 125 089 100	S 220 125 089 100	12,5 mm	89 mm	1,00 mm	1,00 mm	Straight
220 125 089 120	S 220 125 089 120	12,5 mm	89 mm	1,20 mm	1,20 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 120 067 090	S 220 120 067 090	12 mm	67 mm	0,90 mm	0,90 mm	Straight
220 120 067 100	S 220 120 067 100	12 mm	67 mm	1,00 mm	1,00 mm	Straight
220 120 067 120	S 220 120 067 120	12 mm	67 mm	1,20 mm	1,20 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 130 068 090	S 220 130 068 090	13 mm	68 mm	0,90 mm	0,90 mm	Straight
220 130 068 100	S 220 130 068 100	13 mm	68 mm	1,00 mm	1,00 mm	Straight
220 130 068 120	S 220 130 068 120	13 mm	68 mm	1,20 mm	1,20 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 050 018 060	S 220 050 018 060	5 mm	18 mm	0,60 mm	0,60 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 060 015 040	S 220 060 015 040	6 mm	15 mm	0,40 mm	0,40 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 060 030 060	S 220 060 030 060	6 mm	30 mm	0,60 mm	0,60 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 070 015 060	S 220 070 015 060	7 mm	15 mm	0,60 mm	0,60 mm	Straight
X 220 070 015 060	SX 220 070 015 060	7 mm	15 mm	0,60 mm	0,80 mm	Cross



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 080 022 060	S 220 080 022 060	8 mm	22 mm	0,60 mm	0,60 mm	Straight
X 220 080 022 060	SX 220 080 022 060	8 mm	22 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 100 025 060	S 220 100 025 060	10 mm	25 mm	0,60 mm	0,60 mm	Straight
X 220 100 025 060	SX 220 100 025 060	10 mm	25 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
X 220 100 027 040	SX 220 100 027 040	10 mm	27 mm	0,40 mm	0,60 mm	Cross
220 100 027 060	S 220 100 027 060	10 mm	27 mm	0,60 mm	0,60 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 100 030 060	S 220 100 030 060	10 mm	30 mm	0,60 mm	0,60 mm	Straight
X 220 100 030 060	SX 220 100 030 060	10 mm	30 mm	0,60 mm	0,80 mm	Cross



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 100 031 060	S 220 100 031 060	10 mm	31 mm	0,60 mm	0,60 mm	Straight
X 220 100 031 060	SX 220 100 031 060	10 mm	31 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 100 049 060	S 220 100 049 060	10 mm	49 mm	0,60 mm	0,60 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
X 220 100 050 040	SX 220 100 050 040	10 mm	50 mm	0,40 mm	0,60 mm	Cross
220 100 050 060	S 220 100 050 060	10 mm	50 mm	0,60 mm	0,60 mm	Straight
X 220 100 050 060	SX 220 100 050 060	10 mm	50 mm	0,60 mm	0,80 mm	Cross



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 120 022 060	S 220 120 022 060	12 mm	22 mm	0,60 mm	0,60 mm	Straight
X 220 120 022 060	SX 220 120 022 060	12 mm	22 mm	0,60 mm	0,80 mm	Cross



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 120 040 060	S 220 120 040 060	12 mm	40 mm	0,60 mm	0,60 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 140 025 060	S 220 140 025 060	14 mm	25 mm	0,60 mm	0,60 mm	Straight
X 220 140 025 060	SX 220 140 025 060	14 mm	25 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 140 027 060	S 220 140 027 060	14 mm	27 mm	0,60 mm	0,60 mm	Straight



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 140 050 060	S 220 140 050 060	14 mm	50 mm	0,60 mm	0,60 mm	Straight



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 145 050 060	S 220 145 050 060	14,5 mm	50 mm	0,60 mm	0,60 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 198 050 060	S 220 198 050 060	19,8 mm	50 mm	0,60 mm	0,60 mm	Straight



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 200 049 060	S 220 200 049 060	20 mm	49 mm	0,60 mm	0,60 mm	Straight



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
X 220 200 050 040	SX 220 200 050 040	20 mm	50 mm	0,40 mm	0,60 mm	Cross



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 270 050 060	S 220 270 050 060	27 mm	50 mm	0,60 mm	0,60 mm	Straight
X 220 270 050 060	SX 220 270 050 060	27 mm	50 mm	0,60 mm	0,80 mm	Cross




Saw Blades Compatible with **ZIMMER®**




Brand Code **220**


For Reciprocating Saws
223




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 130 090 089	S 220 130 090 089	13 mm	90 mm	0,89 mm	0,89 mm	Straight
220 130 090 090	S 220 130 090 090	13 mm	90 mm	0,90 mm	0,90 mm	Straight
220 130 090 100	S 220 130 090 100	13 mm	90 mm	1,00 mm	1,00 mm	Straight
220 130 090 119	S 220 130 090 119	13 mm	90 mm	1,19 mm	1,19 mm	Straight
220 130 090 120	S 220 130 090 120	13 mm	90 mm	1,20 mm	1,20 mm	Straight
220 130 090 127	S 220 130 090 127	13 mm	90 mm	1,27 mm	1,27 mm	Straight
220 130 090 130	S 220 130 090 130	13 mm	90 mm	1,30 mm	1,30 mm	Straight
220 130 090 137	S 220 130 090 137	13 mm	90 mm	1,37 mm	1,37 mm	Straight
220 130 090 140	S 220 130 090 140	13 mm	90 mm	1,40 mm	1,40 mm	Straight
220 130 090 147	S 220 130 090 147	13 mm	90 mm	1,47 mm	1,47 mm	Straight
220 130 090 150	S 220 130 090 150	13 mm	90 mm	1,50 mm	1,50 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 190 090 089	S 220 190 090 089	19 mm	90 mm	0,89 mm	0,89 mm	Straight
220 190 090 090	S 220 190 090 090	19 mm	90 mm	0,90 mm	0,90 mm	Straight
220 190 090 100	S 220 190 090 100	19 mm	90 mm	1,00 mm	1,00 mm	Straight
220 190 090 119	S 220 190 090 119	19 mm	90 mm	1,19 mm	1,19 mm	Straight
220 190 090 120	S 220 190 090 120	19 mm	90 mm	1,20 mm	1,20 mm	Straight
220 190 090 127	S 220 190 090 127	19 mm	90 mm	1,27 mm	1,27 mm	Straight
220 190 090 130	S 220 190 090 130	19 mm	90 mm	1,30 mm	1,30 mm	Straight
220 190 090 137	S 220 190 090 137	19 mm	90 mm	1,37 mm	1,37 mm	Straight
220 190 090 140	S 220 190 090 140	19 mm	90 mm	1,40 mm	1,40 mm	Straight
220 190 090 147	S 220 190 090 147	19 mm	90 mm	1,47 mm	1,47 mm	Straight
220 190 090 150	S 220 190 090 150	19 mm	90 mm	1,50 mm	1,50 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 190 100 089	S 220 190 100 089	19 mm	100 mm	0,89 mm	0,89 mm	Straight
220 190 100 090	S 220 190 100 090	19 mm	100 mm	0,90 mm	0,90 mm	Straight
220 190 100 100	S 220 190 100 100	19 mm	100 mm	1,00 mm	1,00 mm	Straight
220 190 100 119	S 220 190 100 119	19 mm	100 mm	1,19 mm	1,19 mm	Straight
220 190 100 120	S 220 190 100 120	19 mm	100 mm	1,20 mm	1,20 mm	Straight
220 190 100 127	S 220 190 100 127	19 mm	100 mm	1,27 mm	1,27 mm	Straight
220 190 100 130	S 220 190 100 130	19 mm	100 mm	1,30 mm	1,30 mm	Straight
220 190 100 137	S 220 190 100 137	19 mm	100 mm	1,37 mm	1,37 mm	Straight
220 190 100 140	S 220 190 100 140	19 mm	100 mm	1,40 mm	1,40 mm	Straight
220 190 100 147	S 220 190 100 147	19 mm	100 mm	1,47 mm	1,47 mm	Straight
220 190 100 150	S 220 190 100 150	19 mm	100 mm	1,50 mm	1,50 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 195 090 089	S 220 195 090 089	19,5 mm	90 mm	0,89 mm	0,89 mm	Straight
220 195 090 090	S 220 195 090 090	19,5 mm	90 mm	0,90 mm	0,90 mm	Straight
220 195 090 100	S 220 195 090 100	19,5 mm	90 mm	1,00 mm	1,00 mm	Straight
220 195 090 119	S 220 195 090 119	19,5 mm	90 mm	1,19 mm	1,19 mm	Straight
220 195 090 120	S 220 195 090 120	19,5 mm	90 mm	1,20 mm	1,20 mm	Straight
220 195 090 127	S 220 195 090 127	19,5 mm	90 mm	1,27 mm	1,27 mm	Straight
220 195 090 130	S 220 195 090 130	19,5 mm	90 mm	1,30 mm	1,30 mm	Straight
220 195 090 137	S 220 195 090 137	19,5 mm	90 mm	1,37 mm	1,37 mm	Straight
220 195 090 140	S 220 195 090 140	19,5 mm	90 mm	1,40 mm	1,40 mm	Straight
220 195 090 147	S 220 195 090 147	19,5 mm	90 mm	1,47 mm	1,47 mm	Straight
220 195 090 150	S 220 195 090 150	19,5 mm	90 mm	1,50 mm	1,50 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 222 090 089	S 220 222 090 089	22,2 mm	90 mm	0,89 mm	0,89 mm	Straight
220 222 090 090	S 220 222 090 090	22,2 mm	90 mm	0,90 mm	0,90 mm	Straight
220 222 090 100	S 220 222 090 100	22,2 mm	90 mm	1,00 mm	1,00 mm	Straight
220 222 090 119	S 220 222 090 119	22,2 mm	90 mm	1,19 mm	1,19 mm	Straight
220 222 090 120	S 220 222 090 120	22,2 mm	90 mm	1,20 mm	1,20 mm	Straight
220 222 090 127	S 220 222 090 127	22,2 mm	90 mm	1,27 mm	1,27 mm	Straight
220 222 090 130	S 220 222 090 130	22,2 mm	90 mm	1,30 mm	1,30 mm	Straight
220 222 090 137	S 220 222 090 137	22,2 mm	90 mm	1,37 mm	1,37 mm	Straight
220 222 090 140	S 220 222 090 140	22,2 mm	90 mm	1,40 mm	1,40 mm	Straight
220 222 090 147	S 220 222 090 147	22,2 mm	90 mm	1,47 mm	1,47 mm	Straight
220 222 090 150	S 220 222 090 150	22,2 mm	90 mm	1,50 mm	1,50 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 250 090 089	S 220 250 090 089	25 mm	90 mm	0,89 mm	0,89 mm	Straight
220 250 090 090	S 220 250 090 090	25 mm	90 mm	0,90 mm	0,90 mm	Straight
220 250 090 100	S 220 250 090 100	25 mm	90 mm	1,00 mm	1,00 mm	Straight
220 250 090 119	S 220 250 090 119	25 mm	90 mm	1,19 mm	1,19 mm	Straight
220 250 090 120	S 220 250 090 120	25 mm	90 mm	1,20 mm	1,20 mm	Straight
220 250 090 127	S 220 250 090 127	25 mm	90 mm	1,27 mm	1,27 mm	Straight
220 250 090 130	S 220 250 090 130	25 mm	90 mm	1,30 mm	1,30 mm	Straight
220 250 090 137	S 220 250 090 137	25 mm	90 mm	1,37 mm	1,37 mm	Straight
220 250 090 140	S 220 250 090 140	25 mm	90 mm	1,40 mm	1,40 mm	Straight
220 250 090 147	S 220 250 090 147	25 mm	90 mm	1,47 mm	1,47 mm	Straight
220 250 090 150	S 220 250 090 150	25 mm	90 mm	1,50 mm	1,50 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 254 090 089	S 220 254 090 089	25,4 mm	90 mm	0,89 mm	0,89 mm	Straight
220 254 090 090	S 220 254 090 090	25,4 mm	90 mm	0,90 mm	0,90 mm	Straight
220 254 090 100	S 220 254 090 100	25,4 mm	90 mm	1,00 mm	1,00 mm	Straight
220 254 090 119	S 220 254 090 119	25,4 mm	90 mm	1,19 mm	1,19 mm	Straight
220 254 090 120	S 220 254 090 120	25,4 mm	90 mm	1,20 mm	1,20 mm	Straight
220 254 090 127	S 220 254 090 127	25,4 mm	90 mm	1,27 mm	1,27 mm	Straight
220 254 090 130	S 220 254 090 130	25,4 mm	90 mm	1,30 mm	1,30 mm	Straight
220 254 090 137	S 220 254 090 137	25,4 mm	90 mm	1,37 mm	1,37 mm	Straight
220 254 090 140	S 220 254 090 140	25,4 mm	90 mm	1,40 mm	1,40 mm	Straight
220 254 090 147	S 220 254 090 147	25,4 mm	90 mm	1,47 mm	1,47 mm	Straight
220 254 090 150	S 220 254 090 150	25,4 mm	90 mm	1,50 mm	1,50 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 125 090 127	S 220 125 090 127	12,5 mm	90 mm	1,27 mm	1,27 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 190 095 120	S 220 190 095 120	19 mm	95 mm	1,20 mm	1,20 mm	Straight
220 190 095 127	S 220 190 095 127	19 mm	95 mm	1,27 mm	1,27 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 230 090 127	S 220 230 090 127	23 mm	90 mm	1,27 mm	1,27 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 250 095 119	S 220 250 095 119	25 mm	95 mm	1,19 mm	1,19 mm	Straight
220 250 095 127	S 220 250 095 127	25 mm	95 mm	1,27 mm	1,27 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 090 025 060	S 220 090 025 060	9 mm	25 mm	0,60 mm	0,60 mm	Straight
X 220 090 025 060	SX 220 090 025 060	9 mm	25 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 070 025 060	S 220 070 025 060	7 mm	25 mm	0,60 mm	0,60 mm	Straight



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 100 026 060	S 220 100 026 060	10 mm	26 mm	0,60 mm	0,60 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 140 026 060	S 220 140 026 060	14 mm	26 mm	0,60 mm	0,60 mm	Straight
X 220 140 026 060	SX 220 140 026 060	14 mm	26 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 140 051 060	S 220 140 051 060	14 mm	51 mm	0,60 mm	0,60 mm	Straight
X 220 140 051 060	SX 220 140 051 060	14 mm	51 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 141 050 060	S 220 141 050 060	14,1 mm	50 mm	0,60 mm	0,60 mm	Straight



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 130 045 127	S 220 130 045 127	13 mm	45 mm	1,27 mm	1,27 mm	Straight



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 130 060 090	S 220 130 060 090	13 mm	60 mm	0,90 mm	0,90 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 150 060 090	S 220 150 060 090	15 mm	60 mm	0,90 mm	0,90 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 270 051 060	S 220 270 051 060	27 mm	51 mm	0,60 mm	0,60 mm	Straight
X 220 270 051 060	SX 220 270 051 060	27 mm	51 mm	0,60 mm	0,80 mm	Cross

Zimmer®


Compatible Connections:
Zimmer Universal Power System,
89-8509-450-60




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 140 020 060	S 220 140 020 060	14 mm	20 mm	0,60 mm	0,60 mm	Straight
X 220 140 020 060	SX 220 140 020 060	14 mm	20 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 480 070 060	S 220 480 070 060	48 mm	70 mm	0,60 mm	0,60 mm	Straight
X 220 480 070 060	SX 220 480 070 060	48 mm	70 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
220 490 070 090	S 220 490 070 090	49 mm	70 mm	0,90 mm	0,90 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
223 125 062 090	S 223 125 062 090	12,5 mm	62 mm	0,90 mm	0,90 mm	Straight
223 125 062 100	S 223 125 062 100	12,5 mm	62 mm	1,00 mm	1,00 mm	Straight
223 125 062 120	S 223 125 062 120	12,5 mm	62 mm	1,20 mm	1,20 mm	Straight



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
223 125 073 090	S 223 125 073 090	12,5 mm	73 mm	0,90 mm	0,90 mm	Straight
223 125 073 100	S 223 125 073 100	12,5 mm	73 mm	1,00 mm	1,00 mm	Straight
223 125 073 120	S 223 125 073 120	12,5 mm	73 mm	1,20 mm	1,20 mm	Straight



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
223 125 075 090	S 223 125 075 090	12,5 mm	75 mm	0,90 mm	0,90 mm	Straight
223 125 075 100	S 223 125 075 100	12,5 mm	75 mm	1,00 mm	1,00 mm	Straight
223 125 075 120	S 223 125 075 120	12,5 mm	75 mm	1,20 mm	1,20 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
223 125 089 090	S 223 125 089 090	12,5 mm	89 mm	0,90 mm	0,90 mm	Straight
223 125 089 100	S 223 125 089 100	12,5 mm	89 mm	1,00 mm	1,00 mm	Straight
223 125 089 120	S 223 125 089 120	12,5 mm	89 mm	1,20 mm	1,20 mm	Straight

Zimmer®

Compatible Connections:

Zimmer Universal Power System,
89-8509-451-20




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
223 135 076 090	S 223 135 076 090	13,5 mm	76 mm	0,90 mm	0,90 mm	Straight
223 135 076 100	S 223 135 076 100	13,5 mm	76 mm	1,00 mm	1,00 mm	Straight
223 135 076 120	S 223 135 076 120	13,5 mm	76 mm	1,20 mm	1,20 mm	Straight




**Saw Blades
Compatible with
BRASSELER[®] / TAVA[®]**




Brand Code **260**




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
260 130 090 089	S 260 130 090 089	13 mm	90 mm	0,89 mm	0,89 mm	Straight
260 130 090 090	S 260 130 090 090	13 mm	90 mm	0,90 mm	0,90 mm	Straight
260 130 090 100	S 260 130 090 100	13 mm	90 mm	1,00 mm	1,00 mm	Straight
260 130 090 119	S 260 130 090 119	13 mm	90 mm	1,19 mm	1,19 mm	Straight
260 130 090 120	S 260 130 090 120	13 mm	90 mm	1,20 mm	1,20 mm	Straight
260 130 090 127	S 260 130 090 127	13 mm	90 mm	1,27 mm	1,27 mm	Straight
260 130 090 130	S 260 130 090 130	13 mm	90 mm	1,30 mm	1,30 mm	Straight
260 130 090 137	S 260 130 090 137	13 mm	90 mm	1,37 mm	1,37 mm	Straight
260 130 090 140	S 260 130 090 140	13 mm	90 mm	1,40 mm	1,40 mm	Straight
260 130 090 147	S 260 130 090 147	13 mm	90 mm	1,47 mm	1,47 mm	Straight
260 130 090 150	S 260 130 090 150	13 mm	90 mm	1,50 mm	1,50 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
260 190 090 089	S 260 190 090 089	19 mm	90 mm	0,89 mm	0,89 mm	Straight
260 190 090 090	S 260 190 090 090	19 mm	90 mm	0,90 mm	0,90 mm	Straight
260 190 090 100	S 260 190 090 100	19 mm	90 mm	1,00 mm	1,00 mm	Straight
260 190 090 119	S 260 190 090 119	19 mm	90 mm	1,19 mm	1,19 mm	Straight
260 190 090 120	S 260 190 090 120	19 mm	90 mm	1,20 mm	1,20 mm	Straight
260 190 090 127	S 260 190 090 127	19 mm	90 mm	1,27 mm	1,27 mm	Straight
260 190 090 130	S 260 190 090 130	19 mm	90 mm	1,30 mm	1,30 mm	Straight
260 190 090 137	S 260 190 090 137	19 mm	90 mm	1,37 mm	1,37 mm	Straight
260 190 090 140	S 260 190 090 140	19 mm	90 mm	1,40 mm	1,40 mm	Straight
260 190 090 147	S 260 190 090 147	19 mm	90 mm	1,47 mm	1,47 mm	Straight
260 190 090 150	S 260 190 090 150	19 mm	90 mm	1,50 mm	1,50 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
260 190 100 089	S 260 190 100 089	19 mm	100 mm	0,89 mm	0,89 mm	Straight
260 190 100 090	S 260 190 100 090	19 mm	100 mm	0,90 mm	0,90 mm	Straight
260 190 100 100	S 260 190 100 100	19 mm	100 mm	1,00 mm	1,00 mm	Straight
260 190 100 119	S 260 190 100 119	19 mm	100 mm	1,19 mm	1,19 mm	Straight
260 190 100 120	S 260 190 100 120	19 mm	100 mm	1,20 mm	1,20 mm	Straight
260 190 100 127	S 260 190 100 127	19 mm	100 mm	1,27 mm	1,27 mm	Straight
260 190 100 130	S 260 190 100 130	19 mm	100 mm	1,30 mm	1,30 mm	Straight
260 190 100 137	S 260 190 100 137	19 mm	100 mm	1,37 mm	1,37 mm	Straight
260 190 100 140	S 260 190 100 140	19 mm	100 mm	1,40 mm	1,40 mm	Straight
260 190 100 147	S 260 190 100 147	19 mm	100 mm	1,47 mm	1,47 mm	Straight
260 190 100 150	S 260 190 100 150	19 mm	100 mm	1,50 mm	1,50 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
260 195 090 089	S 260 195 090 089	19,5 mm	90 mm	0,89 mm	0,89 mm	Straight
260 195 090 090	S 260 195 090 090	19,5 mm	90 mm	0,90 mm	0,90 mm	Straight
260 195 090 100	S 260 195 090 100	19,5 mm	90 mm	1,00 mm	1,00 mm	Straight
260 195 090 119	S 260 195 090 119	19,5 mm	90 mm	1,19 mm	1,19 mm	Straight
260 195 090 120	S 260 195 090 120	19,5 mm	90 mm	1,20 mm	1,20 mm	Straight
260 195 090 127	S 260 195 090 127	19,5 mm	90 mm	1,27 mm	1,27 mm	Straight
260 195 090 130	S 260 195 090 130	19,5 mm	90 mm	1,30 mm	1,30 mm	Straight
260 195 090 137	S 260 195 090 137	19,5 mm	90 mm	1,37 mm	1,37 mm	Straight
260 195 090 140	S 260 195 090 140	19,5 mm	90 mm	1,40 mm	1,40 mm	Straight
260 195 090 147	S 260 195 090 147	19,5 mm	90 mm	1,47 mm	1,47 mm	Straight
260 195 090 150	S 260 195 090 150	19,5 mm	90 mm	1,50 mm	1,50 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
260 220 090 089	S 260 220 090 089	22 mm	90 mm	0,89 mm	0,89 mm	Straight
260 220 090 090	S 260 220 090 090	22 mm	90 mm	0,90 mm	0,90 mm	Straight
260 220 090 100	S 260 220 090 100	22 mm	90 mm	1,00 mm	1,00 mm	Straight
260 220 090 119	S 260 220 090 119	22 mm	90 mm	1,19 mm	1,19 mm	Straight
260 220 090 120	S 260 220 090 120	22 mm	90 mm	1,20 mm	1,20 mm	Straight
260 220 090 127	S 260 220 090 127	22 mm	90 mm	1,27 mm	1,27 mm	Straight
260 220 090 130	S 260 220 090 130	22 mm	90 mm	1,30 mm	1,30 mm	Straight
260 220 090 137	S 260 220 090 137	22 mm	90 mm	1,37 mm	1,37 mm	Straight
260 220 090 140	S 260 220 090 140	22 mm	90 mm	1,40 mm	1,40 mm	Straight
260 220 090 147	S 260 220 090 147	22 mm	90 mm	1,47 mm	1,47 mm	Straight
260 220 090 150	S 260 220 090 150	22 mm	90 mm	1,50 mm	1,50 mm	Straight



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
260 250 090 089	S 260 250 090 089	25 mm	90 mm	0,89 mm	0,89 mm	Straight
260 250 090 090	S 260 250 090 090	25 mm	90 mm	0,90 mm	0,90 mm	Straight
260 250 090 100	S 260 250 090 100	25 mm	90 mm	1,00 mm	1,00 mm	Straight
260 250 090 119	S 260 250 090 119	25 mm	90 mm	1,19 mm	1,19 mm	Straight
260 250 090 120	S 260 250 090 120	25 mm	90 mm	1,20 mm	1,20 mm	Straight
260 250 090 127	S 260 250 090 127	25 mm	90 mm	1,27 mm	1,27 mm	Straight
260 250 090 130	S 260 250 090 130	25 mm	90 mm	1,30 mm	1,30 mm	Straight
260 250 090 137	S 260 250 090 137	25 mm	90 mm	1,37 mm	1,37 mm	Straight
260 250 090 140	S 260 250 090 140	25 mm	90 mm	1,40 mm	1,40 mm	Straight
260 250 090 147	S 260 250 090 147	25 mm	90 mm	1,47 mm	1,47 mm	Straight
260 250 090 150	S 260 250 090 150	25 mm	90 mm	1,50 mm	1,50 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
260 254 090 089	S 260 254 090 089	25,4 mm	90 mm	0,89 mm	0,89 mm	Straight
260 254 090 090	S 260 254 090 090	25,4 mm	90 mm	0,90 mm	0,90 mm	Straight
260 254 090 100	S 260 254 090 100	25,4 mm	90 mm	1,00 mm	1,00 mm	Straight
260 254 090 119	S 260 254 090 119	25,4 mm	90 mm	1,19 mm	1,19 mm	Straight
260 254 090 120	S 260 254 090 120	25,4 mm	90 mm	1,20 mm	1,20 mm	Straight
260 254 090 127	S 260 254 090 127	25,4 mm	90 mm	1,27 mm	1,27 mm	Straight
260 254 090 130	S 260 254 090 130	25,4 mm	90 mm	1,30 mm	1,30 mm	Straight
260 254 090 137	S 260 254 090 137	25,4 mm	90 mm	1,37 mm	1,37 mm	Straight
260 254 090 140	S 260 254 090 140	25,4 mm	90 mm	1,40 mm	1,40 mm	Straight
260 254 090 147	S 260 254 090 147	25,4 mm	90 mm	1,47 mm	1,47 mm	Straight
260 254 090 150	S 260 254 090 150	25,4 mm	90 mm	1,50 mm	1,50 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
260 055 025 060	S 260 055 025 060	5,5 mm	25 mm	0,60 mm	0,60 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
260 090 025 060	S 260 090 025 060	9 mm	25 mm	0,60 mm	0,60 mm	Straight
X 260 090 025 060	SX 260 090 025 060	9 mm	25 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
260 070 025 060	S 260 070 025 060	7 mm	25 mm	0,60 mm	0,60 mm	Straight
X 260 070 025 060	SX 260 070 025 060	7 mm	25 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
260 090 040 060	S 260 090 040 060	9 mm	40 mm	0,60 mm	0,60 mm	Straight
X 260 090 040 060	SX 260 090 040 060	9 mm	40 mm	0,60 mm	0,80 mm	Cross



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
260 280 058 060	S 260 280 058 060	28 mm	58 mm	0,60 mm	0,60 mm	Straight
X 260 280 058 060	SX 260 280 058 060	28 mm	58 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
260 520 030 060	S 260 520 030 060	52 mm	30 mm	0,60 mm	0,60 mm	Straight
X 260 520 030 060	SX 260 520 030 060	52 mm	30 mm	0,60 mm	0,80 mm	Cross
260 520 030 090	S 260 520 030 090	52 mm	30 mm	0,90 mm	0,90 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
260 520 040 060	S 260 520 040 060	52 mm	40 mm	0,60 mm	0,60 mm	Straight
X 260 520 040 060	SX 260 520 040 060	52 mm	40 mm	0,60 mm	0,80 mm	Cross
260 520 040 090	S 260 520 040 090	52 mm	40 mm	0,90 mm	0,90 mm	Straight




Saw Blades Compatible with B BRAUN AESCULAP®




Brand Code **320**
330




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
320 130 090 089	S 320 130 090 089	13 mm	90 mm	0,89 mm	0,89 mm	Straight
320 130 090 090	S 320 130 090 090	13 mm	90 mm	0,90 mm	0,90 mm	Straight
320 130 090 100	S 320 130 090 100	13 mm	90 mm	1,00 mm	1,00 mm	Straight
320 130 090 119	S 320 130 090 119	13 mm	90 mm	1,19 mm	1,19 mm	Straight
320 130 090 120	S 320 130 090 120	13 mm	90 mm	1,20 mm	1,20 mm	Straight
320 130 090 127	S 320 130 090 127	13 mm	90 mm	1,27 mm	1,27 mm	Straight
320 130 090 130	S 320 130 090 130	13 mm	90 mm	1,30 mm	1,30 mm	Straight
320 130 090 137	S 320 130 090 137	13 mm	90 mm	1,37 mm	1,37 mm	Straight
320 130 090 140	S 320 130 090 140	13 mm	90 mm	1,40 mm	1,40 mm	Straight
320 130 090 147	S 320 130 090 147	13 mm	90 mm	1,47 mm	1,47 mm	Straight
320 130 090 150	S 320 130 090 150	13 mm	90 mm	1,50 mm	1,50 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
320 190 090 089	S 320 190 090 089	19 mm	90 mm	0,89 mm	0,89 mm	Straight
320 190 090 090	S 320 190 090 090	19 mm	90 mm	0,90 mm	0,90 mm	Straight
320 190 090 100	S 320 190 090 100	19 mm	90 mm	1,00 mm	1,00 mm	Straight
320 190 090 119	S 320 190 090 119	19 mm	90 mm	1,19 mm	1,19 mm	Straight
320 190 090 120	S 320 190 090 120	19 mm	90 mm	1,20 mm	1,20 mm	Straight
320 190 090 127	S 320 190 090 127	19 mm	90 mm	1,27 mm	1,27 mm	Straight
320 190 090 130	S 320 190 090 130	19 mm	90 mm	1,30 mm	1,30 mm	Straight
320 190 090 137	S 320 190 090 137	19 mm	90 mm	1,37 mm	1,37 mm	Straight
320 190 090 140	S 320 190 090 140	19 mm	90 mm	1,40 mm	1,40 mm	Straight
320 190 090 147	S 320 190 090 147	19 mm	90 mm	1,47 mm	1,47 mm	Straight
320 190 090 150	S 320 190 090 150	19 mm	90 mm	1,50 mm	1,50 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
320 190 100 089	S 320 190 100 089	19 mm	100 mm	0,89 mm	0,89 mm	Straight
320 190 100 090	S 320 190 100 090	19 mm	100 mm	0,90 mm	0,90 mm	Straight
320 190 100 100	S 320 190 100 100	19 mm	100 mm	1,00 mm	1,00 mm	Straight
320 190 100 119	S 320 190 100 119	19 mm	100 mm	1,19 mm	1,19 mm	Straight
320 190 100 120	S 320 190 100 120	19 mm	100 mm	1,20 mm	1,20 mm	Straight
320 190 100 127	S 320 190 100 127	19 mm	100 mm	1,27 mm	1,27 mm	Straight
320 190 100 130	S 320 190 100 130	19 mm	100 mm	1,30 mm	1,30 mm	Straight
320 190 100 137	S 320 190 100 137	19 mm	100 mm	1,37 mm	1,37 mm	Straight
320 190 100 140	S 320 190 100 140	19 mm	100 mm	1,40 mm	1,40 mm	Straight
320 190 100 147	S 320 190 100 147	19 mm	100 mm	1,47 mm	1,47 mm	Straight
320 190 100 150	S 320 190 100 150	19 mm	100 mm	1,50 mm	1,50 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
320 195 090 089	S 320 195 090 089	19,5 mm	90 mm	0,89 mm	0,89 mm	Straight
320 195 090 090	S 320 195 090 090	19,5 mm	90 mm	0,90 mm	0,90 mm	Straight
320 195 090 100	S 320 195 090 100	19,5 mm	90 mm	1,00 mm	1,00 mm	Straight
320 195 090 119	S 320 195 090 119	19,5 mm	90 mm	1,19 mm	1,19 mm	Straight
320 195 090 120	S 320 195 090 120	19,5 mm	90 mm	1,20 mm	1,20 mm	Straight
320 195 090 127	S 320 195 090 127	19,5 mm	90 mm	1,27 mm	1,27 mm	Straight
320 195 090 130	S 320 195 090 130	19,5 mm	90 mm	1,30 mm	1,30 mm	Straight
320 195 090 137	S 320 195 090 137	19,5 mm	90 mm	1,37 mm	1,37 mm	Straight
320 195 090 140	S 320 195 090 140	19,5 mm	90 mm	1,40 mm	1,40 mm	Straight
320 195 090 147	S 320 195 090 147	19,5 mm	90 mm	1,47 mm	1,47 mm	Straight
320 195 090 150	S 320 195 090 150	19,5 mm	90 mm	1,50 mm	1,50 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
320 220 090 089	S 320 220 090 089	22 mm	90 mm	0,89 mm	0,89 mm	Straight
320 220 090 090	S 320 220 090 090	22 mm	90 mm	0,90 mm	0,90 mm	Straight
320 220 090 100	S 320 220 090 100	22 mm	90 mm	1,00 mm	1,00 mm	Straight
320 220 090 119	S 320 220 090 119	22 mm	90 mm	1,19 mm	1,19 mm	Straight
320 220 090 120	S 320 220 090 120	22 mm	90 mm	1,20 mm	1,20 mm	Straight
320 220 090 127	S 320 220 090 127	22 mm	90 mm	1,27 mm	1,27 mm	Straight
320 220 090 130	S 320 220 090 130	22 mm	90 mm	1,30 mm	1,30 mm	Straight
320 220 090 137	S 320 220 090 137	22 mm	90 mm	1,37 mm	1,37 mm	Straight
320 220 090 140	S 320 220 090 140	22 mm	90 mm	1,40 mm	1,40 mm	Straight
320 220 090 147	S 320 220 090 147	22 mm	90 mm	1,47 mm	1,47 mm	Straight
320 220 090 150	S 320 220 090 150	22 mm	90 mm	1,50 mm	1,50 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
320 250 090 089	S 320 250 090 089	25 mm	90 mm	0,89 mm	0,89 mm	Straight
320 250 090 090	S 320 250 090 090	25 mm	90 mm	0,90 mm	0,90 mm	Straight
320 250 090 100	S 320 250 090 100	25 mm	90 mm	1,00 mm	1,00 mm	Straight
320 250 090 119	S 320 250 090 119	25 mm	90 mm	1,19 mm	1,19 mm	Straight
320 250 090 120	S 320 250 090 120	25 mm	90 mm	1,20 mm	1,20 mm	Straight
320 250 090 127	S 320 250 090 127	25 mm	90 mm	1,27 mm	1,27 mm	Straight
320 250 090 130	S 320 250 090 130	25 mm	90 mm	1,30 mm	1,30 mm	Straight
320 250 090 137	S 320 250 090 137	25 mm	90 mm	1,37 mm	1,37 mm	Straight
320 250 090 140	S 320 250 090 140	25 mm	90 mm	1,40 mm	1,40 mm	Straight
320 250 090 147	S 320 250 090 147	25 mm	90 mm	1,47 mm	1,47 mm	Straight
320 250 090 150	S 320 250 090 150	25 mm	90 mm	1,50 mm	1,50 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
320 254 090 089	S 320 254 090 089	25,4 mm	90 mm	0,89 mm	0,89 mm	Straight
320 254 090 090	S 320 254 090 090	25,4 mm	90 mm	0,90 mm	0,90 mm	Straight
320 254 090 100	S 320 254 090 100	25,4 mm	90 mm	1,00 mm	1,00 mm	Straight
320 254 090 119	S 320 254 090 119	25,4 mm	90 mm	1,19 mm	1,19 mm	Straight
320 254 090 120	S 320 254 090 120	25,4 mm	90 mm	1,20 mm	1,20 mm	Straight
320 254 090 127	S 320 254 090 127	25,4 mm	90 mm	1,27 mm	1,27 mm	Straight
320 254 090 130	S 320 254 090 130	25,4 mm	90 mm	1,30 mm	1,30 mm	Straight
320 254 090 137	S 320 254 090 137	25,4 mm	90 mm	1,37 mm	1,37 mm	Straight
320 254 090 140	S 320 254 090 140	25,4 mm	90 mm	1,40 mm	1,40 mm	Straight
320 254 090 147	S 320 254 090 147	25,4 mm	90 mm	1,47 mm	1,47 mm	Straight
320 254 090 150	S 320 254 090 150	25,4 mm	90 mm	1,50 mm	1,50 mm	Straight



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
320 130 070 100	S 320 130 070 100	13 mm	70 mm	1,00 mm	1,00 mm	Straight



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
320 090 075 127	S 320 090 075 127	9 mm	75 mm	1,27 mm	1,27 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
320 195 036 090	S 320 195 036 090	19,5 mm	36 mm	0,90 mm	0,90 mm	Straight



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
320 150 050 090	S 320 150 050 090	15 mm	50 mm	0,90 mm	0,90 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
320 250 048 090	S 320 250 048 090	25 mm	48 mm	0,90 mm	0,90 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
320 300 048 090	S 320 300 048 090	30 mm	48 mm	0,90 mm	0,90 mm	Straight


B Braun Aesculap®
Compatible Connections:
GA673,
GA668,
GB660R




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
320 195 065 090	S 320 195 065 090	19,5 mm	65 mm	0,90 mm	0,90 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
320 250 065 090	S 320 250 065 090	25 mm	65 mm	0,90 mm	0,90 mm	Straight



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
320 320 067 090	S 320 320 067 090	32 mm	67 mm	0,90 mm	0,90 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
320 410 067 090	S 320 410 067 090	41 mm	67 mm	0,90 mm	0,90 mm	Straight



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
320 420 050 090	S 320 420 050 090	42 mm	50 mm	0,90 mm	0,90 mm	Straight



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
320 055 025 060	S 320 055 025 060	5,5 mm	25 mm	0,60 mm	0,60 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
320 710 035 090	S 320 710 035 090	71 mm	35 mm	0,90 mm	0,90 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
320 090 025 060	S 320 090 025 060	9 mm	25 mm	0,60 mm	0,60 mm	Straight
X 320 090 025 060	SX 320 090 025 060	9 mm	25 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
320 100 035 060	S 320 100 035 060	10 mm	35 mm	0,60 mm	0,60 mm	Straight
X 320 100 035 060	SX 320 100 035 060	10 mm	35 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
X 320 100 050 050	SX 320 100 050 050	10 mm	50 mm	0,50 mm	0,80 mm	Cross
320 100 050 080	S 320 100 050 080	10 mm	50 mm	0,80 mm	0,80 mm	Straight
320 100 050 090	S 320 100 050 090	10 mm	50 mm	0,90 mm	0,90 mm	Straight



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
320 150 026 060	S 320 150 026 060	15 mm	26 mm	0,60 mm	0,60 mm	Straight
X 320 150 026 060	SX 320 150 026 060	15 mm	26 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
320 150 045 060	S 320 150 045 060	15 mm	45 mm	0,60 mm	0,60 mm	Straight
X 320 150 045 060	SX 320 150 045 060	15 mm	45 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
320 150 050 090	S 320 150 050 090	15 mm	50 mm	0,90 mm	0,90 mm	Straight


B Braun Aesculap®
Compatible Connections:
GA673,
GA668,
GB660R




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
X 320 200 050 070	SX 320 200 050 070	20 mm	50 mm	0,70 mm	0,90 mm	Cross
320 200 050 090	S 320 200 050 090	20 mm	50 mm	0,90 mm	0,90 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
320 155 045 050	S 320 155 045 050	15 mm	45 mm	0,50 mm	0,50 mm	Straight



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
X 320 120 060 070	SX 320 120 060 070	12 mm	60 mm	0,70 mm	0,90 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
X 320 090 075 070	SX 320 090 075 070	9 mm	75 mm	0,70 mm	0,90 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
320 070 035 060	S 320 070 035 060	7 mm	35 mm	0,60 mm	0,60 mm	Straight
X 320 070 035 060	SX 320 070 035 060	7 mm	35 mm	0,60 mm	0,80 mm	Cross

Compatible Connections;


Acculan/Comfort, GA4040,
GA523, GA623, GA633, GA640,
GB125R, GB129, GB431, GB432




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
330 130 090 089	S 330 130 090 089	13 mm	90 mm	0,89 mm	0,89 mm	Straight
330 130 090 090	S 330 130 090 090	13 mm	90 mm	0,90 mm	0,90 mm	Straight
330 130 090 100	S 330 130 090 100	13 mm	90 mm	1,00 mm	1,00 mm	Straight
330 130 090 119	S 330 130 090 119	13 mm	90 mm	1,19 mm	1,19 mm	Straight
330 130 090 120	S 330 130 090 120	13 mm	90 mm	1,20 mm	1,20 mm	Straight
330 130 090 127	S 330 130 090 127	13 mm	90 mm	1,27 mm	1,27 mm	Straight
330 130 090 130	S 330 130 090 130	13 mm	90 mm	1,30 mm	1,30 mm	Straight
330 130 090 137	S 330 130 090 137	13 mm	90 mm	1,37 mm	1,37 mm	Straight
330 130 090 140	S 330 130 090 140	13 mm	90 mm	1,40 mm	1,40 mm	Straight



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
330 190 090 089	S 330 190 090 089	19 mm	90 mm	0,89 mm	0,89 mm	Straight
330 190 090 090	S 330 190 090 090	19 mm	90 mm	0,90 mm	0,90 mm	Straight
330 190 090 100	S 330 190 090 100	19 mm	90 mm	1,00 mm	1,00 mm	Straight
330 190 090 119	S 330 190 090 119	19 mm	90 mm	1,19 mm	1,19 mm	Straight
330 190 090 120	S 330 190 090 120	19 mm	90 mm	1,20 mm	1,20 mm	Straight
330 190 090 127	S 330 190 090 127	19 mm	90 mm	1,27 mm	1,27 mm	Straight
330 190 090 130	S 330 190 090 130	19 mm	90 mm	1,30 mm	1,30 mm	Straight
330 190 090 137	S 330 190 090 137	19 mm	90 mm	1,37 mm	1,37 mm	Straight
330 190 090 140	S 330 190 090 140	19 mm	90 mm	1,40 mm	1,40 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
330 190 100 089	S 330 190 100 089	19 mm	100 mm	0,89 mm	0,89 mm	Straight
330 190 100 090	S 330 190 100 090	19 mm	100 mm	0,90 mm	0,90 mm	Straight
330 190 100 100	S 330 190 100 100	19 mm	100 mm	1,00 mm	1,00 mm	Straight
330 190 100 119	S 330 190 100 119	19 mm	100 mm	1,19 mm	1,19 mm	Straight
330 190 100 120	S 330 190 100 120	19 mm	100 mm	1,20 mm	1,20 mm	Straight
330 190 100 127	S 330 190 100 127	19 mm	100 mm	1,27 mm	1,27 mm	Straight
330 190 100 130	S 330 190 100 130	19 mm	100 mm	1,30 mm	1,30 mm	Straight
330 190 100 137	S 330 190 100 137	19 mm	100 mm	1,37 mm	1,37 mm	Straight
330 190 100 140	S 330 190 100 140	19 mm	100 mm	1,40 mm	1,40 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
330 195 090 089	S 330 195 090 089	19,5 mm	90 mm	0,89 mm	0,89 mm	Straight
330 195 090 090	S 330 195 090 090	19,5 mm	90 mm	0,90 mm	0,90 mm	Straight
330 195 090 100	S 330 195 090 100	19,5 mm	90 mm	1,00 mm	1,00 mm	Straight
330 195 090 119	S 330 195 090 119	19,5 mm	90 mm	1,19 mm	1,19 mm	Straight
330 195 090 120	S 330 195 090 120	19,5 mm	90 mm	1,20 mm	1,20 mm	Straight
330 195 090 127	S 330 195 090 127	19,5 mm	90 mm	1,27 mm	1,27 mm	Straight
330 195 090 130	S 330 195 090 130	19,5 mm	90 mm	1,30 mm	1,30 mm	Straight
330 195 090 137	S 330 195 090 137	19,5 mm	90 mm	1,37 mm	1,37 mm	Straight
330 195 090 140	S 330 195 090 140	19,5 mm	90 mm	1,40 mm	1,40 mm	Straight

Compatible Connections:


Acculan/Comfort, GA4040,
GA523, GA623, GA633, GA640,
GB125R, GB129, GB431, GB432




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
330 220 090 089	S 330 220 090 089	22 mm	90 mm	0,89 mm	0,89 mm	Straight
330 220 090 090	S 330 220 090 090	22 mm	90 mm	0,90 mm	0,90 mm	Straight
330 220 090 100	S 330 220 090 100	22 mm	90 mm	1,00 mm	1,00 mm	Straight
330 220 090 119	S 330 220 090 119	22 mm	90 mm	1,19 mm	1,19 mm	Straight
330 220 090 120	S 330 220 090 120	22 mm	90 mm	1,20 mm	1,20 mm	Straight
330 220 090 127	S 330 220 090 127	22 mm	90 mm	1,27 mm	1,27 mm	Straight
330 220 090 130	S 330 220 090 130	22 mm	90 mm	1,30 mm	1,30 mm	Straight
330 220 090 137	S 330 220 090 137	22 mm	90 mm	1,37 mm	1,37 mm	Straight
330 220 090 140	S 330 220 090 140	22 mm	90 mm	1,40 mm	1,40 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
330 250 090 089	S 330 250 090 089	25 mm	90 mm	0,89 mm	0,89 mm	Straight
330 250 090 090	S 330 250 090 090	25 mm	90 mm	0,90 mm	0,90 mm	Straight
330 250 090 100	S 330 250 090 100	25 mm	90 mm	1,00 mm	1,00 mm	Straight
330 250 090 119	S 330 250 090 119	25 mm	90 mm	1,19 mm	1,19 mm	Straight
330 250 090 120	S 330 250 090 120	25 mm	90 mm	1,20 mm	1,20 mm	Straight
330 250 090 127	S 330 250 090 127	25 mm	90 mm	1,27 mm	1,27 mm	Straight
330 250 090 130	S 330 250 090 130	25 mm	90 mm	1,30 mm	1,30 mm	Straight
330 250 090 137	S 330 250 090 137	25 mm	90 mm	1,37 mm	1,37 mm	Straight
330 250 090 140	S 330 250 090 140	25 mm	90 mm	1,40 mm	1,40 mm	Straight



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
330 254 090 089	S 330 254 090 089	25,4 mm	90 mm	0,89 mm	0,89 mm	Straight
330 254 090 090	S 330 254 090 090	25,4 mm	90 mm	0,90 mm	0,90 mm	Straight
330 254 090 100	S 330 254 090 100	25,4 mm	90 mm	1,00 mm	1,00 mm	Straight
330 254 090 119	S 330 254 090 119	25,4 mm	90 mm	1,19 mm	1,19 mm	Straight
330 254 090 120	S 330 254 090 120	25,4 mm	90 mm	1,20 mm	1,20 mm	Straight
330 254 090 127	S 330 254 090 127	25,4 mm	90 mm	1,27 mm	1,27 mm	Straight
330 254 090 130	S 330 254 090 130	25,4 mm	90 mm	1,30 mm	1,30 mm	Straight
330 254 090 137	S 330 254 090 137	25,4 mm	90 mm	1,37 mm	1,37 mm	Straight
330 254 090 140	S 330 254 090 140	25,4 mm	90 mm	1,40 mm	1,40 mm	Straight



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
330 120 067 090	S 330 120 067 090	12 mm	67 mm	0,90 mm	0,90 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
330 220 046 090	S 330 220 046 090	22 mm	46 mm	0,90 mm	0,90 mm	Straight


Compatible Connections;
Acculan/Comfort, GA4040,
GA523, GA623, GA633, GA640,
GB125R, GB129, GB431, GB432




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
330 250 046 090	S 330 250 046 090	25 mm	46 mm	0,90 mm	0,90 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
330 250 067 090	S 330 250 067 090	25 mm	67 mm	0,90 mm	0,90 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
330 270 046 090	S 330 270 046 090	27 mm	46 mm	0,90 mm	0,90 mm	Straight



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
330 270 048 090	S 330 270 048 090	27 mm	48 mm	0,90 mm	0,90 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
330 380 045 090	S 330 380 045 090	38 mm	45 mm	0,90 mm	0,90 mm	Straight



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
330 095 030 060	S 330 095 030 060	9,5 mm	30 mm	0,60 mm	0,60 mm	Straight
X 330 095 030 060	SX 330 095 030 060	9,5 mm	30 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
330 130 043 060	S 330 130 043 060	13 mm	43 mm	0,60 mm	0,60 mm	Straight
X 330 130 043 060	SX 330 130 043 060	13 mm	43 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
330 130 048 060	S 330 130 048 060	13 mm	48 mm	0,60 mm	0,60 mm	Straight
X 330 130 048 060	SX 330 130 048 060	13 mm	48 mm	0,60 mm	0,80 mm	Cross

Compatible Connections;
Acculan/Comfort, GA4040,
GA523, GA623, GA633, GA640,
GB125R, GB129, GB431, GB432




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
330 140 024 060	S 330 140 024 060	14 mm	24 mm	0,60 mm	0,60 mm	Straight
S 330 140 024 060	SX 330 140 024 060	14 mm	24 mm	0,60 mm	0,80 mm	Cross




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
330 500 067 060	S 330 500 067 060	50 mm	67 mm	0,60 mm	0,60 mm	Straight
X 330 500 067 060	SX 330 500 067 060	50 mm	67 mm	0,60 mm	0,80 mm	Cross




**Saw Blades
Compatible with
BOJIN[®]**




Brand Code 180




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
180 100 090 120	S 180 100 090 120	10 mm	90 mm	1,20 mm	1,20 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
180 110 090 100	S 180 110 090 100	11 mm	90 mm	1,00 mm	1,00 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
180 130 090 090	S 180 130 090 090	13 mm	90 mm	0,90 mm	0,90 mm	Straight
180 130 090 100	S 180 130 090 100	13 mm	90 mm	1,00 mm	1,00 mm	Straight
180 130 090 120	S 180 130 090 120	13 mm	90 mm	1,20 mm	1,20 mm	Straight
180 130 090 127	S 180 130 090 127	13 mm	90 mm	1,27 mm	1,27 mm	Straight
180 130 090 130	S 180 130 090 130	13 mm	90 mm	1,30 mm	1,30 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
180 190 090 090	S 180 190 090 090	19 mm	90 mm	0,90 mm	0,90 mm	Straight
180 190 090 100	S 180 190 090 100	19 mm	90 mm	1,00 mm	1,00 mm	Straight
180 190 090 120	S 180 190 090 120	19 mm	90 mm	1,20 mm	1,20 mm	Straight
180 190 090 127	S 180 190 090 127	19 mm	90 mm	1,27 mm	1,27 mm	Straight
180 190 090 130	S 180 190 090 130	19 mm	90 mm	1,30 mm	1,30 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
180 190 100 100	S 180 190 100 100	19 mm	100 mm	1,00 mm	1,00 mm	Straight
180 190 100 120	S 180 190 100 120	19 mm	100 mm	1,20 mm	1,20 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
180 220 070 100	S 180 220 070 100	22 mm	70 mm	1,00 mm	1,00 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
180 220 080 090	S 180 220 080 090	22 mm	80 mm	0,90 mm	0,90 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
180 220 090 100	S 180 220 090 100	22 mm	90 mm	1,00 mm	1,00 mm	Straight



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
180 250 090 090	S 180 250 090 090	25 mm	90 mm	0,90 mm	0,90 mm	Straight
180 250 090 100	S 180 250 090 100	25 mm	90 mm	1,00 mm	1,00 mm	Straight
180 250 090 120	S 180 250 090 120	25 mm	90 mm	1,20 mm	1,20 mm	Straight
180 250 090 127	S 180 250 090 127	25 mm	90 mm	1,27 mm	1,27 mm	Straight
180 250 090 130	S 180 250 090 130	25 mm	90 mm	1,30 mm	1,30 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
180 070 020 080	S 180 070 020 080	7 mm	20 mm	0,80 mm	0,80 mm	Straight



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
180 070 025 060	S 180 070 025 060	7 mm	25 mm	0,60 mm	0,60 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
180 090 025 060	S 180 090 025 060	9 mm	25 mm	0,60 mm	0,60 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
180 090 080 060	S 180 090 080 060	9 mm	80 mm	0,60 mm	0,60 mm	Straight
180 090 080 090	S 180 090 080 090	9 mm	80 mm	0,90 mm	0,90 mm	Straight




Saw Blades Compatible with Chinese Power Tools




Brand Code **300**



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
300 130 090 090	S 300 130 090 090	13 mm	90 mm	0,90 mm	0,90 mm	Straight
300 130 090 100	S 300 130 090 100	13 mm	90 mm	1,00 mm	1,00 mm	Straight
300 130 090 120	S 300 130 090 120	13 mm	90 mm	1,20 mm	1,20 mm	Straight
300 130 090 130	S 300 130 090 130	13 mm	90 mm	1,30 mm	1,30 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
300 190 090 090	S 300 190 090 090	19 mm	90 mm	0,90 mm	0,90 mm	Straight
300 190 090 100	S 300 190 090 100	19 mm	90 mm	1,00 mm	1,00 mm	Straight
300 190 090 120	S 300 190 090 120	19 mm	90 mm	1,20 mm	1,20 mm	Straight
300 190 090 127	S 300 190 090 127	19 mm	90 mm	1,27 mm	1,27 mm	Straight
300 190 090 130	S 300 190 090 130	19 mm	90 mm	1,30 mm	1,30 mm	Straight




NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
300 190 100 090	S 300 190 100 090	19 mm	100 mm	0,90 mm	0,90 mm	Straight


Chinese Power Tools
Contact customer representative
for compatibility.



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
300 250 090 090	S 300 250 090 090	25 mm	90 mm	0,90 mm	0,90 mm	Straight
300 250 090 100	S 300 250 090 100	25 mm	90 mm	1,00 mm	1,00 mm	Straight
300 250 090 120	S 300 250 090 120	25 mm	90 mm	1,20 mm	1,20 mm	Straight
300 250 090 130	S 300 250 090 130	25 mm	90 mm	1,30 mm	1,30 mm	Straight



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
300 070 020 080	S 300 070 020 080	7 mm	20 mm	0,80 mm	0,80 mm	Straight



NON-STERILE REF	STERILE REF	Cutting Edge	Cutting Depth	Material Thickness	Cutting Thickness	Teeth
300 090 025 060	S 300 090 025 060	9 mm	25 mm	0,60 mm	0,60 mm	Straight



ACF Medikal Urunler Mak. San. ve Tic. Ltd. Sti.

Guzelyurt Mah. 5754 Sok. No: 3/A
Yunusemre / Manisa, Turkey, 45030

+90 236 236 28 20

info@acf.com.tr

www.acf.com.tr